



# SURGICAL TECHNIQUE GUIDE:

TIBIOTALAR/TIBIOTALOCALCANEAL ARTHRODESIS

Exclusively foot & ankle <sup>28</sup>  
**Paragon<sup>®</sup>**



**SILVERBACK<sup>™</sup>**  
ANKLE FUSION PLATING SYSTEM

# SURGICAL TECHNIQUE GUIDE:

## TIBIOTALAR/TIBIOCALCANEAL ARTHRODESIS

### Acknowledgment:

Paragon 28® would like to thank Byron Hutchinson, DPM and Mark Myerson, MD for their contribution to the development of the surgical technique guide.

## PRODUCT DESCRIPTION

The Paragon 28® Silverback™ Ankle Fusion Plating System was designed to give surgeons options for tibiotalar (TT) and tibiotalocalcaneal (TTC) arthrodesis. The hole sizing allows for Ø4.5 mm and Ø5.2 mm screws to be used for the tibia and calcaneus, while the talar screw holes allow for Gorilla® R3CON Ø3.5 mm and Ø4.2 mm screws. A Ø4.7 mm “Compact” screw is available, which was designed with a smaller thread height to help reduce insertion torque in dense bone. Additionally, single lead bone threads result in a decreased pitch differential between the locking screw head and bone threads to reduce the amount of insertion torque required to lock the screw into the plate in areas of dense bone. Precision® Guides are provided to allow for a crossing screw to be inserted outside of the plate while avoiding interference with the on-axis plate screws. Each plate has one Precision® Guide to place a tibiotalar screw, while the lateral TTC plate has two Precision® Guides—one to place a tibiotalar screw and one to place a subtalar screw. The crossing screw works in conjunction with the relatively thinner plate to help evenly distribute force across the construct and help guard against stress shielding during healing.

TABLE OF CONTENTS	PAGE #:
Implants: Plate Families .....	2
Featured Instruments .....	3-5
<b>FEATURED TECHNIQUES:</b>	
Anterior TT Plate .....	6-9
Anterior TTC Plate .....	10-13
Lateral TTC Plate .....	14-17
Lateral TT Plate .....	18
Caddy Layout .....	19
Indications, Contraindications, Warnings .....	20-21



## ANKLE FUSION PLATES

Available in Right (shown) and Left Configurations



### Anterior TT Plates

### Lateral TT Plates

### Lateral TTC Plates








































### Anterior TTC Plates



# SURGICAL TECHNIQUE GUIDE:

## ANKLE FUSION PLATING SYSTEM

### FEATURED INSTRUMENTS

	Ø3.5 mm R3CON Screws	Ø4.2 mm R3CON Screws	Ø4.5 mm SILVERBACK Screws	Ø5.2 mm SILVERBACK Screws	Ø4.7 mm SILVERBACK Compact Screws
Locking:					
Non-locking:					
Screw Lengths:	14 mm - 30 mm in 2 mm increments	14 mm - 50 mm      and      55 mm - 60 mm in 2 mm increments      in 5 mm increments			20 mm - 40 mm in 2 mm increments
Drill Size:	Ø2.4 mm 	Ø2.8 mm 	Ø3.1 mm 	Ø3.6 mm 	Ø3.6 mm 
Driver Size:	HX-10 	HX-10 	HX-15 	HX-15 	HX-15 
Locking Drill Guide Size:	Ø3.5 mm 	Ø4.2 mm 	Ø4.5 mm 	Ø4.7/Ø5.2 mm 	Ø4.7/Ø5.2 mm 
Centering Drill Guide Size:	Ø3.5 mm 	Ø4.2 mm 	Ø4.5 mm 	Ø5.2 mm 	
Compression Slot Drill Guide Size:			Ø4.5 mm 	Ø4.7/Ø5.2 mm 	Ø4.7/Ø5.2 mm 
Cone/Straight Easy Guide Size:	Ø3.5 mm 	Ø4.2 mm 	Ø4.5 mm 	Ø4.7/Ø5.2 mm 	Ø4.7/Ø5.2 mm 

Drill 

Driver 

Locking Drill Guide 

Centering Drill Guide 

Compression Slot Drill Guide 

Cone/Straight Easy Guide 

## SURGICAL TECHNIQUE GUIDE:

### ANKLE FUSION PLATING SYSTEM

#### JOINT PREPARATION INSTRUMENTATION



Bone Fenestration Perforator



Oval Burr



Barrel Burr



Curved Bone Fenestration Chisel



Straight Bone Fenestration Chisel



Cartilage Removal Tool



Curved 3 mm Osteotome



Angled Ring Curette



Straight 6 mm Osteotome



Straight Ring Curette



Curved 6 mm Osteotome



Angled Curette



Straight 12 mm Osteotome




Straight Curette




Curved 12 mm Osteotome

#### OTHER INSTRUMENTATION



Ø2.5 & Ø3.0 mm K-wire  
Hindfoot Distractor




Overdrill Guide



AO Handles




Threaded Plate Bending Bars



Flat Plate Benders



Olive Wire, Short  
(Smooth and Threaded)



Olive Wire, Long  
(Smooth and Threaded)

Ø2.0 x 200 mm K-wire (Smooth and Threaded)

Ø2.5 x 150 mm K-wire (Smooth and Threaded)




## SURGICAL TECHNIQUE GUIDE:


### ANKLE FUSION PLATING SYSTEM

#### PRECISION GUIDE INSTRUMENTATION

##### Anterior TT/TTC Precision Guide



##### Lateral TT/TTC Precision Guide



### INCISION/EXPOSURE

#### Tibiotalar Joint

A longitudinal midline incision is made over the anterior ankle, beginning approximately 10 cm proximal to the ankle joint and terminating just distal to the talonavicular joint. The incision will start approximately 1 cm lateral to the tibial crest and will course just lateral to the tibialis anterior tendon. The initial incision should penetrate skin only, but no direct tension should be placed on the skin margins until full-thickness retraction is possible. Identify the superficial peroneal nerve and retract it laterally. Continue exposure to the extensor retinaculum. Identify the extensor hallucis longus (EHL) tendon below the retinaculum and divide the retinaculum longitudinally over the EHL tendon. Care should be taken to leave the sheath of the tibialis anterior (TA) tendon intact and the retinaculum well preserved for repair at closure.



Retract the EHL tendon laterally and the TA tendon medially. Identify the neurovascular bundle and retract it laterally with the EHL tendon. Continue exposure until the anterior capsule is visualized. Perform an anterior capsulotomy via a longitudinal incision. Elevate the capsule and periosteum over the anterior tibia and talus to expose the anterior ankle joint, the tibial plafond, the medial and lateral gutters and the anterior and dorsal talus. Remove any osteophytes on the tibia and talus to allow for exposure to the ankle joint and facilitate entry of instrumentation for cartilage removal. All osteophytes must be removed from the anterior ankle to facilitate application of the plate.

### JOINT PREPARATION


Preparation of the tibiotalar joint can be performed using the provided joint preparation instrumentation. A pin distractor is provided to allow for space and visualization during joint preparation, to be used with provided Ø2.5 mm K-wires. The pin distractor and K-wires are placed through the lateral aspect of the incision with appropriate soft tissue retraction, in an anterior to posterior direction. Following cartilage removal, it is advised to penetrate the subchondral plate with the subchondral drill, burrs and/or bone fenestration chisels to promote healing.

### PROVISIONAL FIXATION

Align the ankle joint. The foot and ankle should be positioned such that the ankle is neutral with respect to dorsiflexion and plantarflexion. The foot should be in approximately 5-10° of external rotation and 5° of hindfoot valgus. With the foot and ankle held in this alignment, use one or two Ø2.0 mm K-wires to temporarily fixate the tibiotalar joint, per surgeon preference.



### PLATE PLACEMENT




Frequently there is an irregular surface remaining over the anterior joint which requires debridement to a smooth surface with either a rongeur or saw. Retrieve the appropriate anterior TT plate based on the patient's anatomy. To position the plate, palpate the medial and lateral margins of the talus and center the talar portion of the plate. Ensure that the proximal plate is midline or just lateral to midline.



Secure the plate to the tibiotalar joint using a long olive wire in the most proximal circular hole on the tibia and a short olive wire in the medial talar neck screw hole. Confirm plate position using fluoroscopy.

### PERMANENT FIXATION - PLATE SCREWS

**NOTE:** The talar screw holes accept Ø3.5 mm or Ø4.2 mm non-locking and locking screws. Ø4.2 mm screws are recommended for this area, except in the case of a small patient. The use of Ø4.2 mm screws is demonstrated in this technique. When using Ø3.5 mm screws, use the appropriate instrumentation as described on page 3.



Retrieve the Ø4.2 mm threaded drill guide and thread into the lateral talar body screw hole. Drill, using the Ø2.8 mm drill.



Remove the drill guide and measure screw length using the depth gauge. Confirm screw projection and length using the depth gauge under fluoroscopy (not shown). Insert the selected screw size into the plate hole using the provided driver and handle. Do not fully tighten screw until the second talar screw is secure, to prevent toggling of the plate.



Remove the olive wire in the talar neck hole. Insert a second Ø4.2 mm screw into the medial talar body hole using the same procedure previously described. Complete tightening and seating of both talar body screws.



## SURGICAL TECHNIQUE GUIDE:


### ANTERIOR TIBIOTALAR ARTHRODESIS

#### PERMANENT FIXATION - TIBIAL COMPRESSION SCREW

**NOTE:** The tibial screw holes accept Ø4.5 mm or Ø5.2 mm locking or non-locking screws. A laser etched dot on the plate indicates the plate holes that accept the Ø4.5 mm and Ø5.2 mm screws. The technique demonstrates the use of the Ø4.5 mm screws. When using the Ø5.2 mm screws, use the appropriate instrumentation as described on page 3.



Retrieve the Ø4.5 mm oblong compression slot drill guide and insert into the tibial compression slot with the arrow pointing towards the tibiotalar joint. Drill, using a Ø3.1 mm drill through the compression slot drill guide.




Remove the drill guide and measure screw length using the depth gauge. Insert a Ø4.5 mm non-locking screw. Do not fully tighten the screw until the tibiotalar crossing screw has been inserted. This allows for travel of the compression screw distally in the slot during final tightening of the crossing screw. Remove the olive wire from the tibial screw hole and provisional fixation across the tibiotalar joint.

#### PERMANENT FIXATION - PRECISION GUIDED CROSSING SCREW

Retrieve the corresponding anterior TT Precision Guide and set screw. Insert the provided set screw into the hole adjacent to the boss on the undersurface of the Precision Guide arm.

Align the boss on the underside of the Precision Guide into the smaller of the two central holes on the ankle fusion plate to correctly orient the Precision Guide, while aligning the set screw over the larger of the two central holes. Rotate the set screw clockwise to secure the Precision Guide to the ankle fusion plate.

**NOTE:** Alternatively, the Precision Guide can be placed on the plate following plating screw fixation of the plate to bone. Placement of a fully threaded crossing screw would take place following fully seating the tibial non-locking screw in the compression slot.






## SURGICAL TECHNIQUE GUIDE:


### ANTERIOR TIBIOTALAR ARTHRODESIS

#### PERMANENT FIXATION - PRECISION GUIDED CROSSING SCREW



Insert the K-wire sleeve into the one of the two parallel trajectories of the Precision Guide, depending on desired starting point across the arthrodesis site. The posterior trajectory is recommended to provide more even fixation across the joint, if the position is appropriate for the patient.

**NOTE:** K-wire Tubes for the Precision Guide are available in Ø1.6 mm and Ø2.3 mm, allowing for Ø5.5 mm or Ø7.0 mm Monster Screws to be used. Partially threaded and fully threaded screw options are available for each screw diameter, per surgeon preference. The use of Ø7.0 mm Monster Screw is demonstrated in this technique. When using the Ø5.5 mm Monster Screw, use the appropriate instrumentation.



Retrieve the Ø2.3 mm K-wire for a Ø7.0 mm partially threaded Monster Screw. Insert the K-wire through the K-wire sleeve and across the arthrodesis site.

Confirm crossing screw trajectory using fluoroscopy. The ideal position for the trajectory is just proximal to the lateral process of the talus.

Remove the K-wire sleeve from the Ø2.3 mm K-wire. Remove the Precision Guide by turning the set screw counterclockwise to detach the Precision Guide from the plate. Slide the Precision Guide and K-wire sleeve off the K-wire.

Retrieve the countersink for the Ø7.0 mm headed Monster Screw. Rotate the countersink clockwise over the K-wire to remove adequate bone to seat the screw head. Measure screw length using the depth gauge (not shown).

**NOTE:** If soft bone is apparent, it is advised not to countersink in this area to allow for better screw purchase. Monster Screw washers are available; however, the steep screw angle may cause prominence of the washer.

**TIP:** If minimal compression or poor bone purchase occurs with this crossing screw, remove the partially threaded screw and guide wire. Obtain compression by fully tightening the non-locking screw in the compression slot and then place a fully threaded screw via the Precision Guide.



## SURGICAL TECHNIQUE GUIDE:

### ANTERIOR TIBIOTALAR ARTHRODESIS


#### PERMANENT FIXATION - PRECISION GUIDED CROSSING SCREW



Drill over the K-wire using the  $\varnothing 4.6$  mm drill for the  $\varnothing 7.0$  mm headed Monster Screw.




Insert a  $\varnothing 7.0$  mm headed Monster Screw using the provided driver. Confirm screw length and placement using fluoroscopy.




Fully seat the non-locking screw in the tibia compression slot.


#### PERMANENT FIXATION - PLATE SCREWS



Insert two  $\varnothing 4.2$  mm screws into the talar neck holes using the same procedure previously described.



Insert the remaining tibial screws using the technique previously described for  $\varnothing 4.5$  mm screws.



Confirm screw position and placement using fluoroscopy.

#### CLOSURE

Proceed to incision closure or concomitant procedures at this time.



## SURGICAL TECHNIQUE GUIDE:

# ANTERIOR TIBIOTALOCALCANEAL ARTHRODESIS



## INCISION/EXPOSURE

### Tibiotalar Joint

An incision for the tibiotalar joint is made as described on page 6 for the anterior TT arthrodesis surgical technique.

### Subtalar Joint

A longitudinal incision is made starting at the distal aspect of the lateral malleolus over the subtalar joint extending toward the 4<sup>th</sup> metatarsal base and terminating at the calcaneocuboid joint. Continue exposure through the subcutaneous tissue, with care being taken to identify and retract the anterior branch of the sural nerve. The extensor hallucis brevis muscle is reflected distally to expose the sinus tarsi and posterior facet of the subtalar joint. Dissection of the fat pad out of the sinus tarsi should be performed, with reflection of the tissue dorsally.




## JOINT PREPARATION

Preparation of the tibiotalar joint can be performed using the provided joint preparation instrumentation. The medial and lateral gutters should be cleared at this time. Remove cartilage from the posterior, middle and anterior facets of the subtalar joint. A pin distractor is provided to allow for space and visualization during joint preparation, to be used with provided Ø2.5 mm K-wires. It is advised to penetrate the subchondral plate with the subchondral drill, burrs and/or chisels to promote healing.

## PROVISIONAL FIXATION

Align the ankle joint. The foot and ankle should be positioned such that the ankle is neutral with respect to dorsiflexion and plantarflexion. The foot should be in approximately 5-10° of external rotation and 5° of hindfoot valgus. With the foot and ankle held in this alignment, use Ø2.0 mm K-wires to temporarily fixate the tibiotalar joint and subtalar joint, per surgeon preference.




## PLATE POSITIONING

Resection of osteophytes or prominent bone may be necessary to ensure proper plate fit. Retrieve the appropriate anterior TTC plate based on the patient's anatomy.

To position the plate, palpate the medial and lateral margins of the talus and center the talar portion of the plate. Ensure that the proximal plate is midline or just lateral to midline.

Secure the plate to the tibiotalar joint using a long olive wire in the most proximal circular hole on the tibia and a short olive wire in the medial talar neck screw hole. Confirm plate position using fluoroscopy.




## SURGICAL TECHNIQUE GUIDE:

### ANTERIOR TIBIOTALOCALCANEAL ARTHRODESIS

#### PERMANENT FIXATION - TALAR NECK AND TIBIAL COMPRESSION SCREWS



Insert talar neck screws as described on page 7. Insert the tibial compression screw as described on page 8. Remove the K-wire serving as provisional fixation across the tibiotalar joint and the olive wire in the tibia.



#### PERMANENT FIXATION - TALAR NECK AND TIBIAL SCREWS


Assemble and attach the Precision Guide to the anterior TTC plate as described on page 8.



Insert the tibiotalar crossing screw as described on page 9 and 10. Confirm screw length and placement using fluoroscopy. Remove the Ø2.3 mm K-wire.

#### PERMANENT FIXATION - TIBIAL SCREWS

Insert tibial screws as described on page 9.







## SURGICAL TECHNIQUE GUIDE:

### ANTERIOR TIBIOTALOCALCANEAL ARTHRODESIS

#### PERMANENT FIXATION - SUBTALAR JOINT CROSSING SCREW





A Ø7.0 headed Monster Screw is used in this plate hole. Retrieve the subtalar joint wire guide for a Ø7.0 mm headed Monster Screw. Mate the tip of the wire guide into the socket of the plate.




Insert a Ø2.3 mm K-wire through the wire guide across the subtalar joint, aiming for the central posterior aspect of the calcaneus, just above the weight bearing surface.

Remove the K-wire serving as provisional fixation across the subtalar joint.

Measure screw length using the cannulated depth gauge.



Drive the guide wire posteriorly until it exits the skin. Retrieve a parallel K-wire guide from the Monster Screw System instrument caddy. Slide the central, isolated hole of the parallel K-wire guide over the Ø2.3 mm K-wire.




Depending on screw type used (headed vs. headless), place a second Ø2.3 mm K-wire a desired distance proximal to the initial Ø2.3 mm K-wire.

## SURGICAL TECHNIQUE GUIDE:


### ANTERIOR TIBIOTALOCALCANEAL ARTHRODESIS

#### PERMANENT FIXATION – SUBTALAR JOINT CROSSING SCREW

**NOTE:** A partially threaded screw is recommended to be placed in the more proximal position, to allow for the threads of the screw to engage the denser bone in the talus and create compression across the subtalar joint. The second screw placed in the plate is intended to be a fully threaded screw to hold in the compression created from the first screw.




Placement of a partially threaded  $\varnothing 7.0$  mm Monster Screw through the posterior calcaneus into the talus is performed as described on page 10. Remove the  $\varnothing 2.3$  mm K-wire. Confirm screw placement using fluoroscopy.



**NOTE:** Countersinking is not necessary, as the  $\varnothing 7.0$  mm Monster Screw is seated within the plate.

Drill over the K-wire using the  $\varnothing 4.6$  mm drill for the  $\varnothing 7.0$  mm headed Monster Screw.




Insert a  $\varnothing 7.0$  mm headed, fully threaded Monster Screw using the provided driver.

Confirm screw length and placement using fluoroscopy.  
Remove the  $\varnothing 2.3$  mm K-wire.

#### CLOSURE

Proceed to incision closure or concomitant procedures at this time.





## SURGICAL TECHNIQUE GUIDE:

### LATERAL TIBIOTALOCALCANEAL ARTHRODESIS

#### INCISION/EXPOSURE

A lateral incision is made over the posterior half of the fibula, beginning approximately 10 cm proximal to the tip of the fibula extending distally to the plantar aspect of the calcaneus.

Identify the sural nerve and retract it posteriorly. Continue dissection to the fibula. While retracting the peroneal tendons and sural nerve, a transverse fibular osteotomy is performed by beveling the saw from proximal lateral to distal medial to avoid a sharp bony prominence above the plate. Transect the syndesmotic and lateral ankle ligaments to free the fibula from adjacent soft tissues. Resect the fibula and retain for bone graft, if desired. Alternatively, if available, a reamer can be used to remove the fibula, while harvesting the reamed bone for graft material.

Elevate the anterior joint capsule and nearby periosteum to assess the anterior tibiotalar joint articulation. Minimal dissection of the talar neck is recommended to avoid devascularization of this bone. Remove any anterior osteophytes that may interfere with joint reduction. Elevate the posterior soft tissues using a periosteal elevator to allow for retractors to be placed anterior to and posterior to the tibiotalar joint. Elevate the extensor digitorum brevis muscle belly to expose the subtalar joint. Release the lateral ligaments around the subtalar joint including the talocalcaneal intraosseous ligament to allow for appropriate distraction. Dissect the fat pad out of the sinus tarsi to allow for appropriate visualization.




#### JOINT PREPARATION

Preparation of the tibiotalar joint can be performed using the provided joint preparation instrumentation. A pin distractor is provided to allow for space and visualization during joint preparation, to be used with provided Ø2.5 mm K-wires. It is advised to penetrate the subchondral plate with the subchondral drill, burrs and/or chisels to promote healing. A medial arthrotomy may be required and performed to allow for exposure and joint preparation of the medial gutter of the ankle joint.

Remove cartilage from the posterior, middle and anterior facets of the subtalar joint, . Perform subchondral plate penetration to these joints to promote healing.

#### PROVISIONAL FIXATION

Align the ankle joint. The foot and ankle should be positioned such that the ankle is neutral with respect to dorsiflexion and plantarflexion. The foot should be in approximately 5-10° of external rotation and 5° of hindfoot valgus. With the foot and ankle held in this alignment, use multiple Ø2.0 mm K-wires to temporarily fix the tibiotalar and subtalar joints as shown. Provisional fixation wires should be placed in the anterolateral to posteromedial direction for the tibiotalar joint, and posterior to anterior direction for the subtalar joint.



## SURGICAL TECHNIQUE GUIDE:

### LATERAL TIBIOTALOCALCANEAL ARTHRODESIS

#### PLATE POSITIONING AND PROVISIONAL FIXATION



Resection of osteophytes or prominent bone such as the lateral talar process may be necessary to ensure proper plate fit. Retrieve the appropriate lateral TTC plate based on the patient's anatomy. The plate should be positioned such that the proximal aspect is centered from anterior to posterior on the tibia, the central talar hole is centered on the body of the talus and the posterior calcaneal holes align just inferior to the superior surface of the calcaneus. If necessary, a saw can be used to scrape the bone surfaces smooth to fit the contour of the plate.




Secure the plate to the lateral aspect of the tibiotalar joint using a long olive wire in a circular tibial hole and a short olive wire in the talus and calcaneus, as shown. Confirm plate position using fluoroscopy.

**NOTE:** The small hole in the talar portion of the plate is directly in line with the center of the tibial portion of the plate. Talar position can be verified using this hole for guidance. Confirm trajectory of subtalar screw using fluoroscopy.


#### PERMANENT FIXATION - PLATE SCREWS

**NOTE:** The talar screw holes accept Ø3.5 mm or Ø4.2 mm non-locking and locking screws. The use of Ø4.2 mm screws is demonstrated in this technique. When using the Ø3.5 mm screws, use the appropriate instrumentation as described on page 3.




Retrieve the Ø4.2 mm threaded drill guide and thread into a talar body screw hole. Drill, using the Ø2.8 mm drill.

**NOTE:** Screw placement within the talus may be limited in cases of diseased or eroded tali.



Remove the drill guide and measure screw length using the depth gauge (not shown). Insert the selected screw into the plate hole using the provided driver and handle. Do not fully tighten the screw until the second talar screw is secure, to prevent toggling of the plate.



Remove the olive wire from the talus. Insert remaining locking or non-locking talar screws, per surgeon preference. It is recommended to drill using the threaded drill guide or the straight end of the fast guide, in order to achieve on-axis screw trajectories, which allows for placement of a crossing screw via the Precision Guide.



## SURGICAL TECHNIQUE GUIDE:


### LATERAL TIBIOTALOCALCANEAL ARTHRODESIS

#### PERMANENT FIXATION - PLATE SCREWS

**NOTE:** The tibial and calcaneal screw holes accept Ø4.5 or Ø5.2 mm screws. A laser etched dot on the plate indicates the plate holes that accept the Ø4.5 mm and Ø5.2 mm screws. The instructions provided below are for Ø4.5 mm screws. When using the Ø5.2 mm screws, use the appropriate instrumentation as described on page 3.



Retrieve the Ø4.5 mm oblong compression slot drill guide and insert into the tibial compression slot with the arrow pointing toward the tibiotalar joint. Drill, using a Ø3.1 mm drill through the compression slot drill guide.



Remove the drill guide and measure screw length using the depth gauge. Insert a Ø4.5 mm non-locking screw but do not fully seat. Remove the K-wire serving as provisional fixation as well as the tibial olive wire prior to fully seating the Ø4.5 mm screw.




Seat the Ø4.5 mm non-locking screw in the tibia.

#### PERMANENT FIXATION - TIBIOTALAR PRECISION GUIDED CROSSING SCREW

Retrieve the lateral TT precision guide arm and set screw. The Precision Guide is attached to the plate as described on page 8.

**NOTE:** K-wire Tubes for the Precision Guide are available in Ø1.6 mm and Ø2.3 mm, allowing for Ø5.5 mm or Ø7.0 mm Monster Screws to be used. Partially threaded and fully threaded screw options are available for each screw diameter, per surgeon preference. The use of a Ø7.0 mm Monster Screw is demonstrated in this technique. When using the Ø5.5 mm Monster Screw, use the appropriate instrumentation.




Insert a Ø7.0 mm partially threaded headed Monster Screw using the technique described on pages 8 and 9, and remove the provisional fixation wire across the tibiotalar joint and the calcaneal olive wire. Confirm screw length and placement using fluoroscopy. Remove the Ø2.3 mm K-wire.


## SURGICAL TECHNIQUE GUIDE:

### LATERAL TIBIOTALOCALCANEAL ARTHRODESIS

#### PERMANENT FIXATION - SUBTALAR PRECISION GUIDED CROSSING SCREW




Retrieve the K-wire guide and thread into the lateral subtalar Precision Guide arm. Insert a Ø2.3 mm K-wire through the wire guide across the subtalar joint.




Using the method outlined on pages 13 and 14, place a partially threaded Ø7.0 mm Monster Screw across the subtalar joint. Remove provisional fixation across the subtalar joint prior to fully seating the Monster Screw.

#### PERMANENT FIXATION - TIBIA AND CALCANEAL PLATE SCREWS



Screw placement in the remaining tibial and calcaneal plate holes is completed at this time using Ø4.5 mm or Ø5.2 mm locking or non-locking screws, per surgeon preference. It is recommended to drill using the threaded drill guide or the straight end of the fast guide, in order to achieve an on-axis screw trajectory.

**TIP:** If off-axis plate screw trajectories are desired, it is recommended to angle the screws away from the laser markings on the plate, which indicate Precision-guided crossing screw trajectories.



#### CLOSURE

Proceed to incision closure or concomitant procedures at this time.



## SURGICAL TECHNIQUE GUIDE:

### LATERAL TIBIOTALAR ARTHRODESIS

#### INCISION/EXPOSURE

A lateral incision is made over the posterior half of the fibula, beginning approximately 10 cm proximal to the tip of the fibula and curving anterior distally toward the 4<sup>th</sup> metatarsal, just past the tip of the fibula.

Identify the sural nerve and retract it posteriorly. Continue dissection to the fibula. While retracting the peroneal tendons and sural nerve, a transverse fibular osteotomy is performed by beveling the saw from proximal lateral to distal medial to avoid a sharp bony prominence above the plate. Transect the syndesmotomic and lateral ankle ligaments to free the fibula from adjacent soft tissues. Resect the fibula and retain for bone graft, if desired.


Elevate the anterior joint capsule and nearby periosteum to access the anterior tibiotalar joint articulation. Minimal dissection of the talar neck is recommended to avoid devascularization of this bone. Remove any anterior osteophytes that may interfere with joint reduction. Elevate the posterior soft tissues using a periosteal elevator to allow for retractors to be placed anterior to and posterior to the tibiotalar joint.



#### JOINT PREPARATION

Preparation of the tibiotalar joint can be performed as described on page 15, or per surgeon preference using the provided joint preparation instrumentation.

#### PERMANENT FIXATION



Steps for fixation of a lateral TT plate, screws and crossing screw can be performed as described previously for the lateral TTC plate.

#### CLOSURE

Proceed to incision closure or concomitant procedures at this time.

# SILVERBACK™ ANKLE FUSION CADDY


## Silverback™ Anterior Plate Caddy

Anterior TT and TTC plates and corresponding Precision Guides are located within the Anterior Plate Caddy.





## Silverback™ Lateral Plate Caddy

Lateral TT and TTC plates and corresponding Precision Guides are located within the Lateral Plate Caddy.




## Silverback™ K-wire and Olive Wire Caddy

Smooth and threaded K-wires and olive wires and a ruler are located within the K-wire and Olive Wire Caddy.



## Silverback™ Instrument Tray

All drill guides, drills, overdrills, taps, drivers, forceps and a depth gauge are located within the Silverback™ instrument tray.



## Silverback™ Case Base

Handles, plate bending instrumentation and joint preparation instrumentation including curettes, osteotomes, chisels and a cartilage removal tool are located at the bottom of the Silverback™ Instrument Case.

## Silverback™ Screw Caddy

The Silverback™ screw length options for locking and non-locking screws are as follows:

3.5 mm	2 mm increments, 14-30 mm	
4.2 mm	2 mm increments, 14-50 mm	
4.2 mm	5 mm increments, 55-60 mm	
4.5 mm	2 mm increments, 14-50 mm	
4.5 mm	5 mm increments, 55-60 mm	
5.2 mm	2 mm increments, 14-50 mm	
5.2 mm	5 mm increments, 55-60 mm	

The Silverback™ compact screw length options are as follows:

4.7 mm	2 mm increments, 20-40 mm	
--------	---------------------------	--



# SURGICAL TECHNIQUE GUIDE: INDICATIONS, CONTRAINDICATIONS, AND WARNINGS

Refer to [www.paragon28.com/ifus](http://www.paragon28.com/ifus) for the complete and most current instructions for use document.

## INDICATIONS FOR USE (GORILLA®)

The Baby Gorilla®/Gorilla® Bone Plates and Bone Screws of the Baby Gorilla®/Gorilla® Plating System are indicated for use in stabilization and fixation of fractures or osteotomies; intra and extra articular fractures, joint depression, and multi-fragmentary fractures; revision procedures, joint fusion and reconstruction of small bones of the toes, feet and ankles including the distal tibia, talus, and calcaneus, as well as the fingers, hands, and wrists. The system can be used in both adult and pediatric patients. Specific examples include:

### Forefoot:

- Arthrodesis of the first metatarsalcuneiform joint (Lapidus Fusion)
- Metatarsal or phalangeal fractures and osteotomies
- Lesser metatarsal shortening osteotomies (e.g. Weil)
- Fifth metatarsal fractures (e.g. Jones Fracture)

### Mid/Hindfoot:

- LisFranc Arthrodesis and/or Stabilization
- 1<sup>st</sup> (Lapidus), 2<sup>nd</sup>, 3<sup>rd</sup>, 4<sup>th</sup>, and 5<sup>th</sup> Tarsometatarsal (TMT) Fusions
- Intercuneiform Fusions
- Navicular-Cuneiform (NC) Fusion
- Talo-Navicular (TN) Fusion
- Calcaneo-Cuboid (CC) Fusion
- Subtalar Fusion
- Medial Column Fusion
- Cuneiform Fracture
- Cuboid Fracture
- Navicular Fracture

### Ankle:

- Lateral Malleolar Fractures
- Syndesmosis Injuries
- Medial Malleolar Fractures and Osteotomies
- Bi-Malleolar Fractures
- Tri-Malleolar Fractures
- Posterior Malleolar Fractures
- Distal Anterior Tibia Fractures
- Vertical Shear Fractures of the Medial Malleolus
- Pilon Fractures
- Distal Tibia Shaft Fractures
- Distal Fibula Shaft Fractures
- Distal Tibia Periarticular Fractures
- Medial Malleolar Avulsion Fractures
- Lateral Malleolar Avulsion Fractures
- Tibiotalocalcaneal Joint Arthrodesis
- Tibiotalar Joint Arthrodesis
- Tibiocalcaneal Arthrodesis
- Supramalleolar Osteotomy
- Fibular Osteotomy

### First metatarsal osteotomies for hallux valgus correction including:

- Opening base wedge osteotomy
- Closing base wedge osteotomy
- Crescentic Osteotomy
- Proximal Osteotomy (Chevron and Rotational Oblique)
- Distal Osteotomy (Chevron/Austin)

### Arthrodesis of the first metatarsophalangeal joint (MTP) including:

- Primary MTP Fusion due to hallux rigidus and/or hallux valgus
- Revision MTP Fusion
- Revision of failed first MTP Arthroplasty implant

### Flatfoot:

- Lateral Column Lengthening (Evans Osteotomy)
- Plantar Flexion Opening Wedge Osteotomy of the Medial Cuneiform (Cotton Osteotomy)
- Calcaneal Slide Osteotomy

### Charcot:

- Medial column fusion (talus, navicular, cuneiform, metatarsal) for neuropathic osteoarthropathy (Charcot)
- Lateral column fusion (calcaneus, cuboid, metatarsal) for neuropathic osteoarthropathy (Charcot)

In addition, the non-locking, titanium screws and washers are indicated for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, fracture repair and fracture fixation, appropriate for the size of the device.

## CONTRAINDICATIONS

Use of the Baby Gorilla®/Gorilla® Plating System is contraindicated in cases of inflammation, cases of active or suspected sepsis/infection and osteomyelitis; or in patients with certain metabolic diseases.

All applications that are not defined by the indications are contraindicated. In addition, surgical success can be adversely affected by:

- Acute or chronic infections, local or systemic
- Vascular, muscular or neurological pathologies that compromise the concerned extremity
- All concomitant pathologies that could affect the function of the implant
- Osteopathies with reduced bone substance that could affect the function of the implant
- Any mental or neuromuscular disorder that could result in an unacceptable risk of failure at the time of fixation or complications in post-operative treatment
- Known or suspected sensitivity to metal
- Corpulence; an overweight or corpulent patient can strain the implant to such a degree that stabilization or implant failure can occur
- Whenever the use of the implant comes into conflict with the anatomical structures of physiological status

Other medical or surgical pre-conditions that could compromise the potentially beneficial procedure, such as:

- The presence of tumors
- Congenital abnormalities
- Immunosuppressive pathologies
- Increased sedimentation rates that cannot be explained by other pathologies
- Increased leukocyte (WBC) count
- Pronounced left shift in the differential leukocyte count

## POTENTIAL COMPLICATIONS AND ADVERSE REACTIONS

In any surgical procedure, the potential for complications and adverse reactions exist. The risks and complications with these implants include:

- Loosening, deformation or fracture of the implant
- Acute post-operative wound infections and late infections with possible sepsis
- Migration, subluxation of the implant with resulting reduction in range of movement
- Fractures resulting from unilateral joint loading
- Thrombosis and embolism
- Wound hematoma and delayed wound healing
- Temporary and protracted functional neurological perturbation
- Tissue reactions as the result of allergy or foreign body reaction to dislodged particles
- Corrosion with localized tissue reaction and pain
- Pain, a feeling of malaise or abnormal sensations due to the implant used
- Bone loss due to stress shielding

All possible complications listed here are not typical of Paragon 28®, Inc. products but are in principle observed with any implant. Promptly inform Paragon 28®, Inc. as soon as complications occur in connection with the implants or surgical instruments used. In the event of premature failure of an implant in which a causal relationship with its geometry, surface quality or mechanical stability is suspected, please provide Paragon 28®, Inc. with the explant(s) in a cleaned, disinfected and sterile condition. Paragon 28®, Inc. cannot accept any other returns of used implants. The surgeon is held liable for complications associated with inadequate asepsis, inadequate preparation of the osseous implant bed in the case of implants, incorrect indication or surgical technique or incorrect patient information and consequent incorrect patient behavior.

## WARNINGS AND PRECAUTIONS

- Re-operation to remove or replace implants may be required at any time due to medical reasons or device failure. If corrective action is not taken, complications may occur.
- Use of an undersized plate or screw in areas of high functional stresses may lead to implant fracture and failure.
- Plates and screws, wires, or other appliances of dissimilar metals should not be used together in or near the implant site.
- The implants and guide wires are intended for single use only.
- Instruments, guide wires and screws are to be treated as sharps.
- Do not use other manufacturer's instruments or implants in conjunction with the Baby Gorilla®/Gorilla® Plating System.
- If a stainless steel Gorilla® R3LEASE™ Screw is used, it may only be used standalone.
- The device should only be used in pediatric patients where the growth plates have fused or in which active growth plates will not be crossed by the system implants or instrumentation.

## MR SAFETY INFORMATION

The Baby Gorilla®/Gorilla® Plating System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Baby Gorilla®/Gorilla® Plating System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

# SURGICAL TECHNIQUE GUIDE: INDICATIONS, CONTRAINDICATIONS, AND WARNINGS

Refer to [www.paragon28.com/ifus](http://www.paragon28.com/ifus) for the complete and most current instructions for use document.

## INDICATIONS FOR USE (MONSTER®)

The Monster® Screw System is indicated for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, ligament fixation, fracture repair and fracture fixation, appropriate for the size of the device. Specific examples include:

### Fractures and Osteotomies

- Fractures of the tarsals, metatarsals and other fractures of the foot (i.e. LisFranc)
- Avulsion fractures and fractures of the 5th metatarsal (i.e. Jones Fracture)
- Talar fractures
- Ankle fractures
- Navicular fractures
- Fractures of the fibula, malleolus, and calcaneus
- Metatarsal and phalangeal osteotomies
- Weil osteotomy
- Calcaneal osteotomy

### Hallux Valgus Correction

- Fixation of osteotomies (i.e. Akin, Scarf, Chevron)
- Interphalangeal (IP) arthrodesis
- Proximal, midshaft, or distal osteotomy
- Lapidus arthrodesis

### Arthrodesis/Deformity Correction

- 1<sup>st</sup> MTP arthrodesis
- Metatarsal deformity correction
- Tarsometatarsal joint arthrodesis
- Naviculocuneiform joint arthrodesis
- Talonavicular arthrodesis
- Subtalar joint arthrodesis
- Triple arthrodesis
- Medial column arthrodesis
- Subtalar joint distraction arthrodesis
- Ankle arthrodesis
- Lateralizing calcaneal osteotomy
- Lateral column lengthening
- Hammertoe

### Fusion resulting from neuropathic osteoarthopathy (Charcot) such as:

- Medial and lateral column
- Subtalar, talonavicular, and calcaneocuboid

## CONTRAINDICATIONS

Use of the Monster® Screw System is contraindicated in cases of inflammation, cases of active or suspected sepsis / infection and osteomyelitis; or in patients with certain metabolic diseases.

All applications that are not defined by the indications are contraindicated. In addition, surgical success can be adversely affected by:

- Acute or chronic infections, local or systemic
- Vascular, muscular or neurological pathologies that compromise the concerned extremity
- All concomitant pathologies that could affect the function of the implant
- Osteopathies with reduced bone substance that could affect the function of the implant
- Any mental or neuromuscular disorder that could result in an unacceptable risk of failure at the time of fixation or complications in post-operative treatment
- Known or suspected sensitivity to metal
- Corpulence; an overweight or corpulent patient can strain the implant to such a degree that stabilization or implant failure can occur
- Whenever the use of the implant comes into conflict with the anatomical structures of physiological status

Other medical or surgical pre-conditions that could compromise the potentially beneficial procedure, such as:

- The presence of tumors
- Congenital abnormalities
- Immunosuppressive pathologies
- Increased sedimentation rates that cannot be explained by other pathologies
- Increased leukocyte (WBC) count
- Pronounced left shift in the differential leukocyte count

## POTENTIAL COMPLICATIONS AND ADVERSE REACTIONS

In any surgical procedure, the potential for complications and adverse reactions exist. The risks and complications with these implants include:

- Loosening, deformation or fracture of the implant
- Acute post-operative wound infections and late infections with possible sepsis
- Migration, subluxation of the implant with resulting reduction in range of movement
- Fractures resulting from unilateral joint loading
- Thrombosis and embolism
- Wound hematoma and delayed wound healing
- Temporary and protracted functional neurological perturbation
- Tissue reactions as the result of allergy or foreign body reaction to dislodged particles.
- Corrosion with localized tissue reaction and pain
- Pain, a feeling of malaise or abnormal sensations due to the implant used
- Bone loss due to stress shielding

All possible complications listed here are not typical of Paragon 28®, Inc. products but are in principle observed with any implant. Promptly inform Paragon 28® as soon as complications occur in connection with the implants or surgical instruments used. In the event of premature failure of an implant in which a causal relationship with its geometry, surface quality or mechanical stability is suspected, please provide Paragon 28® with the explant(s) in a cleaned, disinfected and sterile condition. Paragon 28® cannot accept any other returns of used implants. The surgeon is held liable for complications associated with inadequate asepsis, inadequate preparation of the osseous implant bed in the case of implants, incorrect indication or surgical technique or incorrect patient information and consequent incorrect patient behavior.

## WARNINGS AND PRECAUTIONS


- Re-operation to remove or replace implants may be required at any time due to medical reasons or device failure. If corrective action is not taken, complications may occur.
- Use of an undersized screw in areas of high functional stresses may lead to implant fracture and failure.
- Plates and screws, wires, or other appliances of dissimilar metals should not be used together in or near the implant site.
- The implants and guide wires are intended for single use only. Re-use may cause product failure and could lead to disease transmission.
- Instruments, guide wires and screws are to be treated as sharps.
- **Do not use other manufacturer's instruments or implants in conjunction with the Monster® Screw System.**

## MR SAFETY INFORMATION

The Monster® Screw System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Monster® Screw System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

## NOTES





**SILVERBACK™**  
ANKLE FUSION PLATING SYSTEM



PATENTED, DESIGNED & EXCLUSIVELY DISTRIBUTED BY

Exclusively foot & ankle **28**  
**Paragon®**


[www.PARAGON28.com](http://www.PARAGON28.com)

**P51-STG-0011 RevA**

™Trademarks and ®Registered Marks of Paragon 28®, Inc.

© Copyright 2019 Paragon 28®, Inc. All rights reserved.

Patents: [www.paragon28.com/patents](http://www.paragon28.com/patents)

Paragon 28, Inc.   
14445 Grasslands Dr.  
Englewood, CO 80112 USA  
(855) 786-2828

**DISCLAIMER**

The purpose of the SILVERBACK™ Ankle Fusion Plating System Surgical Technique Guide is to demonstrate the optionality and functionality of the SILVERBACK™ Ankle Fusion Plating System and Gorilla® R3CON Plating System. Although variations in placement and use of the SILVERBACK™ Ankle Fusion Plating System can be performed, the fixation options demonstrated in this technique were chosen to demonstrate the functionality of the system and for simplicity of explanation. Other uses for the SILVERBACK™ Ankle Fusion Plating System can be employed, appropriate for the size of the device. Federal law (U.S.A.) restricts this device to sale and use by, or on order of, a physician.