Baby Gorilla®/Gorilla® Plating System

BEFORE USING PRODUCT, READ THE FOLLOWING IMPORTANT INFORMATION

A FULL SYMBOLS GLOSSARY CAN BE FOUND AT:
www.paragon28.com/resources

Please check the website, www.paragon28.com/ifus, for the most current instructions for use document.

This booklet is designed to assist in using the Baby Gorilla®/Gorilla® Plating System. It is not a reference for surgical techniques.

CAUTION

Federal Law (USA) restricts this device to sale and use by, or on the order of, a physician.

GENERAL DESCRIPTION

The Baby Gorilla®/Gorilla® Plating System is comprised of plates, screws and washers used for bone fixation and stabilization. The bone plates are available in varying configurations (including, but not limited to, straight, curved, dog bone, rectangular, rhombus, 'T' plates, slanted 'T' plates, 'L' plates and ribbon plates) and varying lengths, which are attached to the bone using screw fixation. These plates are attached to bone using 2.0-4.2 mm diameter titanium self-tapping screws; screw diameter is dependent upon plate thickness and plate option. The screws will be available in both standard (locking) and lag design (non-locking) with a hex drive head feature. The plate screw holes are either threaded or non-threaded. Threaded plate screw holes can accept both standard (locking) screws with threaded screw heads and lag design (non-locking) screws with non-threaded screw head. Non-threaded plate screw holes can accept lag design (non-locking) screws only.

Available plates, screws, washers and instrumentation are packaged as a single system. The system instruments include guide wires, guide wire drill guide, drill bits, drill guides, tissue protectors, depth gauges, distractors, compressors, screwdriver shafts and driver handles. These instruments are used to facilitate the placement of the plates and screws.

INTENDED/ EXPECTED CLINICAL BENEFITS

Intended/ Expected Clinical Benefits of the Baby Gorilla®/ Gorilla® Plating System include:

- Successful Implantation of Baby Gorilla®/ Gorilla® Devices
- Bone healing

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The Summary of Safety and Clinical Performance (SSCP) can be found at: <https://ec.europa.eu/tools/eudamed> by searching for Manufacturer: Paragon 28 (upon Eudamed activation). Until Eudamed is activated, to request a copy of the SSCP, please contact Paragon 28®, Inc. by phone, (+1) (855) 786-2828.

Residual Risks identified within the SSCP are included as warnings and precautions.

IMPLANT MATERIALS

All Baby Gorilla®/Gorilla® Plates, Screws and Washers are made from Commercially Pure (CP Gr3, CP Gr4) Titanium (ASTM F67) and Titanium Alloy (ASTM F136). The maximum possible composition of each element is as follows:

Element	Maximum Composition, % (Mass/Mass)
Nitrogen	0.05
Carbon	0.10
Hydrogen	0.015
Iron	0.50
Oxygen	0.13
Aluminum	6.5
Vanadium	4.5
Titanium	99.185

INSTRUMENT MATERIALS

The instrumentation is made from medical grades of stainless steel, silicone, polymers, anodized aluminum, and Titanium Alloy.

CMR, EDP, PHTHALATES

The Baby Gorilla®/Gorilla® System Implants, guide wires, and drill bits do not contain substances ≥0.1% (w/w) which are classified as carcinogenic, mutagenic or reprotoxic (CMR), substances ≥0.1% (w/w) which have endocrine disrupting properties (EDP) and do not contain phthalates.

INTENDED USE/ INTENDED PURPOSE

The bone plates and bone screws of the Baby Gorilla®/Gorilla® Plating System are intended for use in primary and revision procedures including bone reconstruction/osteotomies, arthrodesis/joint fusion, and fracture repair/fracture fixation of the foot and ankle, appropriate for the size of the device.

INDICATIONS FOR USE

System(s)	Indications for Use	Intended Population(s)
Gorilla® Plating System Plates and Bone Screws in Diameters: 2.7 mm, 3.5 mm, and 4.2 mm	The Gorilla® plates and bone screws of the Gorilla® Plating System are indicated for use in the foot for: - Bone Reconstruction/Osteotomy - Arthrodesis/Joint Fusion - Fracture Repair/Fracture Fixation The Gorilla® plates and bone screws of the Gorilla® Plating System are indicated for use in the ankle for: - Bone Reconstruction/Osteotomy - Fracture Repair/Fracture Fixation	Adult and Pediatric Patients
Baby Gorilla® Plating System Plates and Bone Screws in Diameters: 2.0 mm and 2.5 mm	The Baby Gorilla® plates and bone screws of the Gorilla® Plating System are indicated for use in the foot for: - Bone Reconstruction/Osteotomy - Arthrodesis/Joint Fusion - Fracture Repair/Fracture Fixation	Adult and Pediatric Patients
Silverback™ Plating System Plates and Bone Screws in Diameters: 4.5 mm, 4.7 mm, 5.2 mm	The Silverback™ plates and bone screws of the Gorilla® plating System are indicated for use in the ankle for: - Arthrodesis/joint fusion	Adult Patients
R3LEASE™ Stabilization System Solid Screws and Washers in 3.9 mm Diameter	The R3LEASE™ Stabilization Screws of the Gorilla® Plating System are indicated for use in the ankle for: - Fracture repair/fracture fixation	Adult Patients
Non-locking Screws and Washers Diameters: 2.0 mm, 2.5 mm, 2.7 mm, 3.5 mm, 4.2 mm, 4.5 mm, 4.7 mm, 5.2 mm, and 5.5 mm	The Non-Locking Bone Screws and Washers of the Gorilla® Plating System are indicated for use in the foot and ankle for: - Bone Reconstruction/Osteotomy - Arthrodesis/Joint Fusion - Fracture Repair/Fracture Fixation	Adult and Pediatric Patients

CONTRAINDICATIONS

Use of the Baby Gorilla®/Gorilla® Plating System is contraindicated in cases of inflammation, cases of active or suspected sepsis / infection and osteomyelitis, or in patients with certain metabolic diseases and based on patient need and surgeon judgement.

All applications that are not defined by the indications are contraindicated. In addition, surgical success can be adversely affected by:

- Acute or chronic infections, local or systemic
- Vascular, muscular or neurological pathologies that compromise the concerned extremity
- All concomitant pathologies that could affect the function of the implant
- Osteopathy with reduced bone substance that could affect the function of the implant

- Any mental or neuromuscular disorder that could result in an unacceptable risk of failure at the time of fixation or complications in post-operative treatment
- Known or suspected sensitivity to metal
- Corpulence; an overweight or corpulent patient can strain the implant to such a degree that stabilization or implant failure can occur
- Whenever the use of the implant comes into conflict with the anatomical structures of physiological status
- Indications not included in the **INDICATIONS FOR USE**

Other medical or surgical pre-conditions that could compromise the potentially beneficial procedure, such as:

- The presence of tumors
- Congenital abnormalities
- Immunosuppressive pathologies
- Increased sedimentation rates that cannot be explained by other pathologies
- Increased leukocyte (WBC) count
- Pronounced left shift in the differential leukocyte count

WARNINGS AND PRECAUTIONS

- Re-operation to remove or replace implants may be required at any time due to medical reasons or device failure. If corrective action is not taken, complications may occur.
- Use of an undersized plate or screw in areas of high functional stresses may lead to implant fracture and failure.
- Dissimilar metals should not be used together.
- The implants and guide wires are intended for single use only. Any single use device which has been used may have come in contact with tissue, bone, blood or other bodily fluids. This is considered contaminated and must be discarded for the safety of the patient and other users.
- Instruments, guide wires and screws are to be treated as sharps.
- Patient risk as a result of the Baby Gorilla®/Gorilla® Plating System in the MR environment has been minimized.
- **Do not use other manufacturer's instruments or implants in conjunction with the Baby Gorilla®/Gorilla® Plating System.**
- **Do not implant the instruments.**
- In any surgical procedure, the potential for complications and adverse reactions exist. Adverse events include but are not limited to those described in this document. The residual risks, potential complications and adverse reactions with these implants include:
 - Loosening, deformation or fracture of the implant
 - Acute post-operative wound infections and late infections with possible sepsis
 - Migration, subluxation of the implant with resulting reduction in range of movement
 - Fractures resulting from unilateral joint loading
 - Thrombosis and embolism
 - Wound hematoma and delayed wound healing
 - Temporary or protracted neurological perturbation
 - Tissue reactions as the result of allergy or foreign body reaction to dislodged particles
 - Corrosion with localized tissue reaction and pain
 - Pain, a feeling of malaise or abnormal sensations due to the implant used
 - Delay or failure of bone healing
 - Loss of correction

All possible complications listed here are not typical of Paragon 28®, Inc. products but are in principle observed with any implant. Promptly inform Paragon 28®, Inc. as soon as complications occur in connection with the implants or surgical instruments used. In the event of premature failure of an implant in which a causal relationship with its geometry, surface quality or mechanical stability is suspected, please provide Paragon 28®, Inc. with the explant(s) in a cleaned, disinfected and sterile condition. Paragon 28®, Inc. cannot accept any other returns of used implants. The surgeon is held liable for complications associated with inadequate asepsis, inadequate preparation of the osseous implant bed in the case of implants, incorrect indication or surgical technique or incorrect patient information and consequent incorrect patient behavior.

MR SAFETY INFORMATION

 <p>MRI Safety Information A person with the Baby Gorilla®/Gorilla® Plating System Implant may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.</p>	
Device Name	Baby Gorilla®/Gorilla® Plating System Implant
Static Magnetic Field Strength (B0)	1.5 T or 3 T
Maximum Spatial Field Gradient	30 T/m (3000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Whole body transmit coil, Head RF transmit-receive coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2.0 W/kg
Scan Duration	2.0 W/kg or whole body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)
MR Image Artifact	32 mm

MAINTAINING DEVICE EFFECTIVENESS

- The surgeon should have specific training, experience, and thorough familiarity with the use of screw/plating systems.
- The surgeon must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.
- The Baby Gorilla®/Gorilla® plates and screws are not intended to endure excessive abnormal functional stresses.
- The Baby Gorilla®/Gorilla® plates and screws are intended for temporary fixation only until osteogenesis occurs.
- Failure to use dedicated, unique Baby Gorilla®/Gorilla® Plating System instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Carefully inspect the plates and screws prior to use, inspect the instruments before and after each procedure to assure they are in proper operational condition. Instruments which are faulty, damaged or suspect should not be used.
- Paragon28 non-sterile devices are precision medical devices and must be used and handled with care. Inspect the devices for damage prior to use, and at all stages of handling thereafter. If damage is detected, do not use the device prior to consulting the manufacturer for guidance.
- The surface coating of anodized aluminum devices may degrade due to exposure to highly alkaline cleaning processes and/or inappropriate handling and/or normal use. Such devices may require replacement if they no longer perform as designed.
- If an instrument is damaged, contact Paragon 28® immediately to arrange replacement and/or disposal.
- Devices with cutting functions or sharp points become dull with continuous use. This condition does not indicate a defective device and may indicate normal wear. Any sign of a dull or wearing devices may require replacement if they no longer perform as designed. Inspection prior to use should include verifying the cutting ability and sharpness of these points and edges.
- Paragon 28®, Inc. recommends the use of Paragon 28®, Inc. products in a sterile environment.

HANDLING AND STERILIZATION

Sterile Product

Baby Gorilla®/Gorilla® Plating System implants may be provided sterile. If sterile, product has undergone gamma irradiation. Do not re-sterilize. SINGLE USE ONLY. Any device labeled as Single Use only must never be reused or reprocessed. Any device opened from its packaging is considered used and may have come in contact with tissue, bone, blood or other bodily fluids. This is considered contaminated and must be discarded for the safety of the patient and other users. Do not use implants after expiration date. Packages for implants should be intact upon receipt.

Implants and instruments in sterile packaging should be inspected to ensure that the package has not been damaged or previously opened. If the inner package integrity has been compromised, potential for patient harm could occur, DO NOT USE THE IMPLANT or INSTRUMENT. Contact the manufacturer for further instructions. The implants should be opened using the aseptic technique described in P99-STR-1001. The implant or instrument should only be opened after the correct size has been determined. Once the seal of the product is broken, the product should not be re-sterilized.

Product should be stored in a clean, dry environment.

Non-Sterile Product

Products that are presented in a tray are provided non-sterile. All non-sterile implants and instruments should be cleaned using established hospital methods before sterilization and introduction into the sterile surgical field. Compliance is required with the manufacturer's user instructions and recommendations for chemical detergents. For manual reusable instrument reprocessing instructions, refer to the Baby Gorilla®/Gorilla® Plating System- **Instrument Reprocessing Instructions for Reusable Instruments** (Manual Method P51-CLN-0001) (Automated Method P99-CLN-1001).

Unless specifically labeled sterile, the implants and instruments are supplied NONSTERILE and MUST be sterilized prior to use. Recommended sterilization methods include steam autoclaving after removal of all protective packaging and labeling. Prior to sterilization, verify that all instruments are in their open and unlocked position within the instrument tray(s). The use of 2 layers of sterilization wrap, applied separately to each tray/caddy, is required. Extended dry time parameters are required due to the nature of the load and the materials used. The following minimum validated steam autoclave cycle and dry time is recommended:

Method	Cycle	Temperature	Exposure Time	Dry Time
Steam	Pre-Vacuum	273°F (134°C)	3 minutes	60 minutes within autoclave +45 minutes open door +45 minutes cooling rack

Baby Gorilla, Silverback, and Ankle Fracture Slimline Kits

Method	Cycle	Temperature	Exposure Time	Dry Time
Steam	Pre-Vacuum	273°F (134°C)	3 minutes	30 minutes within autoclave +30 minutes cooling rack

INTENDED USERS

Only surgeons who are fully experienced in the use of such implants and the required specialized surgical techniques should implant the Baby Gorilla®/Gorilla® Plating System. Refer to the Baby Gorilla® (P53-STG-1001) and Gorilla® Plating System Surgical Technique (P51-STG-1001) for complete instructions for use. For product information or to obtain a copy of the surgical technique manual, please contact Paragon 28®, Inc. by phone, (+1) (855) 786-2828

IMPLANT REMOVAL (IF NECESSARY)

- Locate implant with intra-operative imaging.
- Palpate plate and head of screw removing surrounding soft tissue to gain maximum exposure.
- Engage screw head with appropriate driver and rotate counterclockwise until screw is removed from plate. Repeat as many times as necessary.
- OPTION: If screw head is stripped, engage proximal shaft of screw under screw head with a medium sized Kern forceps and continue turning driver shaft and Kern forceps counterclockwise while exerting light pressure upwards with the Kern forceps.
- If screw is integrated into bone, core out with appropriately sized trephine drill.

TARGETED POPULATION

The Baby Gorilla®/Gorilla® Plating System is intended for members of the adult and pediatric general population who are under the care of a doctor who determines the need for bone reconstruction, osteotomy, arthrodesis, joint fusion, fracture repair and fracture fixation of the foot and ankle requiring fixation, appropriate for the size of the device. The population does not include patients that are pregnant except in emergency situations where bone fixation is required; nor does it include patients with conditions listed in the contraindications.

PRODUCT COMPLAINTS

The customer or health care provider should report any dissatisfaction with the product quality, labeling, performance or device malfunction to Paragon 28®, Inc. or its EC Representative immediately by phone or written correspondence. The user and/or patient should report any suspected serious incident related to the device by informing Paragon 28 or the local Paragon 28 distributor and the competent authority, ministry of health, or delegated agency in the country / region where the incident has occurred.

DISPOSAL OF IMPLANT OR INSTRUMENT

The implant and associated instrumentation will be used in a surgical suite, and thus, any disposal of the implant(s), instrument(s), and/or packaging will occur in the surgical suite by the user.

- The sterile packaged implant will be opened in the surgical suite with the packaging disposed of by the user. The outside packaging is not passed onto the operative field and should not be contaminated by human tissue or blood. The outside packaging is disposed of in a standard waste receptacle. If the inner packaging is in contact with human tissue or blood, the inner packaging must be disposed of in a standard biohazard container. If the inner packaging does not come in contact with human tissue or blood, the inner packaging can be disposed of in a standard waste receptacle.
- In the case of an implant removal, or if an implant must be disposed of during surgery, the disposal will occur in the surgical suite by the user. Because the implant has been in contact with human tissue or blood and may have a sharp physical hazard, the implant must be disposed of in a biohazard waste container intended for sharps.
- In the case of instrument disposal, the disposal will occur in the surgical suite by the user. Because the instrument has been in contact with human tissue or blood and may have a sharp physical hazard, the instrument must be disposed of in a biohazard container intended for sharps.

Please contact Paragon 28®, Inc. for product inquiries, cleaning instructions and surgical techniques. To report any adverse event, please contact Paragon 28®, Inc or its EC Representative.

	EC REP	UK Responsible Person	CH REP
 Paragon 28, Inc. 14445 Grasslands Dr., Englewood, CO 80112 USA (+1) (855) 786-2828	Emergo Europe Westervoortsedijk 60 6827 AT Arnhem The Netherlands	MedEnvoy UK Ltd 85 Great Portland St. London, W1W 7LT, UK	Zimmer GmbH Zählerweg 4 6300 Zug Switzerland

 2797	Class Ir, IIa, IIb
	Class I Self Certify

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