



## Paraderm® Dermal Matrix (Paraderm®)

DONATED HUMAN TISSUE

RESTRICTED TO USE BY OR ON THE ORDER OF A LICENSED HEALTHCARE PROFESSIONAL (physician, dentist, podiatrist, optometrist, nurse practitioner or physician assistant). NOT INTENDED FOR VETERINARY USE.

80-573 Rev. 01

### DESCRIPTION

Paraderm® is a freeze-dried decellularized dermal matrix aseptically recovered with consent from qualified donors. Paraderm® is processed using aseptic techniques and freeze-dried. The allograft is aseptically packaged in a tear pouch within a peel pouch and secured in an outer container to ensure allograft integrity. The inner pouch has been backfilled with nitrogen. Paraderm® is available in fenestrated and non-fenestrated options. The non-fenestrated product is available in the following sizes and thicknesses: 4x4cm and 4x8 cm (0.4mm-0.8mm; 1.0mm-2.0mm; and 2.2mm-3.5mm). The fenestrated product is available in the following sizes and thicknesses: 4x4cm and 4x8 cm (0.4mm-0.8mm; 0.9mm-1.4mm).

### INTENDED USE

Paraderm® is intended for use in integumentary augmentation.

### CONTRAINDICATIONS

Paraderm® is contraindicated in patients with known sensitivities or allergies to any agents listed below in the Warnings section.

### DONOR ELIGIBILITY

Paraderm® is recovered from a qualified donor and processed using aseptic techniques in accordance with federal, state, and/or international regulations and to the standards of the American Association of Tissue Banks. The donor is screened and tested for communicable disease risks and other exclusionary medical conditions. The results of the donor screening and testing have been reviewed by the Medical Director of VIVEX Biologics, Inc., and the donor has been deemed suitable for transplantation.

Communicable disease testing is performed by an FDA-registered laboratory certified to perform such testing on human specimens under the Clinical Laboratory Improvement Amendments of 1988 (42 U.S.C. 263a) and 42 CFR Part 493, or that has met equivalent requirements as determined by the Centers for Medicare and Medicaid Services in accordance with those provisions. Results from the following infectious disease tests were found to be nonreactive or negative:

#### Human Immunodeficiency Virus (HIV)

HIV-1/2 Plus O Antibodies (HIV-1/2 Plus O Ab)  
Nucleic Acid Test for HIV-1 RNA (HIV-1 NAT)

#### Hepatitis B Virus (HBV)

HBV Surface Antigen (HBsAg)  
HBV Core Antibody (IgG & IgM) (HbCAb)

Nucleic Acid Test for HBV DNA (HBV NAT)

#### Hepatitis C Virus (HCV)

HCV Antibody (HCVAb)  
Nucleic Acid Test for HCV RNA (HCV NAT)

#### Syphilis\*

Rapid Plasma Reagin (RPR) Screen  
T. pallidum IgG

\*A donor whose blood specimen is unsuitable for the non-treponemal screening assay, such as the RPR test, or with a reactive result from the non-treponemal screening assay, are cleared for transplantation use only when the result from the treponemal-specific (confirmatory) assay is nonreactive.

Screening tests for exposure to other viruses or parasites such as those listed below may have been completed. A negative/nonreactive result may not be required for these tests; however, all results are evaluated on a case-by-case basis by the Medical Director of VIVEX Biologics.

#### Cytomegalovirus\*\*

CMV Ab (IgG & IgM)

#### Epstein Barr Virus

EBV Ab (IgG & IgM)

#### Human T Cell Lymphotropic Virus I/II\*\*

HTLV-I/II (Antibody HTLV-I/II-Ab)

#### Toxoplasma gondii

Toxoplasma Ab (IgG & IgM)

#### Trypanosoma cruzi

T. cruzi Ab (IgG & IgM)

#### West Nile Virus (WNV)

Nucleic Acid Test for WNV RNA (WNV NAT)

#### Zika Virus

Zika Ab (IgM)

Nucleic Acid Test for Zika RNA (Zika NAT)

\*\*A donor with a reactive result for the CMV or HTLV-I/II Antibody test is suitable for use only when the result from a confirmatory assay is nonreactive.

### WARNINGS

The donor of Paraderm® is screened and tested for relevant communicable diseases and disease agents. Paraderm® is processed using aseptic techniques, microbiologically tested and may be exposed to Gentamicin, Vancomycin, peracetic acid, phosphate buffer solution or polyoxyethylene lauryl ether. Although the tissue is rinsed using sterile water or sterile saline during the manufacturing process, trace amounts may remain. Although all efforts have been made to ensure the safety of the allograft, there is no assurance that this allograft is free from all infectious diseases or microbial contamination. Paraderm® may transmit infectious agents.

**DO NOT FREEZE** the allograft by any method.

**FOR USE IN ONE PATIENT, ON A SINGLE OCCASION ONLY.**

**DO NOT STERILIZE** the allograft by any method. Exposure of the allograft and packaging to irradiation, steam, ethylene oxide, or other chemical sterilant may render the allograft unfit for use.

### ADVERSE EVENTS AND REACTIONS

Possible adverse events may include:

- Immunologic response (the possibility that a patient may develop alloantibodies should be considered for any patient who might be a future recipient of allograft tissue or cells)
- Transmission of disease of unknown etiology and transmission of infectious agents including but not limited to: HIV, syphilis, or microbial contaminants
- Infection of soft tissue and/or bone
- Fever

80-573 Rev. 01  
P01-IFU-0009

## STORAGE

Paraderm® must be stored at ambient temperature (2°C to 30°C). It is the responsibility of the tissue dispensing service, tissue distribution intermediary, and/or end-user clinician to maintain the allograft in appropriate storage conditions prior to further distribution and to use and track expiration dates accordingly. Appropriate inventory control should be maintained so that the allograft with the earlier expiration date is preferentially used and expiration is avoided.

## PRECAUTIONS

Paraderm® is processed and packaged using aseptic techniques and must be handled in an aseptic manner to prevent contamination.

Once the allograft container seal has been compromised, the allograft must be transplanted, if appropriate, or otherwise discarded.

Do not use the allograft if the pouch integrity has been compromised.

The outermost pouch is not sterile and should not be placed on an operative field.

The allograft must be reconstituted prior to use. The allograft should be rehydrated until soft and entirely pliable prior to use. Do not use, bend, or manipulate the allograft unless it is rehydrated.

## ALLOGRAFT PREPARATION

Step 1: Remove the pouch containing the allograft from the outer container.

Step 2: Inspect the pouch packaging for any holes, tears, or incomplete seals.


Step 3: Using aseptic technique, peel open the outer peel pouch from the chevron end and present the inner pouch to the operative field, when required.

Step 4: Locate the tear notch on the pouch and tear open. Place the allograft into a sterile basin.

Step 5: Aseptically pour sterile solution into the container until the allograft is completely immersed. Antibiotics of the end-user's preference may be added to the solution if desired. After the graft is completely rehydrated, it is ready for use.

## ALLOGRAFT ORIENTATION

A notch is located on the upper left-hand corner to differentiate between the basement membrane and the dermal sides of the allograft. Orient the allograft according to the figure below.



## RECIPIENT INFORMATION

Patient records must be maintained for the purpose of traceability. It is the responsibility of the end-user clinician to provide VIVEX Biologics with information pertaining to the traceability of the allograft used. For this purpose, the postage paid Tissue Utilization Report (TUR) card is provided with the allograft. Once the allograft is used, peel off the small product labels provided on the product packaging and affix on the TUR card and applicable patient records. Complete the TUR card and mail to VIVEX Biologics, scan, and e-mail to turs@VIVEX.com, or fax to (888) 630-4321.

## ADVERSE REACTION OUTCOME AND COMPLAINT REPORTING

Adverse reaction outcomes potentially attributable to the allograft must be promptly reported to VIVEX Biologics at (888) 684-7783. Any other complaints must be promptly reported to Paragon 28, Inc. at (855) 786-2828.

## RETURNED GOODS POLICY

Due to the delicate biological nature of a processed allograft, it cannot be returned for credit. If for any reason the allograft must be returned, a return authorization is required from Paragon 28, Inc. prior to shipping. It is the responsibility of the healthcare institution returning the allograft to adequately package and label it for return shipment.

The product manufacturer warrants that the allograft will conform to the specifications set forth herein provided that the allograft is handled, stored, and implanted by health care professionals according to the requirements set forth herein or as provided by it in writing. The product manufacturer makes no other warranties regarding the allograft; specifically disclaims any implied or statutory warranties, including any warranty against disease transmission; and makes no representations or warranties concerning the biological properties or biomechanical properties of the allograft.

### Distributed by:

Paragon 28, Inc.  
14445 Grasslands Dr.,  
Englewood, Colorado 80112 USA  
TEL: (855) 786-2828 | FAX: (888) 728-1220  
www.paragon28.com



### Manufactured by:

VIVEX Biologics, Inc.  
2430 NW 116<sup>th</sup> Street  
Miami, Florida 33167 USA

Paraderm™ is a trademark of Paragon 28, Inc.