

# BEAST 100™ Demineralized Bone Matrix



## DONATED HUMAN TISSUE

*RESTRICTED TO USE BY OR ON THE ORDER OF A LICENSED HEALTHCARE PROFESSIONAL (physician, dentist, podiatrist, optometrist, nurse practitioner or physician assistant).*

**80-147 02 05/2017  
P01-IFU-0006 RevB**

BEAST 100™ is a 100% human-derived cortical bone and gel allograft. The allograft is processed using aseptic techniques, exposed to antibiotics (containing Gentamicin and either Vancomycin or Bacitracin), and subjected to a demineralization process by exposure to hydrochloric acid and hydrogen peroxide solutions. Although the allograft was rinsed with sterile saline and (or) sterile water through processing, traces of antibiotics, hydrochloric acid, hydrogen peroxide and phosphate buffer solution may remain in the product. The components of this allograft have positive OsteoInductive (OI) test results. The allograft is packaged into a sterile syringe and sealed in packaging consisting of one tear pouch and one peel pouch. The allograft is secured in an outer container to ensure allograft integrity.

### INTENDED USE

BEAST 100™ is intended for use as a bone void filler.

### CONTRAINDICATIONS

BEAST 100™ is contraindicated in patients with known sensitivities to Gentamicin, Vancomycin, Bacitracin, hydrochloric acid or hydrogen peroxide.

### DONOR ELIGIBILITY

BEAST 100 is recovered from a qualified donor and processed using aseptic techniques in accordance with federal, state, and/or international regulations and to the standards of the American Association of Tissue Banks. Each donor is screened and tested for communicable disease risks and other exclusionary medical conditions. The results of the donor screening and testing have been reviewed by the Medical Director (or licensed physician designee) of UMTB Biomedical, Inc. and the donor has been deemed suitable for transplantation.

Communicable disease testing is performed by an FDA-registered laboratory certified to perform such testing on human specimens under the Clinical Laboratory Improvement Amendments of 1988 (42 U.S.C. 263a) and 42 CFR part 493, or that has met equivalent requirements as determined by the Centers for Medicare and

Medicaid Services in accordance with those provisions. Results from the following infectious disease tests have been found to be nonreactive or negative:

#### **Human Immunodeficiency Virus (HIV)**

HIV-1/2 Antibodies (HIV-1/2-Ab)  
Nucleic Acid Test for HIV-1 RNA (HIV-1 NAT)

#### **Hepatitis B Virus (HBV)**

HBV Surface Antigen (HBsAg)  
HBV Core Antibody (IgG & IgM) (HBcAb)  
Nucleic Acid Test for HBV DNA (if performed) (HBV NAT)

#### **Hepatitis C Virus (HCV)**

HCV Antibody (HCVAb)  
Nucleic Acid Test for HCV RNA (HCV NAT)

#### **Human T Cell Lymphotropic Virus I/II\* (if performed)**

HTLV-I/II (Antibody HTLV-I/II-Ab)

#### **Syphilis\*\***

Rapid Plasma Reagin (RPR) Screen  
T. Pallidum IgG

\*A donor with a reactive result for the HTLV-I/II Antibody test is suitable for use only when the result from a confirmatory assay is nonreactive.

\*\*A donor whose blood specimen is unsuitable for the non-treponemal screening assay, such as the RPR test, or with a reactive result from the non-treponemal screening assay, is suitable for use only when the result from the treponemal-specific (confirmatory) assay is nonreactive.

Screening tests for exposure to other viruses or parasites such as those listed below may have been completed. A negative/nonreactive result is not required for these tests, however, all results are evaluated on a case-by-case basis by the Medical Director (or licensed physician designee).

#### **Cytomegalovirus**

CMV Ab (IgG & IgM)

#### **Epstein Barr Virus**

EBV Ab (IgG & IgM)

#### **Toxoplasma gondii**

Toxoplasma Ab (IgG & IgM)

#### **Trypanosoma cruzi**

T. cruzi Ab (IgG & IgM)

### WARNINGS

The donor of BEAST 100™ is screened and tested for relevant communicable diseases and disease agents in compliance with the FDA regulations relating to human cells, tissues, and cellular and tissue-based products (21 CFR part 1271). BEAST 100™ is processed using aseptic techniques and microbiologically tested. Although all efforts have been made to ensure the safety of the allograft, there is no assurance that this allograft is free from all infectious diseases or microbial contamination.

**DO NOT FREEZE** the allograft by any method.

**FOR USE IN ONE PATIENT, ON A SINGLE OCCASION ONLY.**

**DO NOT STERILIZE** the allograft by any method. Exposure of the allograft and packaging to irradiation, steam, ethylene oxide or other chemical sterilants may render the allograft unfit for use.

**80-147 02 05/2017  
P01-IFU-0006 RevB**

## **PRECAUTIONS**

BEAST 100™ is processed and packaged using aseptic techniques and must be handled in an aseptic manner to prevent contamination.

## **ADVERSE EVENTS**

Allogeneic cells or tissues can induce an immunologic response in the recipient. The possibility that a patient may develop alloantibodies should be considered for any patient who might be a future recipient of allograft tissue or cells.

Possible adverse events may include: immunologic response, transmission of disease of unknown etiology and transmission of infectious agents including but not limited to: HIV, hepatitis, syphilis, or microbial contaminants.

## **STORAGE**

BEAST 100™ must be stored at ambient temperature (2°C - 30°C). It is the responsibility of the Tissue Dispensing Service, Tissue Distribution Intermediary, and/or End-User Clinician to maintain the allograft in appropriate storage conditions prior to further distribution or use and to track expiration dates accordingly. Appropriate inventory control should be maintained so that the allograft with the earlier expiration date is preferentially used and expiration is avoided.

## **ALLOGRAFT PREPARATION**

THE OUTERMOST POUCH IS NOT STERILE AND SHOULD NOT BE PLACED ON AN OPERATIVE FIELD.

ONCE THE TEAR POUCH SEAL HAS BEEN OPENED, the allograft must be used, or otherwise discarded.

ONCE OPENED, the allograft must be used within 24 hours.

DO NOT USE THE ALLOGRAFT if the pouch integrity has been compromised.

**Step 1:** Remove the pouch containing the allograft from the outer packaging.

**Step 2:** Utilizing aseptic technique, peel open the outer peel pouch from the chevron end and present the inner pouch to the operative field.

**Step 3:** Locate the tear notch on the inner pouch and tear open.

**Step 4:** Remove the cap from the syringe tip or remove the syringe end cap completely.

**Step 5:** Apply pressure to the plunger to extrude the allograft.

**Step 6:** Mold the allograft into the desired shape and press it into the targeted area.

## **RECIPIENT INFORMATION**

Patient records must be maintained for the purpose of traceability. It is the responsibility of the End-user or the Clinician to provide UMTB Biomedical, Inc. with information pertaining to the traceability of the allograft used. For this purpose, the postage paid Tissue Utilization Report (TUR) card is provided with the allograft. Once the allograft is used, peel off the small product labels provided on the product packaging and affix on the TUR card and applicable patient records. Complete the TUR card and mail to UMTB Biomedical, Inc., scan and e-mail to turs@umtb.com, or fax to (888) 630-4321.

## **ADVERSE OUTCOME AND COMPLAINT REPORTING**

Adverse outcomes potentially attributable to BEAST 100™ or other complaints should be promptly reported to UMTB Biomedical, Inc. at (888) 684-7783.

## **RETURNED GOODS POLICY**

Due to the delicate biological nature of a processed allograft, it cannot be returned for credit. If for any reason the allograft must be returned, a return authorization is required from Paragon 28® Customer Service prior to shipping. If the outer peel pouch has been opened or compromised, returns will not be accepted. It is the responsibility of the healthcare institution returning the allograft to adequately package and label it for return shipment.



### **Distributed by:**

Paragon 28, Inc.  
4B Inverness Ct. E, Suite 280  
Englewood, Colorado USA  
TEL: (855) 786-2828  
FAX: (888) 728-1220  
www.paragon28.com

### **Manufactured by:**

UMTB Biomedical, Inc.  
1951 N.W. 7<sup>th</sup> Avenue, Suite 200  
Miami, Florida USA 33136

**BEAST 100™ is a trademark of Paragon 28®, Inc.  
UMTB® is a registered trademark of Vivex Biomedical, Inc.**