










## Gorilla® Calc Fracture Plating System




## PRODUCT INFORMATION

	SINUS TARSI PLATE	SINUS TARSI SUPPORT PLATE	PERIMETER PLATE
			
Total Plates	6	6	8
Plate Thickness	1.1 mm	1.1 mm	1.1 mm
Hole Configurations	6-Hole	10-Hole	11-Hole (X-Small), 15-Hole
Size	Small, Medium, Large	Small, Medium, Large	X-Small, Small, Medium, Large
Left and Right Offerings	Yes	Yes	Yes


## SINUS TARSI PLATE AND SINUS SUPPORT PLATES

- ▶ Minimally invasive
- ▶ Low profile
- ▶ Fragment specific
- ▶ Plate-shaped to capture highest quality bone
  - 3 holes under the subtalar joint buttress the posterior facet
  - 3 anterior holes vertically in line with the calcaneocuboid joint designed to capture best bone and allows for stable fracture fixation
- ▶ Titanium (Ti6Al4V-ELI) for increased strength




### Sinus Tarsi Support Plate Features

- ▶ Extended posterior holes allow for additional fixation in the posterior tuberosity
- ▶ Posterior angle of the plate is perpendicular to the subtalar joint and intended to maintain calcaneal length during healing



### Perimeter Plate

- ▶ Low profile
- ▶ Lateral extensile approach allows for increased visibility of fractures and soft tissue bundles
- ▶ Up to 15 holes with off-axis locking capability allow surgeon flexibility in placement of screws to target around comminuted bone
- ▶ Material selection of CP Grade 4 Titanium and open plate design allow for a more malleable plate



## SYSTEM SPECIFIC INSTRUMENTATION

### Sinus Tarsi Incision Guide

- Slides over the threaded drill guides included in this system allowing the surgeon to locate the anterior and posterior hole clusters and expose each with minimally-sized incisions



### Sinus Tarsi Inserter

- Inserter is attached to the plate using the threaded drill guide enabling the surgeon to slide the plate under the soft tissues more easily
  - Increases leverage and grip to facilitate positioning of plate



### Dissection Tool


Blunt tipped instrument designed to allow for controlled dissection of the posterior tuberosity

- Ramped tip sized to help avoid excessive dissection
- Equipped with an AO quick connect



### Schanz Pins and T-Handle


- Threaded 4.0 and 6.0 mm pins designed to allow for increased purchase in poor bone
- Provided to allow the surgeon to manipulate larger fragments of bone
- AO quick connect compatible





## Gorilla® Calc Fracture Plating System

CFPS-01-RevD  
2023-10-09

Paragon 28®, Inc.   
14445 Grasslands Dr.  
Englewood, CO 80112 USA  
(855) 786-2828

Paragon 28® Medical Devices Trading Limited  
First Floor Block 7 Beckett Way  
Park West Business Park  
Dublin 12  
D12 X884  
Ireland  
+353 (0) 1588 0350

™Trademarks and ®Registered Marks of Paragon 28®, Inc.  
© Copyright 2023 Paragon 28®, Inc. All rights reserved.  
Patents: [www.paragon28.com/patents](http://www.paragon28.com/patents)

Exclusively foot & ankle <sup>28</sup>  
**Paragon**®

[www.Paragon28.com](http://www.Paragon28.com)

For the contraindications, potential complications and adverse reactions, warnings and precautions associated with this device, please refer to the device specific instructions for use at <http://www.paragon28.com/ifus>