PRECISION® MIS BUNION SYSTEM

Precision® MIS Bunion System






DESIGN RATIONALE

- ▶ Wide array of screw diameters and lengths to accommodate a variety of patient anatomies
- ▶ Precision™ Outrigger designed to accomplish triplanar correction of hallux valgus
- ▶ Large guide wire diameters for placement of screws with reduced chance of skiving
- ► Angel wing trajectory guide to trial placement of cannulated screws
- ▶ Robust offering of instruments for Precision™ guided or free-hand technique


IMPLANT OFFERINGS



	Ø3.0 mm Chamfer Screws	Ø3.5 mm Chamfer Screws	Ø4.0 mm Chamfer Screws
		<u> </u>	WARRIED
Screw Lengths:	Available lengths: 12 mm–48 mm in 2 mm increments	Available lengths: 12 mm-70 mm in 2 mm increments	Available lengths: 16 mm–70 mm in 2 mm increments



OUTRIGGER FEATURES





OUTRIGGER FEATURES

Distal Metatarsal Head Translation


➤ Designed to translate the distal metatarsal fragment up to 100%





Metatarsal Head Derotation

▶ Derotate the metatarsal head up to 48°



Distal Metatarsal Articular Angle Correction

▶ Introduce up to 25° of DMAA adjustment





Angel Wing Guide Wire Trajectory Trial


► Visualize screw entry point, trajectory, and position in metatarsal head





FEATURED INSTRUMENTATION

Ø3.5 mm and Ø4.0 mm sizes)



Distal Screw Drill Guide Tubes*

(Available in Ø3.0 mm, Ø3.5 mm and Ø4.0 mm sizes)

TECHNIQUE SUMMARY





Incision & Osteotomy

Make incision and perform the burr osteotomy





Outrigger Guide Attachment

Insert the Distal Hook into the metatarsal and place K-wires through the Outrigger into the navicular





Deformity Correction—IM

Translate the 1st metatarsal until sufficient IM correction is achieved





Deformity Correction—Derotation


Swivel the distal K-wire dorsally until sufficient derotation is achieved





Deformity Correction—DMAA

Rotate the DMAA Knob on the Rotation Fan clockwise until sufficient DMAA correction is achieved





Chamfer Screw Targeting

Attach the Angel Wing and Guide Tubes to visualize wire and screw trajectories

TECHNIQUE SUMMARY



Wire Insertion

Insert guide wires through Guide Tubes



Proximal Screw Insertion

Drill and insert the proximal screw through the Guide Tubes



Outrigger Removal

Withdraw the wires places through the Outrigger and remove





Use the Drill Guide/Soft Tissue Protector to drill and insert the distal screw

PRECISION® MIS BUNION SYSTEM

Precision® MIS Bunion System

P28-PMISB-01-RevA 2024-01-29

Paragon 28°, Inc. 4 14445 Grasslands Dr. Englewood, CO 80112 USA (855) 786-2828 Paragon 28° Medical Devices Trading Limited First Floor Block 7 Beckett Way Park West Business Park Dublin 12 D12 X884 Ireland +353 (0) 1588 0350

™Trademarks and *Registered Marks of Paragon 28*, Inc. © Copyright 2024 Paragon 28*, Inc. All rights reserved. Patents: www.paragon28.com/patents



www.Paragon28.com

For the contraindications, potential complications and adverse reactions, warnings and precautions associated with this device, please refer to the device specific instructions for use at http://www.paragon28.com/ifus