

WEIGHT BEARING CT PROTOCOL GUIDE—CURVEBEAM HIRISE

MAVEN® PSI System—APEX 3D™ Total Ankle Replacement



Disclaimer: MAVEN® Patient-Specific Guides and Surgical Planning Case Reports were designed to effectively expedite alignment and positioning of the APEX 3D™ Total Ankle Replacement System. The contoured guides are based on individual patient anatomic structures, are aligned to the mechanical axis and utilize a CT based coordinate system to aid in implant placement. Adhering to this CT scan protocol is critical.

Preparation

- ▶ Patient should be in shorts for scan.
- ▶ No socks or any jewelry—ex: no ankle chains, toe rings, magnets, socks with copper fiber, etc.
- ▶ All scans for the MAVEN® PSI System are weight-bearing. Patients must be able to stand still for approximately 5 minutes for the entire patient exam.
- ▶ Bilateral, LFOV (Large Field of View) scan of both knees, shins, and ankles is required.

Procedure Selection & Patient Positioning:

- ▶ On the Patient Tab, select the “**CT P28 MAVEN TAR Protocol**”
- ▶ On the Protocol Tab, select the “**CT P28 MAVEN TAR LFOV**”
- ▶ Bilateral, LFOV positioning of the feet and knees (Figure 1 and Figure 3)
- ▶ Align Ankle Maleoli ideally along the green line that horizontally cut the blue circle (Figure 2)

Patient Positioning Images




Figure 1: Green box represents scanned area of patient




Figure 2: Alignment of feet within the CurveBeam HiRise with ankle maleoli aligned perpendicular to the green line




Figure 3: Lateral view of patient with stabilizer

Scan Requirements & DICOM Output:

- ▶ 5cm proximal from the joint line of the knee through the bottom of foot – must include metatarsals
- ▶ Scans with patient movement will need to be repeated
- ▶ **DICOM Output:**
 - Send **axial RECON of SV** (Super Volume) from patient tab
 - Images should be **0.6mm slice thickness**
 - There should be no reformatting, no use of VOI (Volume of Interest) Tool. Original **RECON should be exported**

DICOM Export/Sending:

- ▶ Locate the RECON study with “SV” at the start of the series description within the patient list tab in CubeVue
- ▶ Right click on the file and select DICOM export, and the Study Exporter window will open:



- Do NOT edit identity. This will eliminate patient information needed by Paragon 28®
- Deselect the box for “Compressed Format”
- Check the box for “**Generate zip file output**”
- PREFERABLY, Upload through the secure, rapid electronic transfer system: www.apexankle.com
- Burned DICOM CD can be mailed to:

Attn: MAVEN PSI
Paragon 28®
14445 Grasslands Dr.
Englewood, CO 80112



Contact for Assistance:

Paragon 28® MAVEN® PSI Support Team

Phone: 720-253-0327

E-mail: MAVENorders@paragon28.com

P10-CTS-0003 Rev A [11/07/23]

CurveBeam AI Technical Support

Phone: 267-483-8081, Option 1 - (USA)

E-mail: techsupport@curvebeam.com