Phantom® Hindfoot TTC/TC Nail System

BEFORE USING PRODUCT, READ THE FOLLOWING IMPORTANT INFORMATION

A FULL SYMBOLS GLOSSARY CAN BE FOUND AT: www.paragon28.com/resources

Please check the website, www.paragon28.com/ifus, for the most current instructions for use document.

This booklet is designed to assist in using the Phantom® Hindfoot TTC/TC Nail System. It is not a reference for surgical techniques.

CAUTION

Federal Law (USA) restricts this device to sale and use by, or on the order of, a physician.

GENERAL DESCRIPTION

The Paragon 28® Phantom® Hindfoot TTC/TC Nail System is comprised of the Phantom® TTC Nail, and the Phantom® ActivCore™ Nail. These devices include intramedullary nails, screws, and end cap components. The Phantom® nails are offered in a variety of lengths to accommodate variations in patient anatomy. The Phantom® screws insert through the intramedullary nail to secure the construct. These are offered in varying lengths to accommodate the anatomical fixation required.

The system instruments include, but are not limited to, guide wires, drills/reamers, guides, jigs, gauges, and drivers. The instruments are used to place and/or remove the implants.

INTENDED/EXPECTED CLINICAL BENEFITS

The Intended clinical benefits of the Phantom® Hindfoot TTC/TC Nail System are:

- Successful Implantation of Phantom® Hindfoot TTC/TC Nail System devices
- Bone Healing
- Improvement in Pain
- Improvement in Function

Summary of Safety and Clinical Performance (SSCP) can be found at:

<https://ec.europa.eu/tools/eudamed>

Residual Risks identified within the SSCP are included as warnings and precautions.

IMPLANT MATERIALS

The implants of the Paragon 28® Phantom® Hindfoot TTC/TC Nail System are made from Titanium Alloy (ASTM F136). The maximum possible composition of each element is as follows:

Element	Maximum Composition, % (Mass/Mass)
Nitrogen	0.05
Carbon	0.08
Hydrogen	0.015
Iron	0.25
Oxygen	0.13
Aluminium	6.5
Vanadium	4.5
Titanium	90.475

INSTRUMENT MATERIALS

The instrumentation is manufactured from medical grades of titanium alloy, stainless steel, anodized aluminum, and polymer.

CMR, EDP, PHTHALATES

The Paragon 28® Phantom® Hindfoot TTC/TC Nail System does not contain substances ≥0.1% (w/w) which are classified as carcinogenic, mutagenic, or toxic to reproduction (CMR), substances ≥0.1% (w/w) which have endocrine disrupting properties (EDP), and does not contain phthalates.

INTENDED USE/INTENDED PURPOSE

The Phantom® Hindfoot TTC/TC Nail System is intended for tibiototalcanal and tibiototalcanal arthrodesis (fusion) and to provide stabilization of the hindfoot and ankle.

INDICATIONS FOR USE

The Phantom® Hindfoot TTC/TC Nail System is intended for primary and revision tibiototalcanal and tibiototalcanal arthrodesis (fusion) and to provide stabilization of the hindfoot and ankle.

Examples of specific indications include:

- Post-traumatic or degenerative arthritis
- Previously infected arthrosis
- Revision of failed ankle arthrodesis
- Revision of failed total ankle arthroplasty
- Talar deficiency conditions such as avascular necrosis of the talus (requiring tibiototalcanal arthrodesis)
- Rheumatoid arthritis
- Osteoarthritis
- Nonunions or pseudarthrosis of hindfoot and distal tibia
- Trauma (severe or malunited tibial pilon fracture)
- Charcot foot (neuroarthropathy)
- Severe end-stage degenerative arthritis
- Severe foot/ankle deformity
- Pantalar arthrodesis

CONTRAINDICATIONS

The Phantom® Hindfoot TTC/TC Nail System implants are not designed or sold for any use except as indicated. Use of the Phantom® Hindfoot TTC/TC Nail System is contraindicated in the following situations:

- Active, suspected, or latent infection in the affected area
- Patients who are physiologically or psychologically inadequate
- Patients previously sensitized to titanium
- Longitudinal splits, fractures, or deformities
- Insufficient quantity or quality of bone to permit stabilization, conditions that retard healing (not including pathological fractures) and conditions causing poor blood supply
- Corpulence; an overweight or corpulent patient can strain the implant to such a degree that stabilization or implant failure can occur
- Open epiphyseal plates
- Patients with an insufficient plantar fat pad
- Patients with an intact asymptomatic subtalar joint
- Patients with significant tibial malalignment (>10 degrees in either sagittal or coronal plane)
- In patients where there is a possibility for conservative treatment
- Indications not included in the **INDICATIONS FOR USE**

WARNINGS AND PRECAUTIONS

- Re-operation to remove or replace implants may be required at any time due to medical reasons or device failure. If corrective action is not taken, complications may occur.
- Use of an undersized implant in areas of high functional stresses may lead to implant fracture and failure.
- Implants, wires, or other appliances of dissimilar metals should not be used together in or near the implant site.
- The implants and guide wires are intended for single use only.
- Instruments, guide wires and screws are to be treated as sharps.
- **Do not use other manufacturer's instruments or implants in conjunction with the Phantom® Hindfoot TTC/TC Nail System**
- **Do not re-sterilize** The Phantom® Hindfoot TTC/TC Nail System implants
- In any surgical procedure, the potential for complications and adverse reactions exist. The risks and complications with these implants include:
 - Loosening, deformation, or fracture of the implant
 - Acute post-operative infections and late infections with possible sepsis
 - Migration, subluxation of the implant with resulting reduction in range of movement
 - Fractures resulting from unilateral joint loading
 - Thrombosis and embolism
 - Wound hematoma and delayed wound healing
 - Temporary and protracted functional neurological perturbation
 - Tissue reactions as a result of allergy or foreign body reaction to dislodged particles
 - Corrosion with localized reaction and pain
 - Pain, a feeling of malaise or abnormal sensations due to the implant used
 - Bone loss due to stress shielding
 - Loss of fixation in bone attributable to nonunion, osteoporosis and/or markedly unstable comminuted fractures
 - Nonunion or malunion with rotation or angulation resulting in limb shortening or loss of anatomic positioning
 - Irritation of soft tissues, including impingement syndrome
 - Adverse events include but are not limited to those described in this document

All possible complications listed here are not typical of Paragon 28®, Inc. products but are in principle observed with any implant. Promptly inform Paragon 28®, Inc. as soon as complications occur in connection with the implants or surgical instruments used. In the event of premature failure of an implant in which a causal relationship with its geometry, surface quality or mechanical stability is suspected, please provide Paragon 28®, Inc. with the explant(s) in a cleaned, disinfected, and sterile condition. Paragon 28®, Inc. cannot accept any other returns of used implants. The

surgeon is held liable for complications associated with inadequate asepsis, inadequate preparation of the osseous implant bed in the case of implants, incorrect indication or surgical technique or incorrect patient information and consequent incorrect patient behavior.

MR SAFETY INFORMATION

MRI Safety Information	
	A patient with the Paragon 28® Phantom® Hindfoot TTC/TC Nail System may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.
Name/Identification of device	Paragon 28® Phantom® Hindfoot TTC/TC Nail System
Nominal value(s) of Static Magnetic Field [T]	1.5 T or 3 T
Maximum Spatial Field Gradient [T/m and gauss/cm]	30 T/m (3000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Whole body transmit coil, Head RF transmit-receive coil
Maximum Whole Body SAR [W/kg]	2.0 W/kg
Limits on Scan Duration	All anatomical regions can be safely scanned under the following conditions: 1.0 W/kg whole body average SAR for 40 minutes of continuous RF (a sequence or back to back series/scan without breaks) with a 20 minute cooling period for an hour long scanning session. Scanning of the knees and all anatomy superior to the knees can be safely scanned under the following conditions: 2.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)
MR Image Artifact	The presence of this implant may produce an image artifact of 23 mm.
If information about a specific parameter is not included, there are no conditions associated with that parameter.	

MAINTAINING DEVICE EFFECTIVENESS

- The surgeon should have specific training, experience, and thorough familiarity with the use of fixation devices.
- The surgeon must exercise reasonable judgment when deciding which implant type to use for specific indications.
- The Phantom® Hindfoot TTC/TC Nail System implants are not intended to endure excessive abnormal functional stresses.
- The Phantom® Hindfoot TTC/TC Nail System implants are intended for temporary fixation only until osteogenesis occurs.
- Failure to use dedicated, unique Phantom® Hindfoot TTC/TC Nail System instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Carefully inspect the implants prior to use. Do not modify the implants. Inspect the instruments before and after each procedure to assure they are in proper operational condition. Instruments which are faulty, damaged, or suspect should not be used.
- The surface coating of anodized aluminum devices may degrade due to exposure to highly alkaline cleaning processes and/or inappropriate handling and/or normal use. Such devices may require replacement if they no longer perform as designed.
- If an instrument is damaged, contact Paragon 28® immediately to arrange replacement and/or disposal.
- Paragon 28® non-sterile devices are precision medical devices and must be used and handled with care. Inspect the devices for damage prior to use, and at all stages of handling thereafter. If damage is detected, do not use the device prior to consulting the manufacturer for guidance.
- Devices with cutting functions or sharp points become dull with continuous use. This condition does not indicate a defective device and may indicate normal wear. Any sign of dull or wearing devices may require replacement if they no longer perform as designed. Inspection prior to use should include verifying the cutting ability and sharpness of these points and edges.
- Paragon 28®, Inc. recommends the use of Paragon 28®, Inc. products in a sterile environment.

HANDLING AND STERILIZATION

All sterile products should be stored in a clean, dry environment.

Sterile Product:

Paragon 28® Phantom® Hindfoot TTC/TC Nail System Nails are provided sterile using gamma irradiation. DO NOT re-sterilize. SINGLE USE ONLY. The risk of re-use of device includes potential for patient to develop infection. DO NOT use sterile product after expiration date. Packages for implants should be intact upon receipt.

The sterile packaging should be inspected to ensure that the package has not been damaged or previously opened. If the inner package integrity has been compromised, DO NOT USE THE PRODUCT. Contact the manufacturer for further instructions. The packages should be opened using aseptic technique. Once the seal of the product is broken, the product should not be re-sterilized.

Non-Sterile Product:

Product that is presented in a tray is provided non-sterile. All non-sterile implants and instruments should be cleaned using established hospital methods before sterilization and introduction into the sterile surgical field. Compliance is required with the manufacturer’s user instructions and recommendations for chemical detergents. Refer to the Paragon 28®, Inc. *Instrument Reprocessing Instructions for Reusable Instruments* (Manual Method P99-CLN-0001) (Automated Method P99-CLN-1001). For product information or to obtain a copy of the surgical technique manual and/or the Summary of Safety and Clinical Performance, please contact Paragon 28®, Inc by calling (+1) (855) 786-2828.

Unless specifically labeled sterile, the implants and instruments are supplied NONSTERILE and MUST be sterilized prior to use. Recommended sterilization methods include steam autoclaving after removal of all protective packaging and labeling. Prior to sterilization, verify that all instruments are in their open and unlocked position within the instrument tray(s). The use of 2 layers of sterilization wrap is recommended. The following validated steam autoclave cycle is recommended:

Method	Cycle	Minimum Temperature	Minimum Exposure Time	Dry Time
Steam	Pre-Vacuum	270° F (132° C)	4 Minutes	30 Minutes

INTENDED USERS

Only surgeons who are fully experienced in the use of such implants and the required specialized surgical techniques should implant the Paragon 28® Phantom® Hindfoot TTC/TC Nail System. Refer to the product-specific Phantom® Hindfoot TTC/TC Nail System Surgical Technique Guide (P31-STG-1001) or Phantom® ActivCore™ Surgical Technique Guide (P31-STG-1002) for complete instructions for use. For product information or to obtain a copy of the product-specific surgical technique guide please contact Paragon 28®, Inc. by phone, +1 (855) 786-2828.

IMPLANT REMOVAL (IF NECESSARY)

- Locate the implant with intra-operative imaging
- Attach jig to the implant construct
- Remove all screws from the construct and pass from the operative field
- Use Instrumentation provided to back the nail out of the foot until completely removed.
- See the product-specific surgical technique for detailed implant removal instructions.
- Removal instructions are provided in the Phantom® Hindfoot TTC/TC Nail System Surgical Technique Guide (P31-STG-1001) or ActivCore™ Surgical Technique (P31-STG-1002).

TARGETED POPULATION

The Phantom® Hindfoot TTC/TC Nail System is intended for members of the adult general population requiring tibiototalcaneal or tibioalcaneal arthrodesis (fusion). The population does not include patients that are pregnant nor any others with conditions listed in the contraindications.

PRODUCT COMPLAINTS

The customer or health care provider should report any dissatisfaction with the product quality, labeling, or performance to Paragon 28®, Inc. immediately. Paragon 28®, Inc. should be notified immediately of any product malfunction by telephone or written correspondence. When filing a complaint, the name, part number and lot number of the part should be provided along with the name and address of the person filing the complaint.

DISPOSAL OF IMPLANT OR INSTRUMENT

The implant and associated instrumentation will be used in a surgical suite, and thus, any disposal of the implant(s), instrument(s), and/or packaging will occur in the surgical suite by the user.

- The sterile packaged implant will be opened in the surgical suite with the packaging disposed of by the user. The outside packaging is not passed onto the operative field and should not be contaminated by human tissue or blood. The outside packaging is disposed of in a standard waste receptacle. If the inner packaging is in contact with human tissue or blood, the inner packaging must be disposed of in a standard biowaste container. If the inner packaging does not come in contact with human tissue or blood, the inner packaging can be disposed of in a standard waste receptacle.
- In the case of an implant removal, or if an implant must be disposed of during surgery, the disposal will occur in the surgical suite by the user. Because the implant has been in contact

- with human tissue or blood and may have a sharp physical hazard, the implant must be disposed of in a biohazard waste container intended for sharps.
- In the case of instrument disposal, the disposal will occur in the surgical suite by the user. Because the instrument has been in contact with human tissue or blood and may have a sharp physical hazard, the instrument must be disposed of in a biohazard container intended for sharps.

Please contact Paragon 28®, Inc. for product inquiries, cleaning instructions and surgical techniques. To report any adverse event, please contact Paragon 28®, Inc. or its EC Representative.

		UK Responsible Person	
Paragon 28, Inc. 14445 Grasslands Dr., Englewood, CO 80112 USA (+1) (855) 786 -2828	Emergo Europe Westervoortsedijk 60 6827 AT Arnhem The Netherlands	MedEnvoy UK Ltd 85 Great Portland St. London, W1W 7LT, UK	MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland

 2797	Class Ir, IIa, IIb
	Class I Self Certify

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- Paragon 28®, Inc. anbefaler at Paragon 28®, Inc.-produkter brukes i et sterilt miljø.