Grappler™ Interference Screw System

BEFORE USING PRODUCT, READ THE FOLLOWING IMPORTANT INFORMATION

A FULL SYMBOLS GLOSSARY CAN BE FOUND AT:
www.paragon28.com/resources

Please check the website, www.paragon28.com/ifu, for the most current instructions for use document.

This booklet is designed to assist in using the Grappler™ Interference Screw System. It is not a reference for surgical techniques.

CAUTION

Federal Law (USA) restricts this device to sale and use by, or on the order of, a physician.

GENERAL DESCRIPTION

The Grappler™ Interference Screw System is comprised of PEEK screws that are intended for soft tissue reattachment. The Grappler™ Interference Screw System implant is offered in various sizes to accommodate patient anatomy. The Grappler™ Interference Screw System instruments are used to facilitate the implantation of the Grappler™ Interference Screw.

INTENDED/EXPECTED CLINICAL BENEFITS

Intended/Expected Clinical Benefits of the Grappler™ Interference Screw System include:

- Deformity Correction
- Successful Implantation of Interference Screws
- Bone healing and soft tissue repair
- Acceptable Foot and Ankle Ability Post-Operative
- Acceptable Pain Score Post-Operative
- Patient Satisfaction

Summary of Safety and Clinical Performance (SSCP) can be found at:
<https://ec.europa.eu/tools/eudamed>.

Residual Risks identified within the SSCP are included as warnings and precautions.

IMPLANT MATERIALS

The implants of the Grappler™ Interference Screw System are manufactured from polyetheretherketone (PEEK) per ASTM F2026. The instrumentation is manufactured from medical grades of stainless steel, nitinol, anodized aluminum, polymers and silicone.

CMR, EDP, PHTHALATES

The Grappler™ Interference Screw System Implants, drill bits and k-wires do not contain substances ≥0.1% (w/w) which are classified as carcinogenic, mutagenic or toxic to reproduction (CMR), substances ≥0.1% (w/w) which have endocrine disrupting properties (EDP) and do not contain phthalates.

INTENDED USE / INTENDED PURPOSE

The Grappler™ Interference Screw System is intended for soft tissue to bone repair including ligament and tendon reconstruction in foot and ankle surgery. The Grappler™ Interference Screw System instruments are used to facilitate the implantation of the Grappler™ Interference Screw.

INDICATIONS FOR USE

The Grappler™ Interference Screw System is indicated for:

Lateral Ankle Instability

- Lateral Ankle Stabilization Procedure
- Brostrom Procedure
- Brostrom Gould Procedure
- Syndesmosis Repair Procedure

Achilles Tendon Repair/Reconstruction

- Flexor Hallucis Longus Transfer Procedure

Flatfoot Deformity ("Adult Acquired Flatfoot Deformity" or "Progressive Collapsing Foot Deformity") with one or more of the following clinical findings – deformity, synovitis, tendinosis, posterior tibial tendon dysfunction, deltoid insufficiency, and valgus deformity

- Deltoid Ligament Reconstruction Procedure
- Flexor Digitorum Longus Transfer Procedure
- Spring Ligament Repair Procedure
- Flexor Hallucis Longus Transfer Procedure

Neurogenic Dropfoot

- Posterior Tibial Tendon Transfer Procedure

Pes Cavus Deformity

- Lateral Ankle Stabilization Procedure
- Tarsal Coalition**
- Flexor Digitorum Longus Transfer Procedure

CONTRAINDICATIONS

The Grappler™ Interference Screw System implants are not designed or sold for any use except as indicated. Use of the Grappler™ Interference Screw System is contraindicated in the following situations:

- Active, suspected or latent infection in the affected area
- In patients who are physiologically or psychologically inadequate
- Insufficient quantity or quality of bone to permit stabilization, conditions that retard healing (not including pathological fractures) and conditions causing poor blood supply
- In patients where there is a possibility for conservative treatment
- Placement across open growth plates in skeletally immature patients
- In patients who are pregnant
- In locations where the tendon is irreparable
- Corpulence: an overweight or corpulent patient can strain the implant to such a degree that stabilization or implant failure can occur
- Indications not included in the **INDICATIONS FOR USE**

Other medical or surgical pre-conditions that could compromise the potentially beneficial procedure, such as:

- The presence of tumors
- Congenital abnormalities
- Immunosuppressive pathologies
- Increased sedimentation rates that cannot be explained by other pathologies
- Increased leukocyte (WBC) count
- Pronounced left shift in the differential leukocyte count

WARNINGS AND PRECAUTIONS

- Re-operation to remove or replace implants may be required at any time due to medical reasons or device failure. If corrective action is not taken, complications may occur.
- Use of an undersized implant in areas of high functional stresses may lead to implant fracture and failure.
- The implants and k-wires (single use instruments) are intended for single use only. Any single use device which has been used may have come in contact with tissue, bone, blood or other bodily fluids. This is considered contaminated and must be discarded for safety of the patient and other users.
- Instruments and implants are to be treated as sharps.
- **Do not use other manufacturer's instruments or implants in conjunction with the Grappler™ Interference Screw System.**
- **Do not re-sterilize the Grappler™ Interference Screw System Implants.**
- In any surgical procedure, the potential for complications and adverse reactions exists. The residual risks, potential complications and adverse reactions with these implants include:
 - Loosening, deformation or fracture of the implant
 - Acute post-operative infections and late infections with possible sepsis
 - Migration, subluxation of the implant
 - Wound hematoma and delayed wound healing
 - Temporary and protracted functional neurological perturbation
 - Tissue reactions as a result of allergy or foreign body reaction to dislodged particles or implant material
 - Pain, a feeling of malaise or abnormal sensations due to the implant used
 - Bone resorption or over-production
 - Migration of particle wear debris possibly resulting in a bodily response
 - Adverse events include but are not limited to those described in this document

All possible complications listed here are not typical of Paragon 28®, Inc. products but are in principle observed with any implant. Promptly inform Paragon 28®, Inc. as soon as complications occur in connection with the implants or surgical instruments used. In the event of premature failure of an implant in which a causal relationship with its geometry, surface quality or mechanical stability is suspected, please provide Paragon 28®, Inc. with the explant(s) in a cleaned, disinfected and sterile condition. Paragon 28®, Inc. cannot accept any other returns of used implants. The surgeon is held liable for complications associated with inadequate asepsis, inadequate preparation of the osseous implant bed in the case of implants, incorrect indication or surgical technique or incorrect patient information and consequent incorrect patient behavior.

MR SAFETY INFORMATION

The Grappler™ Interference Screw System has been evaluated for safety and compatibility in the magnetic resonance (MR) environment. It has not been tested for heating, migration, or image artifact in the MR environment. Since the Grappler devices are composed of nonconductive, nonmetallic, and nonmagnetic materials, there are no known hazards resulting from exposure of these devices to any MR environment. The Grappler™ Interference Screw System implants are MR Safe.

MAINTAINING DEVICE EFFECTIVENESS

- The surgeon should have specific training, experience, and thorough familiarity with the use of fixation devices.
- The surgeon must exercise reasonable judgment when deciding which implant type to use for specific indications.
- The Grappler™ Interference Screw System implants are not intended to endure excessive abnormal functional stresses.
- Failure to use dedicated, unique Grappler™ Interference Screw System instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Carefully inspect the implants prior to use, inspect the instruments before and after each procedure to assure they are in proper operational condition. Instruments which are faulty, damaged or suspect should not be used.
- Paragon28 non-sterile devices are precision medical devices and must be used and handled with care. Inspect the devices for damage prior to use, and at all stages of handling thereafter. If damage is detected, do not use the device prior to consulting the manufacturer for guidance.
- Devices with cutting functions or sharp points become dull with continuous use. This condition does not indicate a defective device and may indicate normal wear. Any sign of a dull or wearing devices may require replacement if they no longer perform as designed. Inspection prior to use should include verifying the cutting ability and sharpness of these points and edges.
- Disposal of any implants and instruments should be considered sharps and disposed of according to a healthcare facilities policy.
- Paragon 28®, Inc. recommends the use of Paragon 28®, Inc. products in a sterile environment.

HANDLING AND STERILIZATION

Sterile Product

Grappler™ Interference Screw System implants are provided sterile using gamma irradiation. Instruments may be provided sterile. Do not re-sterilize. SINGLE USE ONLY. Any device labeled as Single Use only must never be reused or reprocessed. Any device opened from its packaging is considered used and may have come in contact with tissue, bone, blood or other bodily fluids. This is considered contaminated and must be discarded for the safety of the patient and other users. Do not use sterile packed implants or instruments after expiration date. Packages for implants and instruments should be intact upon receipt.

Implants and instruments in sterile packaging should be inspected to ensure that the package has not been damaged or previously opened. If the package integrity has been compromised, potential for patient harm could occur, DO NOT USE THE IMPLANT or INSTRUMENT. Contact the manufacturer for further instructions. The implants and instruments should be opened and transferred to the sterile surgical field using the aseptic technique described in document P99-STR-1001. The implant or instrument should only be opened after the correct size has been determined. Once the seal of the product is broken, the product should not be re-sterilized.

All implants and instruments should be stored in a clean, dry environment.

Non-Sterile Product

Instruments that are presented in a case are provided non-sterile. All non-sterile instruments should be cleaned using established hospital methods before sterilization and introduction into the sterile surgical field. Compliance is required with the manufacturer's user instructions and recommendations for chemical detergents. Refer to the Paragon 28®, Inc. **Instrument Reprocessing Instructions for Reusable Instruments** (Manual Method P99-CLN-0001) (Automated Method P99-CLN-1001). This is also available by calling +1 (855) 786-2828.

Unless specifically labeled sterile, the instruments are supplied NONSTERILE and MUST be sterilized prior to use. Recommended sterilization methods include steam autoclaving after removal of all protective packaging and labeling. Prior to sterilization, verify that all instruments are in their open and unlocked position within the instrument tray(s). The use of 2 layers of sterilization wrap is recommended. Using one of the following validated steam autoclave cycles is recommended:

Method	Cycle	Minimum Temperature	Minimum Exposure Time	Dry Time
Steam	Pre-vacuum	270°F (132°C)	4 Minutes	30 Minutes
Steam	Pre-vacuum	273°F (134°C)	3 Minutes	30 Minutes

INTENDED USERS

Only surgeons who are fully experienced in the use of such implants and the required specialized surgical techniques should implant the Grappler™ Interference Screw System. Refer to the Grappler™ Interference Screw System Surgical Technique P41-STG-1001 for complete instructions for use. For product information or to obtain a copy of the surgical technique manual please contact Paragon 28®, Inc. by phone, +1 (855) 786-2828.

IMPLANT REMOVAL (IF NECESSARY)

- Instrumentation can be provided for implant removal.
- Removal instructions are provided in the Grappler™ Interference Screw System Surgical Technique (P41-STG-1001).

TARGETED POPULATION

The Grappler™ Interference Screw System is intended for members of the adult general population suffering from soft tissue injuries or deformities in the foot/ankle. The population does not include patients that are pregnant nor any others with conditions listed in the contraindications

PRODUCT COMPLAINTS

The customer or health care provider should report any dissatisfaction with the product quality, labeling, or performance to Paragon 28®, Inc. immediately. Paragon 28®, Inc. or its EC Representative should be notified immediately of any product malfunction by telephone or written correspondence. When filing a complaint, the name, part number and lot number of the part should be provided along with the name and address of the person filing the complaint.

DISPOSAL OF IMPLANT OR INSTRUMENT

The implant and associated instrumentation will be used in a surgical suite, and thus, any disposal of the implant(s), instrument(s), and/or packaging will occur in the surgical suite by the user.

- The sterile packaged implant will be opened in the surgical suite with the packaging disposed of by the user. The outside packaging is not passed onto the operative field and should not be contaminated by human tissue or blood. The outside packaging is disposed of in a standard waste receptacle. If the inner packaging is in contact with human tissue or blood, the inner packaging must be disposed of in a standard biowaste container. If the inner packaging does not come in contact with human tissue or blood, the inner packaging can be disposed of in a standard waste receptacle.
- In the case of an implant removal, or if an implant must be disposed of during surgery, the disposal will occur in the surgical suite by the user. Because the implant has been in contact with human tissue or blood and may have a sharp physical hazard, the implant must be disposed of in a biohazard waste container intended for sharps.
- In the case of instrument disposal, the disposal will occur in the surgical suite by the user. Because the instrument has been in contact with human tissue or blood and may have a sharp physical hazard, the instrument must be disposed of in a biohazard container intended for sharps.

Please contact Paragon 28®, Inc. for product inquiries, cleaning instructions and surgical techniques. To report any adverse event, please contact Paragon 28®, Inc. or its EC Representative.

	EC REP	UK Responsible Person	CH REP
Paragon 28, Inc. 14445 Grasslands Dr., Englewood, CO 80112 USA (+1) (855) 786 -2828	Emergo Europe Westervoortsedijk 60 6827 AT Arnhem The Netherlands	MedEnvoy UK Ltd 85 Great Portland St. London, W1W 7LT, UK	MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland

 2797	Class Ir, IIa, IIb
	Class I Self Certify

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- Paragon 28®, Inc. anbefaler at Paragon 28®, Inc.-produkter brukes i et sterilt miljø.