

Paragon²⁸[®]

SURGICAL TECHNIQUE GUIDE

Gorilla[®] Pilon Fusion Plating System

GORILLA[®]
ANKLE FRACTURE 360[™]



PRODUCT DESCRIPTION:

Paragon 28® designed the Gorilla® Pilon Fusion Plating System to provide surgeons versatility in fixation selection for pilon fusion procedures.. The system has 6 plate options for anterior and medial fixation of the tibia and talus for a pilon fusion. All circular plate holes accept Ø2.7 mm, Ø3.5 mm and Ø4.2 mm locking and non-locking screws.

A Ø3.5 mm Compact Screw is available in locking and non-locking screws in 10-40 mm lengths to address dense bone in the proximal tibia. The Compact Screw was designed with a smaller thread height to help reduce insertion torque in dense bone. Additionally, single lead bone threads result in a decreased pitch differential between the locking screw head and bone threads to reduce the amount of insertion torque required to lock the screw into the plate in areas of dense bone.

Instrumentation is included in the Gorilla® Ankle Fracture Plating System that facilitates reduction and fixation of ankle fractures.

CONTENTS:

| | |
|------------------|--|
| SECTION 1 | GORILLA® PILON FUSION PLATING SYSTEM OVERVIEW |
| | IMPLANT OFFERING.....3-4 |
| | FEATURED INSTRUMENTATION.....5-6 |
| SECTION 2 | SURGICAL TECHNIQUE |
| | ANTEROLATERAL PLATE.....7-15 |
| | MEDIAL PLATE.....16-19 |
| SECTION 3 | CADDY AND SAFETY INFORMATION |
| | CADDY CONTENTS20-25 |
| | INDICATIONS, CONTRAINDICATIONS, WARNINGS26-28 |


ACKNOWLEDGMENT:

Paragon28® would like to thank Clayton Bettin, M.D., John Kwon, M.D., Eric Moghadamian, M.D., and Aaron Perdue, M.D. for their contribution to the development of the surgical technique guide.

IMPLANT OFFERING:

ANTEROLATERAL PILON FUSION PLATE


- ▶ Offered in 3 sizes –Short, Medium, Long
- ▶ Left/Right Specific
- ▶ Designed to contour to anterolateral face of distal tibia
- ▶ Features optional compression slot and allows for placement of a Ø5.5 Monster® screw through the plate and across the tibiotalar joint
- ▶ Compatible with an optional Backpack Guide and Drill Guides that allows nominal trajectory screw placement in the distal cluster



IMPLANT OFFERING:

MEDIAL PILON FUSION PLATE- Supplemental (Use in conjunction with Anterolateral Plate for Pilon Fusion)

- ▶ Offered in 3 sizes –Short, Medium, Long
- ▶ Universal for left or right
- ▶ Designed to contour to medial malleolus and medial tibia
- ▶ Compatible with provided plate bending instrumentation for patient specific side contouring along the length of the medial tibia
- ▶ Distal portion of the plates allows for optional screw placement from the medial malleolus into the talus



FEATURED INSTRUMENTATION:



Proximal Trocar



Proximal Drill Guide



Assembled Proximal Trocar and Drill Guide



Proximal Screw Driver



Depth Gauge



Reciprocating Rasp (Compatible with Stryker Power System)



King Cobra™ Cartilage Removal Tool



Monster K-wire Guide



K-wire and Drill Guide Assembly



Monster Drill Guide



Backpack Drill Guide



Optional* Backpack Guide



Thread Attachment Rod



Plate Handle



Bone Fenestration Perforator




























































Curved Bone Fenestration Chisel



NOTE: Additional wires, olive wires, fracture reduction instrumentation, joint preparation instrumentation, plate bending instrumentation and drills/drill guides are offered in the AF360 instrumentation case.

SCREW OFFERING AND INSTRUMENTATION MATRIX

| | Ø2.7 mm R3CON Screws | Ø3.5 mm R3CON Screws | Ø4.2 mm R3CON Screws | Ø3.5 mm Compact Screws |
|--|--|---|--|--|
| Locking: |  |  |  |  |
| Non-locking: |  |  |  |  |
| Screw Lengths: | 8 mm - 20 mm in 1 mm increments 20 - 54 mm in 2 mm increments 56 mm - 70 mm in 2 mm increments (Non-locking ONLY) 70 mm - 75 mm in 5 mm increments (Non-locking ONLY) | 8 mm - 60 mm in 2 mm increments 60 mm - 70 mm in 5 mm increments 75 mm - 100 mm in 5 mm increments (Non-locking ONLY) | 8 mm - 60 mm in 2 mm increments 60 mm - 70 mm in 5 mm increments | 10 mm - 40 mm in 2 mm increments |
| Drill Size: |  Ø2.0 mm |  Ø2.4 mm |  Ø2.8 mm |  Ø2.8 mm |
| Driver Size: |  HX-10 |  HX-10 |  HX-10 |  HX-10 |
| Locking Drill Guide Size: |  Ø2.7 mm |  Ø3.5 mm |  Ø3.5 mm C / Ø4.2 mm |  Ø3.5 mm C / Ø4.2 mm |
| Centering Drill Guide Size: |  Ø2.7 mm |  Ø3.5 mm |  Ø4.2 mm |  Ø3.5 mm |
| Compression Slot Drill Guide Size: |  Ø2.7 mm |  Ø3.5 mm |  Ø3.5 mm C / Ø4.2 mm |  Ø3.5 mm C / Ø4.2 mm |
| Cone/Straight Easy Guide Size: |  Ø2.7 mm |  Ø3.5 mm |  Ø3.5 mm C / Ø4.2 mm |  Ø3.5 mm C / Ø4.2 mm |
| Tap Size: |  Ø2.7 mm |  Ø3.5 mm |  Ø4.2 mm |  Ø3.5 mm C |
| OverDrill Size: |  Ø2.7 mm |  Ø3.5 mm |  Ø4.2 mm |  Ø3.5 mm C |
| Double Ended Drill / Over Drill Guides: |  Ø2.0 mm |  Ø2.4 mm |  Ø2.8 mm |  Ø2.8 mm |
| Drill Sleeve (for use with Double Ended Guide): |  Ø2.0 mm Drill / Ø2.7 mm Over Drill |  Ø2.4 mm Drill / Ø3.5 mm Over Drill |  Ø2.8 mm Drill / Ø4.2 mm Over Drill |  Ø2.8 mm Drill / Ø4.2 mm Over Drill |

| | Ø4.0 mm Mini-Monster Cannulated Screws |
|---------------------------------|--|
| Headed, Long Thread: |  |
| Screw Lengths: | 24 mm - 50 mm in 2 mm increments 50 mm - 60 mm in 5 mm increments |
| Drill Size: |  Ø2.6 mm |
| Driver Size: |  HX-10 |
| Drill Guide Size: |  Ø4.0 mm |
| Headed Countersink Size: |  Ø4.0 mm |
| Tap Size: |  Ø4.0 mm |
| Over Drill Size: |  Ø4.0 mm |
| Over Drill Guide Size: |  Ø4.0 mm |
| K-wire Size: |  Ø1.2 mm x 15 cm |

ANTEROLATERAL PILON FUSION PLATE

INCISION



Pre-operative planning should be performed prior to the surgery, including evaluation of soft tissue condition, and review of radiographs and/or advanced imaging to determine approach and internal fixation needs. An anterior (shown) or anterolateral approach is recommended with supine patient positioning, although an anteromedial approach can be utilized depending on fracture morphology. Instrumentation is provided to assist with cartilage removal from the talus and tibia and a fresh scalpel can be used to remove cartilage from small loose articular fracture fragments.

For pilon fractures, it is recommended to restore large metadiaphyseal pieces together with provisional screw and/or wire fixation per surgeon preference.



Reciprocating Rasp



King Cobra™ Cartilage Removal Tool



NOTE: Autograft, allograft, or other biologics (e.g. V92) per surgeon preference are recommended to be used before plate placement.

ANTEROLATERAL PILON FUSION PLATE




Plate Handle



Thread Attachment

Use the provided instrumentation (Lagenbeck Elevator) to elevate the soft tissues along the length of the anterolateral tibial crest for the intended plate to be used. Retrieve the appropriate length Anterolateral Plate. Attach the handle onto the plate by threading the thread attachment rod through the handle and into the plate. Depending on the approach type, plate size, and/or patient anatomy the handle can be used on the medial or lateral side of the anatomy. When placing the plate, ensure the plate is contacting both the proximal lateral and distal anterior surface of the tibia and is not in floating in the soft tissue. Confirm the plate is positioned appropriately on the talus. Confirm plate placement using fluoroscopy. If provisional fixation of the plate to bone is preferred, olive wires are available to place within the screw holes.



Threaded Hole



ANTEROLATERAL PILON FUSION PLATE



After temporarily securing the plate position with olive wires and confirming position, unthread and remove the handle from the plate. Continue with drilling and placing appropriately sized screws according to surgeon preference in the distal portion of the plate.



ANTEROLATERAL PILON FUSION PLATE

OPTIONAL BACKPACK GUIDE USE

If using the optional Backpack Guide, it can be attached to the plate using the same threaded attachment rod as used with the handle. Place the Backpack drill guide into the Backpack Guide to guide drills at a nominal trajectory with respect to the plate using the backpack drill guides.




Backpack Guide




Thread Attachment



Backpack Drill Guide



NOTE: When measuring screw length through the backpack guide, remove the drill guides and use provided pilon fusion depth gauge.



Monster Screw Hole



ANTEROLATERAL PILON FUSION PLATE



Monster K-Wire Guide



Monster Drill Guide



Assembled Monster K-Wire and Drill Guide




Ø1.60 mm x 15 cm Smooth K-Wire


For placement of a Monster® screw through the plate and across the tibiotalar joint, assemble the Monster K-wire and Drill Guide. Place the assembled guide into the Monster screw hole in the Backpack Guide on the plate. Place a Ø1.6 K-wire through the guide and advance the wire to the desired termination point based on surgeon preference. If the K-wire trajectory is not appropriate for patient anatomy, remove the wire and guides and continue with freehand placement of the wire and screw using Monster® instrumentation.



ANTEROLATERAL PILON FUSION PLATE




Remove the K-wire guide and drill over the K-wire with the Ø3.5 mm cannulated drill or remove the K-wire and K-wire Guide and drill with a solid 3.5 mm drill to the intended termination point.






Remove the drill and drill guide and place the appropriate length Monster screw (partial or fully threaded per surgeon preference) through the plate and in to the talus. Confirm screw position and trajectory using fluoroscopy

ANTEROLATERAL PILON FUSION PLATE





The drill guides should press fit into the backpack guide. Drill through the drill guides with the appropriately sized drills for the intended screw diameter to be used. Measure screw length based on the markings on the drill or measure using the provided depth gauge.



Remove the drill guide and place the selected screw into the drilled hole. Repeat the steps of drilling and screw placement for remaining plate holes based on surgeon preference.


ANTEROLATERAL PILON FUSION PLATE

For placement of the proximal lateral screws on the plate, use fluoroscopy to locate the plate screw holes and make a stab or open incision at the intended locations. A trocar and drill guide are available to help assist with minimizing soft tissue disruption for screw placement. After using the trocar to reach bone, and thread the drill guide into the plate.




ANTEROLATERAL PILON FUSION PLATE

If using the trocar with the drill guide, remove the trocar and drill through the drill guide using the calibrated drill for the appropriately sized screw to be used.




NOTE: When measuring through the drill guide with an AF360 drill, subtract 40 mm from the measured length.

Remove the drill guide and continue with placement of remaining proximal lateral screws according to surgeon preference.



MEDIAL PILON FUSION PLATE



A Medial Plate is offered to provide additional fixation and strength to the plating construct. **It is recommended to use the medial pilon plate along with the anterolateral plate and not as the only plate used.**

Make a longitudinal incision over the central aspect of the medial malleolus, appropriately sized for the fracture and plate length. Continue soft tissue dissection until the fusion site is visible. Use the provided Langenbeck elevator to elevate soft tissues along the length of the medial tibia for the intended plate being used.



Retrieve the appropriate length Medial Plate. Slide the plate under the soft tissues and along the bone to the appropriate position. It may be necessary to use the provided plate bending instrumentation to adjust the bend of the proximal portion of the plate to better fit patient anatomy. If provisional fixation of the plate to bone is preferred, olive wires are available to place within the screw holes or within the K-wire/Olive Wire holes in the plate. Confirm plate placement using fluoroscopy. Proceed with drilling and placement of screws along length of the plate according to surgeon preference and fracture pattern. .



NOTE: The plates should not extend to the same spot, but overlap to minimize stress risers proximally.




NOTE: The distal portion of the plate can be bent with provided instrumentation to better fit the shape of the medial malleolus and/or tibia. A slight external rotation bend from proximal to distal should help the plate sit along the length of the tibia.


MEDIAL PILON FUSION PLATE




All plate holes accept Ø2.7 mm, Ø3.5 mm or Ø4.2 mm locking or non-locking screws. Drill through the provided drill guides using the drill sized for the desired screw diameter.



Locking screws have the ability to be placed off-axis 15° in any direction. The cone end of the Easy Cone Drill Guide can be used to limit drilling to 15° in any direction. Drill in desired direction.




Measure screw length using a depth gauge. Insert the selected screw into the drilled hole in the tibia. Confirm screw position and length using fluoroscopy.



NOTE: Using a washer with a 3.5 compact or a 4.2 non locking screw may allow additional compression of the plate down to bone.

Refer to pages 14-15 for use of the drill guide and trocar for placement of proximal screws where the plate is under the soft tissue.

MEDIAL PILON FUSION PLATE




Continue with drilling and placement of screws along the length of the plate according to surgeon preference. If additional compression through the plate and into the talus is preferred, the screws in the two most distal rows can be angled to be placed into the talus.



MEDIAL PILON FUSION PLATE

Confirm plate and screw placement using fluoroscopy. Proceed to incision closure.





REMOVAL

Use provided drivers to remove all plate screws and any crossing screws. Remove the plate and confirm removal of all plates and screws under fluoroscopy.

INSTRUMENT CASE:

1. GORILLA® PILON FUSION CADDY

The Gorilla® Ankle Pilon Fusion Caddy contains the right and left anterolateral plates, Medial plates, Trocar, Drill Guide, Monster Drill and K-wire Guides, Backpack Guide, Backpack Drill Guides, Thread Attachment, Handle, Drill Guide, Curved Bone Fenestration Chisel, King Cobra™ Cartilage Removal Tool, Reciprocating Rasp, and Bone Fenestration Perforator are also included in this caddy.




2. ANKLE FRACTURE 360™ SCREW CADDY

The Ankle Fracture 360™ Screw Caddy accommodates Ø2.7 mm, Ø3.5 mm and Ø4.2 mm R3CON Locking and Non-locking Screws. As well as Gorilla R3CON Ø3.5 mm Compact Locking and Non-locking screws and Ø4.0 mm Mini Monster Screws.



3. GORILLA® ANKLE FRACTURE 360™ INSTRUMENT TRAYS

Drills, drill guides, centering guides, taps, drivers, plate bending instrumentation, K-wires, olive wires, handles and Ø4.0 Mini Monster Instrumentation.




4. ADDITIONAL GORILLA® CADDIES

Drills, drill guides, centering guides, taps, drivers, plate bending instrumentation, K-wires, olive wires, handles and Ø4.0 Mini Monster Instrumentation.



5. ANKLE FRACTURE 360™ INSTRUMENT CASE

Reduction clamps, joint preparation instrumentation and retractors are located at the bottom of the Ankle Fracture 360 Instrument Case Base.



INSTRUMENT CASE: _____

| Part # | Description | Quantity | Use |
|--------------|--|----------|------------|
| P53-440-L001 | Pilon Fusion Plate, Anterolateral, Short, Left | 1 | Single-use |
| P53-440-L002 | Pilon Fusion Plate, Anterolateral, Medium, Left | 1 | Single-use |
| P53-440-L003 | Pilon Fusion Plate, Anterolateral, Long, Left | 1 | Single-use |
| P53-440-R001 | Pilon Fusion Plate, Anterolateral, Short, Right | 1 | Single-use |
| P53-440-R002 | Pilon Fusion Plate, Anterolateral, Medium, Right | 1 | Single-use |
| P53-440-R003 | Pilon Fusion Plate, Anterolateral, Long, Right | 1 | Single-use |
| P53-441-0001 | Pilon Fusion Plate, Medial, Short | 1 | Single-use |
| P53-441-0002 | Pilon Fusion Plate, Medial, Medium | 1 | Single-use |
| P53-441-0003 | Pilon Fusion Plate, Medial, Long | 1 | Single-use |
| P51-926-1001 | Trocar | 1 | Reusable |
| P51-926-1003 | Drill Guide | 1 | Reusable |
| P51-926-2001 | Monster Drill Guide | 1 | Reusable |
| P51-926-1005 | Proximal Driver | 1 | Reusable |
| P51-926-2002 | Monster K-Wire Guide | 1 | Reusable |
| P51-926-4001 | King Cobra Cartilage Removal Tool | 1 | Reusable |
| P51-926-5001 | Backpack Guide, Right | 1 | Reusable |
| P51-926-5002 | Thread Attachment | 1 | Reusable |
| P51-926-5003 | Backpack Guide, Left | 1 | Reusable |
| P51-926-5004 | Backpack Drill Guide | 3 | Reusable |
| P51-926-6001 | Insertion Handle | 1 | Reusable |
| P51-926-7001 | Reciprocating Rasp | 1 | Single-use |
| P51-926-8001 | Depth Gage | 1 | Reusable |
| P99-150-0135 | Bone Fenestration Chisel, Curved | 1 | Reusable |
| P99-100-2010 | Bone Fenestration Perforator | 1 | Single-Use |

GORILLA ANKLE FRACTURE 360 PLATING SYSTEM - SCREW CADDY

| Part # | Description | Use |
|-------------------|---|------------|
| P50-353-27[08-54] | R3CON Locking Screw, Ø2.70 x 8-54 mm | Single-use |
| P50-453-27[08-75] | R3CON Non-Locking Screw, Ø2.70 x 8-75 mm | Single-use |
| P50-353-35[08-70] | R3CON Locking Screw, Ø3.50 x 8-70 mm | Single-use |
| P50-453-35[08-00] | R3CON Non-Locking Screw, Ø3.50 x 8-100 mm | Single-use |
| P50-353-42[08-70] | R3CON Locking Screw, Ø4.20 x 8-70 mm | Single-use |
| P50-453-42[08-70] | R3CON Non-Locking Screw, Ø4.20 x 8-70 mm | Single-use |
| P50-853-35[10-40] | Gorilla, R3CON, Non-Locking Screw, Compact Thread, Ø3.50 x 10-40 mm | Single-use |
| P50-753-35[10-40] | Gorilla, R3CON, Locking Screw, Compact Thread, Ø3.50 x 10-40 mm | Single-use |
| P20-140-0[24-60]L | Mini Monster Ø4.0 X 24-60 mm Headed Screw, Long Thread | Single-use |
| P20-040-DW00 | Ø4.0 Flat Washer, Mini Monster | Single-use |

*See page 6 for detailed screw size offerings

GORILLA ANKLE FRACTURE 360 PLATING SYSTEM - WIRE CADDY

| Part # | Description | Use |
|--------------|--|------------|
| P99-150-0005 | P28, K-Wire Hole Gauge & Ruler, SS | Reusable |
| P99-192-1615 | Ø1.60 mm x 15 cm Smooth, K-Wire | Single-use |
| P99-192-2015 | Ø2.0 mm x 15 cm Smooth, K-Wire | Single-use |
| P99-192-2020 | Ø2.0 mm x 20 cm Smooth, K-Wire | Single-use |
| P99-193-1615 | Ø1.60 mm x 15 cm Threaded, K-Wire | Single-use |
| P99-193-2015 | Ø2.0 mm x 15 cm Threaded, K-Wire | Single-use |
| P99-193-2020 | Ø2.0 mm x 20 cm Threaded K-Wire | Single-use |
| P99-200-1406 | Ø1.40 mm x 6 cm Olive Wire, Smooth | Single-use |
| P99-201-1406 | Ø1.40 mm x 6 cm Olive Wire, Threaded | Single-use |
| P99-250-1608 | Ø1.60 mm X 8 cm Olive Wire, Smooth, 316 LVM | Single-use |
| P99-250-1610 | Ø1.60 mm X 10 cm Olive Wire, Smooth, 316 LVM | Single-use |
| P99-251-1608 | Ø1.60 mm X 8 cm Olive Wire, Threaded, 316 LVM | Single-use |
| P99-251-1610 | Ø1.60 mm X 10 cm Olive Wire, Threaded, 316 LVM | Single-use |

GORILLA ANKLE FRACTURE 360 PLATING SYSTEM - INSTRUMENT TRAYS

| Part # | Description | Use |
|--------------|--|----------|
| P51-900-1004 | Gorilla R3CON, Drill Guide, Ø2.7 Threaded, Locking, SS, Long | Reusable |
| P51-900-1005 | Gorilla R3CON, Drill Guide, Ø3.5 Threaded, Locking, SS, Long | Reusable |
| P51-900-1006 | Gorilla R3CON, Drill Guide, Ø4.2 Threaded, Locking, SS, Long | Reusable |
| P51-901-5027 | Gorilla R3CON, Drill Guide, EZ Cone Ø2.7 mm, Long | Reusable |
| P51-901-5035 | Gorilla R3CON, Drill Guide, EZ Cone Ø3.5 mm, Long | Reusable |
| P51-901-5042 | Gorilla R3CON, Drill guide, EZ Cone Ø4.2 mm, Long | Reusable |
| P51-902-1004 | Gorilla R3CON, Drill Guide, Centering, Ø2.7 mm, SS, Long | Reusable |
| P51-902-1005 | Gorilla R3CON, Drill Guide, Centering, Ø3.5 mm, SS, Long | Reusable |
| P51-902-1006 | Gorilla R3CON, Drill Guide, Centering, Ø4.2 mm, SS, Long | Reusable |

GORILLA ANKLE FRACTURE 360 PLATING SYSTEM - INSTRUMENT TRAYS

| Part # | Description | Use |
|-----------------|--|----------|
| P51-902-1008 | Gorilla R3CON, Drill Guide, Centering, Ø3.5 mm Compact, SS, Long | Reusable |
| P51-903-3004 | Gorilla, R3CON, Drill Guide, Compression Slot, Ø2.7, Long | Reusable |
| P51-903-3005 | Gorilla, R3CON, Drill Guide, Compression Slot, Ø3.5, Long | Reusable |
| P51-903-3006 | Gorilla, R3CON, Drill Guide, Compression Slot, Ø4.2, Long | Reusable |
| P51-905-2027 | Gorilla R3CON, Drill/Over Drill Guide, Double Ended, Ø2.0 / Ø2.7, SS | Reusable |
| P51-905-2027-03 | Gorilla R3CON, Ø2.7 mm Drill Sleeve | Reusable |
| P51-905-2435 | Gorilla R3CON, Drill/Over Drill Guide, Double Ended, Ø2.4 / Ø3.5, SS | Reusable |
| P51-905-2435-03 | Gorilla R3CON, Ø3.5 mm Drill Sleeve | Reusable |
| P51-905-2842 | Gorilla R3CON, Drill / Over Drill Guide, Double Ended, Ø2.8 / Ø4.2, SS | Reusable |
| P51-905-2842-03 | Gorilla R3CON, Ø4.2 mm Drill Sleeve | Reusable |
| P51-913-2027 | Gorilla R3CON, Drill/Over Drill Guide, Double Ended, Ø2.0 / Ø2.7, SS | Reusable |
| P51-913-2027-03 | Gorilla R3CON, Ø2.7 mm Drill Sleeve | Reusable |
| P51-913-2435 | Gorilla R3CON, Drill/Over Drill Guide, Double Ended, Ø2.4 / Ø3.5, SS | Reusable |
| P51-913-2435-03 | Gorilla R3CON, Ø3.5 mm Drill Sleeve | Reusable |
| P51-913-2842 | Gorilla R3CON, Drill / Over Drill Guide, Double Ended, Ø2.8 / Ø4.2, SS | Reusable |
| P51-913-2842-03 | Gorilla R3CON, Ø4.2 mm Drill Sleeve | Reusable |
| P51-900-2002 | Gorilla, Ankle Fracture, Threaded Positioning Tower | Reusable |
| P51-910-1001 | Gorilla, Plate Bender, Threaded Bending Bar, SS | Reusable |
| P51-910-1002 | Gorilla, Plate Bender, Bending Irons, SS | Reusable |
| P51-910-1003 | Gorilla, Plate Bending Pliers | Reusable |
| P51-911-2700 | Gorilla Tap, Double Lead, Ø2.7 mm | Reusable |
| P51-911-3500 | Gorilla Tap, Double Lead, Ø3.5 mm | Reusable |
| P51-911-3501 | Gorilla Tap, Compact Thread, Ø3.5 mm | Reusable |
| P51-911-4200 | Gorilla Tap, Double Lead, Ø4.2 mm | Reusable |
| P99-100-2013 | P28, Drill, Gorilla Ø2.0 x 13 cm, Solid, A/O, SS | Reusable |
| P99-100-2416 | P28, Drill, Gorilla Ø2.4 x 16 cm, Solid A/O, SS | Reusable |
| P99-100-2418 | P28, Drill, Gorilla, Ø2.8 x 18 cm, Solid, A/O, SS | Reusable |
| P99-100-2713 | P28, Drill, Gorilla Ø2.7 x 110 mm, Solid, A/O, Overdrill, SS | Reusable |
| P99-100-2816 | P28, Drill, Gorilla, Ø2.8 x 16 cm, Solid, A/O, SS | Reusable |
| P99-100-3513 | P28, Drill, Gorilla Ø3.5 x 110 mm, Solid, A/O, Overdrill, SS | Reusable |
| P99-100-4214 | P28, Drill, Gorilla Ø4.2 x 120 mm, Solid, A/O, Overdrill, SS | Reusable |
| P99-191-AF10 | Gorilla, R3CON, Long, Driver, Solid, HX10 x 83 mm, SS | Reusable |
| P99-000-AOLG | Mini-AO, Ratchet Handle, Purple, Large | Reusable |
| P99-000-AOMN | Mini-AO, Ratchet Handle, Purple | Reusable |
| P99-150-0001 | P28, Screw Forceps, Ti | Reusable |
| P99-150-0013 | Sharp Tip Bone Reduction Clamp | Reusable |
| P99-150-0017 | Lobster Claw Clamp | Reusable |
| P99-150-0051 | 8 mm, Mini Hohmann Retractor | Reusable |
| P99-150-0052 | 16 mm, Mini Hohmann Retractor | Reusable |
| P99-150-0083 | Webber Spin Down Clamp | Reusable |

GORILLA ANKLE FRACTURE 360 PLATING SYSTEM - INSTRUMENT TRAYS

| Part # | Description | Use |
|--------------|---|----------|
| P99-150-0111 | P28, Depth Gauge Plate Screw, 90 mm, SS | Reusable |
| P99-150-0112 | P28, Percutaneous Depth Gauge, SS | Reusable |
| P99-150-0113 | Large, Lobster Claw Clamp | Reusable |
| P99-150-0114 | Lamina Spreader | Reusable |
| P99-150-0133 | Syndesmotic Clamp | Reusable |
| P99-150-0134 | Malleolar Reduction Clamp | Reusable |
| P99-150-1013 | (Long) Bone Holding Forceps | Reusable |
| P99-150-2003 | Langenbeck Elevator | Reusable |
| P99-150-2004 | Periosteal Elevator | Reusable |
| P99-157-0810 | Ribbon Retractor 1.0" X 8" | Reusable |
| P99-157-0850 | Ribbon Retractor 0.5" X 8" | Reusable |
| P99-191-SL10 | Screw Retention Sleeve | Reusable |
| P99-150-0038 | Tenaculum Clamp | Reusable |

Indications, Contraindications and Warnings

INDICATIONS FOR USE (GORILLA®)

The Baby Gorilla®/Gorilla® Bone Plates and Bone Screws of the Baby Gorilla®/Gorilla® Plating System are indicated for use in stabilization and fixation of fractures or osteotomies; intra and extra articular fractures, joint depression, and multi-fragmentary fractures; revision procedures, joint fusion and reconstruction of small bones of the toes, feet and ankles including the distal tibia, talus, and calcaneus, as well as the fingers, hands, and wrists. The system can be used in both adult and pediatric patients. Specific examples include:

Forefoot:

- Arthrodesis of the first metatarsalcuneiform joint (Lapidus Fusion)
- Metatarsal or phalangeal fractures and osteotomies
- Lesser metatarsal shortening osteotomies (e.g. Weil)
- Fifth metatarsal fractures (e.g. Jones Fracture)

Mid/Hindfoot:

- LisFranc Arthrodesis and/or Stabilization
- 1st (Lapidus), 2nd, 3rd, 4th, and 5th Tarsometatarsal (TMT) Fusions
- Intercuneiform Fusions
- Navicular-Cuneiform (NC) Fusion
- Talo-Navicular (TN) Fusion
- Calcaneo-Cuboid (CC) Fusion
- Subtalar Fusion
- Medial Column Fusion
- Cuneiform Fracture
- Cuboid Fracture
- Navicular Fracture

Ankle:

- Lateral Malleolar Fractures
- Syndesmosis Injuries
- Medial Malleolar Fractures and Osteotomies
- Bi-Malleolar Fractures
- Tri-Malleolar Fractures
- Posterior Malleolar Fractures
- Distal Anterior Tibia Fractures
- Vertical Shear Fractures of the Medial Malleolus
- Pilon Fractures
- Distal Tibia Shaft Fractures
- Distal Fibula Shaft Fractures
- Distal Tibia Periarticular Fractures
- Medial Malleolar Avulsion Fractures
- Lateral Malleolar Avulsion Fractures
- Tibiotalocalcaneal Joint Arthrodesis
- Tibiotalar Joint Arthrodesis
- Tibiocalcaneal Arthrodesis
- Supramalleolar Osteotomy
- Fibular Osteotomy

First metatarsal osteotomies for hallux valgus correction including:

- Opening base wedge osteotomy

- Closing base wedge osteotomy
- Crescentic Osteotomy
- Proximal Osteotomy (Chevron and Rotational Oblique)
- Distal Osteotomy (Chevron/Austin)

Arthrodesis of the first metatarsophalangeal joint (MTP) including:

- Primary MTP Fusion due to hallux rigidus and/or hallux valgus
- Revision MTP Fusion
- Revision of failed first MTP Arthroplasty implant

Flatfoot:

- Lateral Column Lengthening (Evans Osteotomy)
- Plantar Flexion Opening Wedge Osteotomy of the Medial Cuneiform (Cotton Osteotomy)
- Calcaneal Slide Osteotomy

Charcot:

- Medial column fusion (talus, navicular, cuneiform, metatarsal) for neuropathic osteoarthropathy (Charcot)
- Lateral column fusion (calcaneus, cuboid, metatarsal) for neuropathic osteoarthropathy (Charcot)

In addition, the non-locking, titanium screws and washers are indicated for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, fracture repair and fracture fixation, appropriate for the size of the device.

CONTRAINDICATIONS

Use of the Baby Gorilla®/Gorilla® Plating System is contraindicated in cases of inflammation, cases of active or suspected sepsis/infection and osteomyelitis; or in patients with certain metabolic diseases.

All applications that are not defined by the indications are contraindicated. In addition, surgical success can be adversely affected by:

- Acute or chronic infections, local or systemic
- Vascular, muscular or neurological pathologies that compromise the concerned extremity
- All concomitant pathologies that could affect the function of the implant
- Osteopathies with reduced bone substance that could affect the function of the implant
- Any mental or neuromuscular disorder that could result in an unacceptable risk of failure at the time of fixation or complications in post-operative treatment
- Known or suspected sensitivity to metal
- Corpulence; an overweight or corpulent patient can strain the implant to such a degree that stabilization or implant failure can occur
- Whenever the use of the implant comes into conflict with the anatomical structures of physiological status

Other medical or surgical pre-conditions that could compromise the potentially beneficial procedure, such as:

- The presence of tumors
- Congenital abnormalities
- Immunosuppressive pathologies
- Increased sedimentation rates that cannot be explained by other pathologies

Increased leukocyte (WBC) count

- Pronounced left shift in the differential leukocyte count

POTENTIAL COMPLICATIONS AND ADVERSE REACTIONS

In any surgical procedure, the potential for complications and adverse reactions exist. The risks and complications with these implants include:

- Loosening, deformation or fracture of the implant
- Acute post-operative wound infections and late infections with possible sepsis
- Migration, subluxation of the implant with resulting reduction in range of movement
- Fractures resulting from unilateral joint loading
- Thrombosis and embolism
- Wound hematoma and delayed wound healing
- Temporary and protracted functional neurological perturbation
- Tissue reactions as the result of allergy or foreign body reaction to dislodged particles
- Corrosion with localized tissue reaction and pain
- Pain, a feeling of malaise or abnormal sensations due to the implant used
- Bone loss due to stress shielding

All possible complications listed here are not typical of Paragon 28®, Inc. products but are in principle observed with any implant. Promptly inform Paragon 28®, Inc. as soon as complications occur in connection with the implants or surgical instruments used. In the event of premature failure of an implant in which a causal relationship with its geometry, surface quality or mechanical stability is suspected, please provide Paragon 28®, Inc. with the explant(s) in a cleaned, disinfected and sterile condition. Paragon 28®, Inc. cannot accept any other returns of used implants. The surgeon is held liable for complications associated with inadequate asepsis, inadequate preparation of the osseous implant bed in the case of implants, incorrect indication or surgical technique or incorrect patient information and consequent incorrect patient behavior.

WARNINGS AND PRECAUTIONS

- Re-operation to remove or replace implants may be required at any time due to medical reasons or device failure. If corrective action is not taken, complications may occur.
- Use of an undersized plate or screw in areas of high functional stresses may lead to implant fracture and failure.
- Plates and screws, wires, or other appliances of dissimilar metals should not be used together in or near the implant site.

- The implants and guide wires are intended for single use only.
- Instruments, guide wires and screws are to be treated as sharps.
- Do not use other manufacturer's instruments or implants in conjunction with the Baby Gorilla®/Gorilla® Plating System.
- If a stainless steel Gorilla® R3LEASE™ Screw is used, it may only be used standalone.
- The device should only be used in pediatric patients where the growth plates have fused or in which active growth plates will not be crossed by the system implants or instrumentation.

MR SAFETY INFORMATION

The Baby Gorilla®/Gorilla® Plating System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Baby Gorilla®/Gorilla® Plating System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Indications, Contraindications and Warnings

INDICATIONS FOR USE (MONSTER®)

The Monster® Screw System is indicated for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, ligament fixation, fracture repair and fracture fixation, appropriate for the size of the device. Specific examples include:

Fractures and Osteotomies

- Fractures of the tarsals, metatarsals and other fractures of the foot (i.e. LisFranc)
- Avulsion fractures and fractures of the 5th metatarsal (i.e. Jones Fracture)
- Talar fractures
- Ankle fractures
- Navicular fractures
- Fractures of the fibula, malleolus, and calcaneus
- Metatarsal and phalangeal osteotomies
- Weil osteotomy
- Calcaneal osteotomy

Hallux Valgus Correction

- Fixation of osteotomies (i.e. Akin, Scarf, Chevron)
- Interphalangeal (IP) arthrodesis
- Proximal, midshaft, or distal osteotomy
- Lapidus arthrodesis

CONTRAINDICATIONS (CONTINUED)

All applications that are not defined by the indications are contraindicated. In addition, surgical success can be adversely affected by:

- Acute or chronic infections, local or systemic
- Vascular, muscular or neurological pathologies that compromise the concerned extremity
- All concomitant pathologies that could affect the function of the implant
- Osteopathies with reduced bone substance that could affect the function of the implant
- Any mental or neuromuscular disorder that could result in an unacceptable risk of failure at the time of fixation or complications in post-operative treatment
- Known or suspected sensitivity to metal
- Corpulence; an overweight or corpulent patient can strain the implant to such a degree that stabilization or implant failure can occur
- Whenever the use of the implant comes into conflict with the anatomical structures of physiological status

Other medical or surgical pre-conditions that could compromise the potentially beneficial procedure, such as:

- The presence of tumors
- Congenital abnormalities
- Immunosuppressive pathologies
- Increased sedimentation rates that cannot be explained by other pathologies
- Increased leukocyte (WBC) count
- Pronounced left shift in the differential leukocyte count

POTENTIAL COMPLICATIONS AND ADVERSE REACTIONS

In any surgical procedure, the potential for complications and adverse reactions exist. The risks and complications with these implants include:

- Loosening, deformation or fracture of the implant
- Acute post-operative wound infections and late infections with possible sepsis
- Migration, subluxation of the implant with resulting reduction in range of movement
- Fractures resulting from unilateral joint loading
- Thrombosis and embolism
- Wound hematoma and delayed wound healing
- Temporary and protracted functional neurological perturbation
- Tissue reactions as the result of allergy or foreign body reaction to dislodged particles.
- Corrosion with localized tissue reaction and pain
- Pain, a feeling of malaise or abnormal sensations due to the implant used
- Bone loss due to stress shielding

All possible complications listed here are not typical of Paragon 28®, Inc. products but are in principle observed with any implant. Promptly inform Paragon 28® as soon as complications occur in connection with the implants or surgical instruments used. In the event of premature failure of an implant in which a causal relationship with its geometry, surface quality or mechanical stability is suspected, please provide Paragon 28® with the explant(s) in a cleaned, disinfected and sterile condition. Paragon 28® cannot accept any other returns of used implants. The surgeon is held liable for complications associated with inadequate asepsis, inadequate preparation of the osseous implant bed in the case of implants, incorrect indication or surgical technique or incorrect patient information and consequent incorrect patient behavior.

WARNINGS AND PRECAUTIONS

- Re-operation to remove or replace implants may be required at any time due to medical reasons or device failure. If corrective action is not taken, complications may occur.
- Use of an undersized screw in areas of high functional stresses may lead to implant fracture and failure.
- Plates and screws, wires, or other appliances of dissimilar metals should not be used together in or near the implant site.
- The implants and guide wires are intended for single use only. Re-use may cause product failure and could lead to disease transmission.
- Instruments, guide wires and screws are to be treated as sharps.
- Do not use other manufacturer's instruments or implants in conjunction with the Monster® Screw System.

MR SAFETY INFORMATION

The Monster® Screw System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Monster® Screw System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

NOTES:

This image shows a single sheet of white paper with horizontal blue or grey ruling lines, typical of notebook paper. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



SURGICAL TECHNIQUE GUIDE

Gorilla® Pilon Fusion System

P51-STG-0015 RevA [2025-07-28]

™Trademarks and ®Registered Marks of Paragon 28®, Inc.
© Copyright 2025 Paragon 28®, Inc. All rights reserved.
Patents: www.paragon28.com/patents

Paragon 28®, Inc.
14445 Grasslands Dr.
Englewood, CO 80112 USA
(855) 786-2828 📞



Disclaimer:

The purpose of the Gorilla® Pilon Fusion Plating System Surgical Technique Guide is to demonstrate use of the Gorilla® Pilon Fusion in the Gorilla® Pilon Fusion Plating System. Although various methods can be employed for this procedure, the fixation options demonstrated were chosen for simplicity of explanation and demonstration of the unique features of our device. Federal law (U.S.A.) restricts this device to sale and use by, or on order of, a physician.