




SURGICAL TECHNIQUE GUIDE

MTP Arthrodesis

Gorilla[®] MTP Plating System



PRODUCT DESCRIPTION

The Paragon 28® Gorilla® MTP Plating System was designed to provide surgeons versatility in plate selection for MTP Arthrodesis procedures, assistance with joint preparation and the ability to position the hallux in all 3 planes during temporary fixation. The system has 32 anatomically contoured plating options to address primary, revision and graft-spanning arthrodesis. The primary and short arthrodesis plates are offered in 0°, 5° and 10° of dorsiflexion to account for differences in anatomy. The revision arthrodesis plates are thicker and provide more holes proximally to avoid previous screw placement. The graft-spanning plates allow for placement of the Paragon 28® PRESERVE™ MTP Length Restoring Graft. All Gorilla® MTP plate holes accommodate Gorilla® R3CON™ 2.7, 3.5 and 4.2 mm locking and non-locking screws.


The instrumentation provided in the Gorilla® MTP Plating System was designed to address joint preparation while facilitating compression and/or stability at the arthrodesis site. The patented Precision® Guide mates with the desired MTP plate to provide five trajectories of guide wire paths to allow for insertion of a 3.0 or 3.5 mm cannulated Mini-Monster® crossing screw across the arthrodesis site, while avoiding on-axis locking and non-locking plate screws within the construct. Spin Guard Reamers are available in the set to assist the surgeon in cartilage removal and joint preparation prior to plate and screw fixation.

PLATE OFFERING

Primary Plates

- Left and Right Side Specific

- 1.3 mm thickness
- Available in small, medium and large
- 3 holes in proximal phalanx



- Available Plate Angles:



Short Plate

- Left and Right Side Specific

- 1.3 mm thickness
- 2 holes in proximal phalanx




- Available Plate Angles:



Revision Plates

- Left and Right Side Specific

- 1.6 mm thickness
- Available in small, medium and large



- Plate Angle:



Graft-Spanning Plate

- Left and Right Side Specific

- 1.6 mm thickness
- Built to span PRESERVE™ MTP Graft



- Plate Angle:



FEATURED INSTRUMENTS



Patented Female Spin Guard



Patent Pending Male Spin Guard



Female Reamer



Male Reamer

- Available in 17 mm, 19 mm, 21 mm and 23 mm diameters
- Reamers allow for cup and cone engagement of the 1st MTP, providing the ability to adjust hallux position in any plane without re-cutting
 - Helps to remove less bone than flat cuts
- Spin Guard designed to protect adjacent bones while reaming



1.6 x 150 mm K-wire

Precision Guide MTP



1.0 mm Precision Guide K-wire Sleeve



Mini-Monster® 3.0 mm screw



1.2 mm Precision Guide K-wire Sleeve



Mini-Monster® 3.5 mm screw



Precision Guide Set Screw



Precision Guide Arm

- Patented
- Allows a lag screw to cross the joint while avoiding the locking and non-locking plate screws in the construct

ANCILLARY IMPLANTS


PRESERVE™ MTP Length Restoring Graft

- Can help to restore length in cases of a short 1st metatarsal or revision procedure
- Available in 19 mm diameter with lengths of 5, 8, 10, 15 and 20 mm
- Available in 21 mm diameter with lengths of 5, 8 and 10 mm



INCISION/EXPOSURE

Supine patient positioning with fluoroscopy available is recommended for this procedure. A dorsomedial incision over the 1st metatarso-phalangeal joint is recommended. Soft tissue dissection is continued to expose the 1st metatarsophalangeal joint. Release the soft tissue to obtain exposure of the articular surfaces of the 1st metatarsal head and hallux proximal phalanx base.



JOINT PREPARATION

After exposing the joint surfaces, remove any osteophytes surrounding the joint using a saw or osteotome. The method of cartilage resection is according to surgeon preference. The method for using the Paragon 28 MTP reamers is described below:

PREPARATION OF THE 1st METATARSAL HEAD:


A 1.6 mm K-wire is inserted down the central shaft of the 1st metatarsal head, from distal to proximal. Select the female reamer size based on the 1st metatarsal head size, with the reamer as wide or slightly wider than the diameter of the cartilage covering the 1st metatarsal head.



Connect the matching female spin guard to the end of the female reamer. Attach the construct to a powered driver and slide over the 1.6 mm K-wire. Begin motion of the reamer prior to making contact with the 1st metatarsal head.

Remove all of the cartilage on the 1st metatarsal head, using the reamer in a pulsing motion to facilitate cartilage removal, if necessary. If the outer cartilage remains, go up a reamer size. If the reamer is too large, go down a size. Remove cartilage until bleeding subchondral bone is observed and take care not to over-shorten the 1st metatarsal.

Take note of the last reamer size used, as this will be the size designated for reaming the proximal phalanx of the hallux. When finished, remove the K-wire from the central shaft of the 1st metatarsal. A rongeur or curette can be employed to resect any cartilage or rough edges from the 1st metatarsal head joint surface.



JOINT PREPARATION


PREPARATION OF THE PROXIMAL PHALANX BASE:

Insert the same 1.6 mm K-wire down the central shaft of the hallux proximal phalanx, from proximal to distal. The convex male reamer is selected to match the last size of reamer used on the 1st metatarsal head.




Attach the matching male spin guard to the end of the male reamer. The male reamer is placed over the K-wire and is secured to the powered driver. Begin motion of the reamer prior to contact with the proximal phalanx.

Remove all of the cartilage from the base of the hallux proximal phalanx until bleeding subchondral bone is observed. Remove the K-wire from the base of the hallux proximal phalanx. A rongeur or curette can be employed to resect any additional cartilage from the edges or joint surface. Subchondral bone preparation can be performed following reaming using the Paragon 28 subchondral drill or surgeon's preferred technique. If used, bone grafting material or a PRESERVE MTP Length Restoring Graft can be inserted at this time.



TEMPORARY FIXATION




The correct alignment of the hallux can be determined with the assistance of the Paragon 28 Foot Plate. A K-wire is recommended to be inserted from proximal medial to distal lateral to serve as temporary fixation.

PLATE SELECTION

An appropriately sized "Left" or "Right" MTP arthrodesis plate is selected at this time. The laser markings at the central aspect of the plate should align with the joint.


TIP: If a thicker plate is desired for a patient that may have difficulty with restricted weight-bearing following the procedure, a small revision plate (thickness 1.6 mm) may be used in lieu of a medium primary plate (thickness 1.3 mm), as the length is comparable.



PERMANENT FIXATION


Insert the Precision Guide set screw into the corresponding “R” or “L” hole of the Precision Guide Arm such that the set screw is aligned with the larger of the two central holes of the MTP Plate and the Precision Guide Arm is oriented medially. The small peg on the underside of the Precision Guide mates with the smaller of the two central holes of the MTP Plate. Rotate the knob clockwise to secure the Precision Guide set screw to the MTP plate.

NOTE: Use of the Precision Guide MTP is optional. It may be used as demonstrated below, or it may be attached to the plate after plate screw fixation. In the latter case, a fully threaded crossing screw is recommended.



Secure the MTP plate (with the attached Precision Guide) to the bone using olive wires. One olive wire should be placed in the proximal aspect of the compression slot. Insert the guide wire sleeve for the selected screw size into the Precision Guide Arm. If soft tissue dissection is not performed around the area of the guide wire sleeve, a stab incision with blunt soft tissue separation can be made in the skin prior to driving the guide wire through the arthrodesis site.


TIP: The central hole of the Precision Guide generally works for most patients. For larger patients, the plantar hole may work best. The inner and outer holes may help prevent skiving if that is a concern.




Insert the guide wire for the selected crossing screw size into the guide wire sleeve such that it crosses the arthrodesis site at a desired location.

The Precision Guide can be used during initial permanent fixation. In this case (as shown), no further compression across the arthrodesis site should be attempted following crossing screw placement, as this will result in loosening of the crossing screw.

PERMANENT FIXATION




The position and length of the guide wire is confirmed using fluoroscopy. When correct, remove the Precision Guide from the plate by rotating the set screw in a counter-clockwise manner and sliding the Precision Guide Arm off of the guide wire.




The drill and drill guide for the selected crossing screw diameter are slid over the guide wire and drilling is performed.



Countersinking for the headed screw is performed. If using a headless screw, countersink after measuring. The depth gauge is used to determine screw length.



The selected crossing screw is inserted over the guide wire into the proximal phalanx.



Complete screw insertion.



TIP: If a surgeon prefers to not use a crossing screw, or if a crossing screw is unable to be placed, eccentric drilling of the compression slot can be performed and compression can be achieved in this manner or by external compression of the joint prior to locking screw placement.

SURGICAL TECHNIQUE GUIDE:


MTP ARTHRODESIS

PERMANENT FIXATION


Remove the guide wire.



Insert a threaded drill tower into one of the distal screw holes corresponding to desired plate screw diameter.



Drill using the drill corresponding to the desired plate screw diameter.




Screw length can be measured using the provided depth gauge or by measuring off of the drill using the drill guide. Insert the plate screw using the provided driver. It is advised to avoid final tightening of the distal plate screws into the locked position until proximal plate screws are inserted.




Insert an oblong drill guide into the plate compression slot. The oblong drill guide can be reversed with the arrow pointing away from the joint (as shown) to insert a non-locking screw into the compression slot without creating compression.

TIP: If the surgeon prefers to create compression via the compression slot, insert distal plate (screw)s first. Use the compression drill guide with the arrow pointing towards the joint to eccentrically drill. Place a non-locking screw to achieve compression. The cross screw would be placed last in this scenario.




Drill using the drill corresponding to the desired screw diameter.



Insert the screw using the driver provided.

Insert screws in the remaining screw holes, as desired. Confirm final screw length and placement using fluoroscopy.



CLOSURE

Proceed to incision closure or concomitant procedures at this time.

MTP Caddy

The Gorilla® MTP Plating Caddy includes 32 varieties of Gorilla® MTP plates, MTP reamers in 4 sizes and 1.6 mm K-wires for use with the MTP reamers. The Precision® Guide MTP is also located in this caddy.



Gorilla® R3CON Instrument Caddy

Drills, drill guides, centering guides, olive wires, plate benders, drivers, K-wires and a depth gauge are located in the Gorilla® R3CON Instrument Caddy.

Additional Gorilla® Caddies

The Gorilla® Case has room for additional Gorilla® Plate Caddies or PRESERVE™ Allograft caddies that may be needed for procedures performed in addition to an MTP Arthrodesis case.




Mini-Monster® Screw Caddy

The Gorilla® Case can accommodate one Mini-Monster® Screw Caddy for a 3.0 mm or 3.5 mm Mini-Monster® screw to be used during an MTP Arthrodesis procedure.



Gorilla® Case



Gorilla® R3CON Screw Optionality

The Gorilla® R3CON screw length options for both locking and non-locking screws are as follows:

2.7 mm	1 mm increments, 8-20 mm	
2.7 mm	2 mm increments, 22-40 mm	
3.5 mm	2 mm increments, 10-50 mm	
4.2 mm	2 mm increments, 10-50 mm	
4.2 mm	5 mm increments, 55-70 mm	

Gorilla® R3CON Instruments

The Caspar Compression/Distraction device, osteotomes, baby Bennet retractors, bone reduction clamps, periosteal elevator, cartilage removal device, pin distractor and handles are located at the bottom of the Gorilla® Case.

GORILLA R3CON PLATING SYSTEM - CASE INSTRUMENTS:

Part #	Description	Use
P99-000-AOMN	Mini-AO, Ratchet Handle, Purple	Reusable
P99-000-AOSS	Mini-AO, Non-Ratcheting Handle Streamline, Cannulated	Reusable
P99-150-0004	San Gio™ Baby Bennetts (2 Each)	Reusable
P99-150-0013	Sharp Tip Bone Reduction Clamp	Reusable
P99-150-0015	Medium HINTERMANN Distractor	Reusable
P99-150-0017	Lobster Claw	Reusable
P99-150-0018	Periosteal Elevator	Reusable
P99-150-0019	Lapidus Ronjeur	Reusable
P99-150-0034	The Honey Badger™ (Cartilage Removal Tool)	Reusable
P99-150-0040	Osteotome, Straight, 6mm	Reusable
P99-150-0041	Osteotome, Straight, 13mm	Reusable
P99-150-0042	Osteotome, Straight, 19mm	Reusable
P99-150-0140	Osteotome, Curved, 6mm	Reusable
P99-150-0141	Osteotome, Curved, 13mm	Reusable
P99-150-0142	Osteotome, Curved, 19mm	Reusable
P99-900-1005	Foot Plate, Plantigrade Guide	Reusable
P99-150-0011R	Caspar Type Compression/Distraction Device	Reusable

GORILLA R3CON PLATING SYSTEM - INSTRUMENT CADDY:

Part #	Description	Use
P51-900-1001	Gorilla R3CON, Drill Guide, 2.7 Threaded, Locking, SS	Reusable
P51-900-1002	Gorilla R3CON, Drill Guide, 3.5 Threaded, Locking, SS	Reusable
P51-900-1003	Gorilla R3CON, Drill Guide, 4.2 Threaded, Locking, SS	Reusable
P51-901-3101	Gorilla R3CON, Drill Guide, 2.7 / 3.5 Easy Guide, SS	Reusable
P51-901-3102	Gorilla R3CON, Drill Guide, 2.7 / 3.5 Cone Guide, SS	Reusable
P51-901-3142	Gorilla R3CON, Drill Guide, 4.2 Easy / Cone Guide, SS	Reusable
P51-902-1001	Gorilla R3CON, Drill Guide, Centering, 2.7mm, SS	Reusable
P51-902-1002	Gorilla R3CON, Drill Guide, Centering, 3.5mm, SS	Reusable
P51-902-1003	Gorilla R3CON, Drill Guide, Centering, 4.2mm, SS	Reusable
P51-903-3001	Gorilla R3CON, Drill Guide, Compression Slot, 2.7mm, SS	Reusable
P51-903-3002	Gorilla R3CON, Drill Guide, Compression Slot, 3.5mm, SS	Reusable
P51-903-3003	Gorilla R3CON, Drill Guide, Compression Slot, 4.2mm, SS	Reusable
P51-904-3001	Gorilla R3CON, Countersink / Punch, SS	Reusable
P51-905-1001	Gorilla, Over Drill Guide 2.7 / 3.5-4.2, SS	Reusable
P51-910-1001	Gorilla, Plate Bender, Threaded Bending Bar, SS	Reusable
P51-910-1002	Gorilla, Plate Bender, Bending Irons, SS	Reusable
P51-955-0001	Gorilla, Caspar Insert, Plate Engaging Compression Tool, SS	Reusable
P99-100-2009	P28, Drill, Bone Fenestration Perforator, A/O, SS	Reusable
P99-100-2011	P28, Drill, Gorilla 2.0 x 110 mm, Solid, A/O, SS	Reusable

GORILLA R3CON PLATING SYSTEM - INSTRUMENT CADDY:

Part #	Description	Use
P99-100-2414	P28, Drill, Gorilla 2.4 x 140 mm, Solid, A/O, SS	Reusable
P99-100-2814	P28, Drill, Gorilla 2.8 x 140 mm, Solid, A/O, SS	Reusable
P99-100-2711	P28, Drill, Gorilla 2.7 x 110 mm, Solid, A/O, Overdrill, SS	Reusable
P99-100-3511	P28, Drill, Gorilla 3.5 x 110 mm, Solid, A/O, Overdrill, SS	Reusable
P99-100-4212	P28, Drill, Gorilla 4.2 x 120 mm, Solid, A/O, Overdrill, SS	Reusable
P99-150-0001	P28, Screw Forceps, Ti	Reusable
P99-150-0005	P28, K-Wire Hole Gauge & Ruler, SS	Reusable
P99-150-0014	P28, Depth Gauge Plate Screw, 70mm, SS	Reusable
P99-150-0035	P28, Chisel, Bone Fenestration Chisel, SS	Reusable
P99-191-TR10	P28, Driver, HX10 Gorilla R3CON Screw, Solid, AO, SS	Reusable
P99-192-2015	P28, K-wire, 2.0 x 150 mm, Smooth, SS	Single-use
P99-192-2315	P28, K-wire, 2.3 x 150 mm, Smooth, SS	Single-use
P99-200-1406	P28, Olive Wire, 1.4 mm, Smooth, SS	Single-use
P99-201-1406	P28, Olive Wire, 1.4 mm, Threaded, SS	Single-use
P99-202-1810	P28, Olive Wire, 1.8 mm, Short Tip, SS	Single-use

GORILLA R3CON PLATING SYSTEM - R3CON SCREW CADDY:

Part #	Description	Use
P50-153-WS00	Screw Flat/Dome Washer, 2.7, 3.5 & 4.2mm	Single-use
P50-353-27[08-40]	R3CON™ Locking Plate Screw, 2.7 x 8-40mm	Single-use
P50-353-35[10-50]	R3CON™ Locking Plate Screw, 3.5 x 10-50mm	Single-use
P50-353-42[10-70]	R3CON™ Locking Plate Screw, 4.2 x 10-70mm	Single-use
P50-453-27[08-40]	R3CON™ Non-Locking Plate Screw, 2.7 x 8-40mm	Single-use
P50-453-35[10-50]	R3CON™ Non-Locking Plate Screw, 3.5 x 10-50mm	Single-use
P50-453-42[10-70]	R3CON™ Non-Locking Plate Screw, 4.2 x 10-70mm	Single-use

*See page 10 for screw size increments

Refer to www.paragon28.com/ifus for the complete and most current instructions for use document.

INSTRUCTIONS FOR USE: MONSTER®/MINI-MONSTER® SCREW SYSTEM

Indications, Contraindications, Warnings and Precautions relevant to the Monster® Screw System are contained in the Instructions for Use document of the Monster® Screw System P20-IFU-3001.

	A patient with the Paragon 28® Monster® Screw System may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Name/Identification of Device	Paragon 28® Monster® Screw System	
Nominal value(s) of Static Magnetic Field [T]	1.5 T or 3 T	
Maximum Spacial Field Gradient [T/m and gauss/cm]	30 T/m (3000 gauss/cm)	
RF Excitation	Circularly Polarized (CP)	
RF Transmit Coil Type	Whole body transmit coil, Head RF transmit-receive coil	
Maximum Whole Body SAR [W/kg]	2.0 W/kg (Normal Operating Mode)	
Limits on Scan Duration	All anatomical regions can be safely scanned under the following conditions: 2.0 W/kg whole body average SAR for 5 minutes of continuous RF (a sequence or back to back series/scan without breaks) with a 20 minute cooling period between scans for an hour long scanning session	
	Scanning of the knees and all anatomy superior to the knees can be safely scanned under the following conditions: 2.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)	
MR Image Artifact	The presence of this implant may produce an image artifact of 20 mm	
If information about a specific parameter is not included, there are no conditions associated with that parameter.		

Refer to www.paragon28.com/ifus for the complete and most current instructions for use document.

INSTRUCTIONS FOR USE: BABY GORILLA®/GORILLA® PLATING SYSTEM

Indications, Contraindications, Warnings and Precautions relevant to the Baby Gorilla®/Gorilla® Plating System are contained in the Instructions for Use document of the Baby Gorilla®/Gorilla® Plating System P51-IFU-3001.

 A patient with the Paragon 28® Baby Gorilla®/Gorilla® Plating System may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Name/Identification of Device	Paragon 28® Baby Gorilla®/Gorilla® Plating System
Nominal value(s) of Static Magnetic Field [T]	1.5 T or 3 T
Maximum Spacial Field Gradient [T/m and gauss/cm]	30 T/m (3000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Whole body transmit coil, Head RF transmit-receive coil
Operating Mode	Normal Operating Mode
Maximum Whole Body SAR [W/kg]	2.0 W/kg (Normal Operating Mode)
Limits on Scan Duration	2.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)
MR Image Artifact	The presence of this implant may produce an image artifact
If information about a specific parameter is not included, there are no conditions associated with that parameter.	

NOTES:

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

NOTES:

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

NOTES:

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Gorilla[®] MTP Plating System

PATENTED, DESIGNED & EXCLUSIVELY DISTRIBUTED BY



P51-STG-3002 RevC [2026-06-25]

© Copyright 2026 Zimmer Biomet[®]. All rights reserved.
Patents: www.paragon28.com/patents

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

Paragon 28, Inc. 
14445 Grasslands Dr.
Englewood, CO 80112 USA
(855) 786-2828

Australian Sponsor
Zimmer Biomet Ltd
(A.C.N: 096 480 992)(A.B.N: 96 096 480 992)
Level 4, 12 Narabang Way
BELROSE NSW 2085
AUSTRALIA
Phone: (02) 9483 5400

DISCLAIMER

Zimmer Biomet does not practice medicine. This technique was developed in conjunction with healthcare professionals. This document is intended for surgeons and is not intended for laypersons. Each surgeon should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training surgeons have received. As with all surgical procedures, the product/technique used in each case will depend on the surgeon's medical judgment as the best treatment for each patient. Results will vary based on health, weight, activity and other variables. Not all patients are candidates for this product and/or procedure. Caution: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only. ©2025 Zimmer Biomet.

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the instructions for use or contact your local representative; visit www.paragon28.com for additional product information. This material is intended for health care professionals. Distribution to any other recipient is prohibited.

www.Paragon28.com