



# Paragon<sup>20</sup>

a  ZIMMER BIOMET company



 **GORILLA**<sup>®</sup>  
R3CON PLATING SYSTEM

## SURGICAL TECHNIQUE GUIDE

Gorilla<sup>®</sup> Evans Calcaneal Osteotomy

## PRODUCT DESCRIPTION:

The Paragon 28® Gorilla® Plating System provides two distinct plating options for Evans calcaneal osteotomies:

### **HEvans™ Plates and BOW & ARROW™ Plates**

**HEvans plates** use a 3 screw fixation construct to support the graft placed into the osteotomy and helps to prevent graft subsidence with subsequent loss of calcaneal length seen in the early post-operative period due to graft resorption. It is low profile (1.1 mm), with a tapering of the thickness to 0.5 mm toward the end to prevent soft tissue irritation.

The design allows for a single proximal screw to be placed in the calcaneus with two distal screws placed in the anterior fragment of the calcaneus. Dual point fixation at the distal fragment mitigates the prospect of rotation of the distal fragment with the calcaneocuboid anatomy. The posterior aspect of the plate is tapered to prevent peroneal tendon irritation in the location of the plate.

**BOW & ARROW plates** have a tapered “BOW” at the back of the plate helps to evenly distribute vertical pressure along the graft due to the anatomically matched taper, with the point decreasing risk of tearing or attenuation of the lateral aspect of the long plantar ligament.

Anatomic considerations were at the forefront of the Gorilla® BOW & ARROW® Evans Plate design. The low profile (thickness 1.15 mm) nature of the plate helps to minimize risk of soft tissue irritation. The ramped plantar aspect of the plate helps to minimize interference with the peroneal tendons and sural nerve. Locking screw technology allows for screw heads that are flush or below flush to the plate.

## CONTENTS:


|                  |                                                           |
|------------------|-----------------------------------------------------------|
| <b>Section 1</b> | <b>GORILLA® HEVANS™ AND BOW &amp; ARROW™ PLATES</b>       |
|                  | <b>PLATE OFFERINGS.....3</b>                              |
|                  | <b>PLATE FEATURES .....3-4</b>                            |
| <b>Section 2</b> | <b>REPAIR TECHNIQUE</b>                                   |
|                  | <b>INCISION/EXPOSURE .....5</b>                           |
|                  | <b>OSTEOTOMY AND JOINT DISTRACTION..... 5-7</b>           |
|                  | <b>REPAIR USING HEVANS™ PLATE.....8-9</b>                 |
|                  | <b>REPAIR USING BOW AND ARROW™ PLATE.....10-12</b>        |
| <b>Section 3</b> | <b>CADDY, INDICATIONS, AND WARNINGS</b>                   |
|                  | <b>CADDY LAYOUT .....13-15</b>                            |
|                  | <b>INDICATIONS, CONTRAINDICATIONS, WARNINGS.....16-17</b> |

## PLATE OFFERING:

### HEVANS™ PLATE OFFERINGS

#### Left and Right Side Specific


- ▶ Available in small, medium and large
  - Small: used with a 4 or 6 mm Graft
  - Medium: used with an 8 mm Graft
  - Large: used with a 10 or 12 mm Graft
- ▶ Located in the Universal Plate Caddy
- ▶ Compatible with the Ø2.7 mm, Ø3.5 mm and Ø4.2 mm GORILLA® Screws in Locking or Non-locking Varieties



### BOW & ARROW™ EVANS PLATE OFFERING

#### Universal for Right and Left

- ▶ Available in 4 Sizes: Small (for 6 mm grafts), Medium (8 mm), Large (10 mm), and X-Large (12 mm)
- ▶ Located in the BOW & ARROW™ Plate Caddy
- ▶ Compatible with the Ø2.7 mm, Ø3.5 mm and Ø4.2 mm GORILLA® Screws in Locking or Non-locking Varieties



## PLATE FEATURES:

### HEVANS™ PLATE FEATURES

#### Three Screw Fixation

- ▶ Supports Graft
- ▶ Helps to prevent Graft subsidence with subsequent loss of calcaneal length seen in the early post-operative period due to Graft resorption

#### Single Proximal Screw

- ▶ Helps prevent peroneal tendon irritation



#### Two Distal Screws

- ▶ Helps to avert subluxation of the calcaneocuboid joint with dorsal displacement of the distal calcaneal fragment
- ▶ Dual point fixation at the distal fragment mitigates the prospect of rotation of the distal fragment with the calcaneocuboid joint anatomy

#### Low Profile

- ▶ 1.1 mm plate thickness tapering down to 0.5 mm at the end of the plate to help prevent soft tissue, peroneal tendon, and sural nerve irritation caused by the plate and screws
- ▶ Locking screw technology allows for screw heads that are flush to or below flush to the plate
- ▶ Ramped proximal arm designed specifically to help prevent peroneus brevis tendon irritation such that the tendon can glide over the plate



# System Overview - The Gorilla® Evans Calcaneal Osteotomy System

## PLATE FEATURES:

### BOW & ARROW™ PLATE FEATURES


#### Four Screw Fixation

- ▶ Helps to prevent subluxation of the calcaneocuboid joint with dorsal displacement of the distal calcaneal fragment
- ▶ Helps to prevent graft subsidence with subsequent loss of length of calcaneus seen in the early post-operative period due to graft resorption




#### Tapered “Bow”

- ▶ Mimics the shape of the PRESERVE™ Evans Lateral Column Lengthening Grafts based on plate size
- ▶ Helps evenly distribute vertical pressure along the graft due to the anatomically matched taper
- ▶ Reduced plantar point loading decreases risk of tearing or attenuation of the lateral aspect of the long plantar ligament
- ▶ Arrow on “bow” latches onto the cortex to prevent lateral dislocation of Plate once inserted



#### Low Profile

- ▶ Thickness of 1.15 mm tapering to 0.5 mm
- ▶ Ramped plantar aspect of Plate minimizes risk of soft tissue irritation
- ▶ Locking screw technology allows for screw heads that are flush to or below flush to the plate



## PLATE SCREW COMPATIBILITY:

|                | Ø2.7 mm R3CON Screws:                                                               |                                 | Ø3.5 mm R3CON Screws:                                                               | Ø4.2 mm R3CON Screws:                                                                 |                                     |
|----------------|-------------------------------------------------------------------------------------|---------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-------------------------------------|
| Locking:       |  |                                 |  |  |                                     |
| Non-locking:   |  |                                 |  |  |                                     |
| Screw Lengths: | 8 mm - 20 mm<br>in 1 mm increments                                                  | 22 -40 mm<br>in 2 mm increments | 10 mm - 50 mm<br>in 2 mm increments                                                 | 10 mm - 50 mm<br>in 2 mm increments                                                   | 55 mm - 70 mm<br>in 5 mm increments |

*The surgical technique begins the same way regardless of which type of plate is being used.*

*Following the osteotomy, joint distraction, and placement of graft, **individual surgical techniques will be outlined for each plate, starting on page 8.***


## INCISION/EXPOSURE:



This procedure may be done as a sole procedure or combined with various procedures to address Stage II Posterior Tibial Tendon Dysfunction or a pediatric flatfoot at the discretion of the surgeon. This technique will describe only the Evans Calcaneal Osteotomy (Lateral Column Lengthening) procedure for flatfoot correction.

Patient positioning in a lateral decubitus or supine position with fluoroscopy available is recommended for this procedure. A standard lateral incision is shown but can be varied according to surgeon preference. Dissection is carried down to the lateral wall of the calcaneus without entering the capsule of the calcaneocuboid joint.


## CALCANEAL OSTEOTOMY:



An osteotomy is made in the calcaneus according to surgeon's desired method and technique, generally 1.0 cm - 1.5 cm proximal to and parallel to the calcaneocuboid joint.

# System Overview - The Gorilla® Evans Calcaneal Osteotomy System


## JOINT DISTRACTION:



A K-wire can be placed across the calcaneocuboid joint at this time to prevent dorsal dislocation of the distal fragment of the calcaneus.

There are multiple options available in the Gorilla Plating System for distraction of the calcaneal osteotomy and for visualizing length to achieve desired correction. These options include curved and straight osteotomes, the Paragon 28 Caspar Device or a Hintermann distractor (shown). Other bone spreaders and distraction devices are available upon request.


**TIP:** Although surgeons have their preference on distraction, it is recommended to use a device that is external to the osteotomy such as a Hintermann distractor or the Paragon 28 Caspar Device.





### Hintermann distractor:

- Hold distraction once the appropriate length is achieved and do not prohibit insertion of the wedge into the osteotomy.
- The smaller, inside hole accepts up to a Ø1.6 mm K-wire (more suitable for bending wires out of the way) and the larger, outer hole accepts up to a Ø2.3 mm K-wire.

The Hintermann distractor is secured to either side of the osteotomy using two K-wires.





## JOINT DISTRACTION:



Distract the Hintermann distractor until adequate correction is achieved in desired planes. Confirm reduction of deformity using fluoroscopy, if desired.


## GRAFT INSERTION:



With the osteotomy maintained using a distractor, insert an Evans bone graft properly sized and shaped per surgeon preference, into the osteotomy site.


Once the graft position and size is confirmed using fluoroscopy, remove all instrument distractors.

### HEVANS™ PLATE PLATE PLACEMENT:




The locking Drill Guide is screwed into the dorsal distal hole on the HEvans Plate and can be used as a lever for Plate insertion.

**TIP:** It is advised to use a Drill Guide with both locking and non-locking Screws to avoid Screw head prominence by off-axis drilling due to the thinness of the HEvans™ Plate.




One olive wire can be placed on either side of the osteotomy, if desired. The position of the HEvans Plate should be more dorsal than plantar on the lateral wall of the calcaneus to avoid the peroneal tendons. Confirm the Plate position using fluoroscopy, if desired.

### PERMANENT FIXATION:



Using the Ø2.0 mm Drill, drill a hole through the locking guide.


**TIP:** Locking Screws are recommended for use with the HEvans Plate due to the cancellous nature of bone in this area as well as unicortical fixation. The Plate and Screw construct is designed to provide better protection against loss of correction during Graft resorption and subsequent remodeling.




Remove the locking guide. Measure the Screw length using the depth gauge. Alternatively, the surgeon may read the Screw length from the numerical markings on the Drill off the Drill Guide, if desired.

Insert the selected locking Screw using the provided Screw Driver but do not fully seat the Screw.

## PERMANENT FIXATION:




Remove the proximal olive wire from the HEvans Plate.




Using the technique described above, insert the proximal Screw. It is recommended to complete final tightening of the proximal Screw at this time.

Remove the remaining olive wire. Insert the final distal Screw using the same technique.

Complete final tightening of the two distal Screws using a two finger tightening technique and remove the k-wire crossing the calcaneocuboid joint.



Confirm Plate placement and Screw length using fluoroscopy, if desired.



## CLOSURE:

Proceed to incision closure or concomitant procedures at this point.

### BOW AND ARROW™ PLATE

#### PLATE PLACEMENT:


If used with the BOW & ARROW Plate, the graft is traditionally seated a few millimeters deep to the cortex to allow room for the Plate.



#### PLATE PLACEMENT:




The locking Drill Guide is screwed into one of proximal locking holes on the BOW & ARROW Plate and can be used as a lever for Plate insertion.




It is recommended to insert the Plate from a more dorsal lateral approach, allowing the Plate to be slid down the osteotomy and sit in the position where the “BOW” of the Plate matches the contour of the graft.

## PERMANENT FIXATION:




Because the “arrows” on the “bow” of the Plate latch onto the inside of the cortex to help and prevent lateral displacement of the Plate, temporary fixation using olive wires can be performed but is not necessary.

Confirm the Plate position using fluoroscopy, if desired.




Drill a hole through the locking guide using the Ø2.0 mm Drill. Measure Screw length using a depth gauge or measure off of the Drill using the Drill Guide.




Insert the selected screw size using the Screw Driver but do not completely tighten.

### PERMANENT FIXATION:

---



Insert the final three screws using the technique just described. Tighten all of the screws using a two finger tightening technique, returning to the first screw for tightening.



Remove the K-wire across the calcaneocuboid joint and confirm plate placement and screw length using fluoroscopy, if desired.


### CLOSURE:

---

Proceed to incision closure or concomitant procedures at this time.

## GORILLA® UNIVERSAL PLATE CADDY

The Gorilla® Universal Plate Caddy contains many of the Paragon 28® Gorilla® plates, including the HEvans™ plate discussed in this surgical technique guide.



## GORILLA® BOW & ARROW™ PLATE CADDY

The Gorilla® Plate Caddy contains the BOW & ARROW™ Base Wedge, Cotton, and Evans plates. Additionally, instrumentation is included that is complementary to these plates comprising the tear drop spreader, osteotomes, and the BOW & ARROW™ sizing tool.



## GORILLA® R3CON INSTRUMENT CADDY


Drills, Drill Guides, centering guides, olive wires, Plate benders, Drivers, K-wires and a depth gauge are located in the Gorilla® R3CON Instrument Caddy.



## GORILLA® SCREW CADDY



## GORILLA® CASE



## MINI-MONSTER® SCREW CADDY

The Gorilla® Case can accommodate one Mini-Monster® Screw Caddy if a Ø2.0 mm, Ø2.5 mm, Ø3.0 mm Ø3.5 mm or Ø4.0 mm cannulated Screw is needed during a case.

## GORILLA® R3CON INSTRUMENTS

The Caspar Compression/Distraction device, osteotomes, baby Bennet retractors, bone reduction clamps, periosteal elevator, cartilage removal device, pin distractor and handles are located at the bottom of the Gorilla® Case.

## GORILLA® UNIVERSAL PLATE CADDY:

| Part #       | Description                                      | Quantity |
|--------------|--------------------------------------------------|----------|
| P53-101-1012 | 28 Plate™, 20 mm, 2-Hole Compression             | 1        |
| P53-101-1021 | 28 Plate™, 15 mm, 2-Hole Compression             | 1        |
| P53-101-1022 | 28 Plate™, 20 mm, 2-Hole Compression             | 1        |
| P53-101-1023 | 28 Plate™, 25 mm, 2-Hole Compression             | 1        |
| P53-101-2011 | 28 Plate™, 15 mm, 2-Hole Locking                 | 1        |
| P53-101-2012 | 28 Plate™, 20 mm, 2-Hole Locking                 | 1        |
| P53-101-2013 | 28 Plate™, 25 mm, 2-Hole Locking                 | 1        |
| P53-101-2021 | 28 Plate™, 15 mm, 2-Hole Locking                 | 1        |
| P53-101-2022 | 28 Plate™, 20 mm, 2-Hole Locking                 | 1        |
| P53-101-2023 | 28 Plate™, 25 mm, 2-Hole Locking                 | 1        |
| P53-101-3011 | 28 Plate™, 25 mm, 3-Hole Compression             | 1        |
| P53-101-3012 | 28 Plate™, 30 mm, 3-Hole Compression             | 1        |
| P53-101-3021 | 28 Plate™, 25 mm, 3-Hole Compression             | 1        |
| P53-101-3022 | 28 Plate™, 30 mm, 3-Hole Compression             | 1        |
| P53-101-3023 | 28 Plate™, 35 mm, 3-Hole Compression             | 1        |
| P53-101-4021 | 28 Plate™, 4-Hole Straight, Compression          | 1        |
| P53-101-4121 | 28 Plate™, 4-Hole Straight, Locking              | 1        |
| P53-101-5021 | 28 Plate™, 5-Hole Straight, Compression          | 1        |
| P53-105-L001 | HEvans® Plate, Small, Left, 16 mm                | 1        |
| P53-105-L002 | HEvans® Plate, Medium, Left, 20 mm               | 1        |
| P53-105-L003 | HEvans® Plate, Long, Left, 24 mm                 | 1        |
| P53-105-R001 | HEvans® Plate, Small, Right, 16 mm               | 1        |
| P53-105-R002 | HEvans® Plate, Medium, Right, 20 mm              | 1        |
| P53-105-R003 | HEvans® Plate, Long, Right, 24 mm                | 1        |
| P53-151-0001 | Dogbone™ Plate, 10 mm                            | 1        |
| P53-151-0002 | Dogbone™ Plate, 15 mm                            | 1        |
| P53-151-0003 | Dogbone™ Plate, Compression, 20 mm               | 1        |
| P53-151-0004 | Dogbone™ Plate, Compression, 25 mm               | 1        |
| P53-151-0005 | Dogbone™ Plate, Compression, 30 mm               | 1        |
| P53-151-1002 | Dogbone™ Plate, 15 mm Slanted, Left              | 1        |
| P53-151-1003 | Dogbone™ Plate, 20 mm Slanted, Left              | 1        |
| P53-151-1004 | Dogbone™ Plate, Compression, 25 mm Slanted, Left | 1        |
| P53-152-0000 | T-Plate, Compression, 3-Hole                     | 1        |
| P53-152-0001 | T-Plate, Compression, 4-Hole                     | 1        |
| P53-152-0002 | T-Plate, Compression, 5-Hole                     | 1        |
| P53-153-0001 | Trapezoid Plate, Small (12 mm)                   | 1        |
| P53-153-0002 | Trapezoid Plate, Medium (16 mm)                  | 1        |
| P53-153-0003 | Trapezoid Plate, Long (20 mm)                    | 1        |
| P53-153-0004 | Trapezoid Plate, X-Long (24 mm)                  | 1        |
| P53-154-0001 | Teddy Bear™ Plate, Small                         | 1        |
| P53-154-0002 | Teddy Bear™ Plate, Medium                        | 1        |
| P53-154-0003 | Teddy Bear™ Plate, Large                         | 1        |

## GORILLA® BOW &amp; ARROW PLATE CADDY:

| Part #         | Description                                  | Quantity |
|----------------|----------------------------------------------|----------|
| P50-900-0001   | BOW & Arrow Sizing Tool, Base & Cotton Small | 1        |
| P50-900-BOWC-1 | BOW Implant Caddy, Gorilla                   | 1        |
| P50-900-BOWC-2 | BOW Implant caddy Lid, Gorilla               | 1        |
| P51-920-1001   | Spreader, TearDrop, BOW & Arrow              | 1        |
| P51-920-1002   | BOW & Arrow Sizing Tool, Evans Large         | 1        |
| P51-920-1003   | BOW Osteotome, 6 mm, Dual Taper              | 2        |
| P51-920-1004   | BOW Osteotome, 12 mm, Single Taper           | 1        |
| P53-104-0000   | BOW & Arrow, Base Wedge Plate, 0.0 mm Wedge  | 1        |
| P53-104-0020   | BOW & Arrow, Base Wedge Plate, 2.0 mm Wedge  | 1        |
| P53-104-0025   | BOW & Arrow, Base Wedge Plate, 2.5 mm Wedge  | 1        |
| P53-104-0030   | BOW & Arrow, Base Wedge Plate, 3.0 mm Wedge  | 1        |
| P53-104-0035   | BOW & Arrow, Base Wedge Plate, 3.5 mm Wedge  | 1        |
| P53-104-0040   | BOW & Arrow, Base Wedge Plate, 4.0 mm Wedge  | 1        |
| P53-104-0045   | BOW & Arrow, Base Wedge Plate, 4.5 mm Wedge  | 1        |
| P53-104-0150   | BOW & Arrow Cotton Plate, 5 mm Wedge         | 1        |
| P53-104-0160   | BOW & Arrow Cotton Plate, 6 mm Wedge         | 1        |
| P53-104-0170   | BOW & Arrow Cotton Plate, 7 mm Wedge         | 1        |
| P53-104-0180   | BOW & Arrow Cotton Plate, 8 mm Wedge         | 1        |
| P53-104-0206   | BOW & Arrow, Evans Plate, 6 mm Wedge         | 1        |
| P53-104-0208   | BOW & Arrow, Evans Plate, 8 mm Wedge         | 1        |
| P53-104-0210   | BOW & Arrow, Evans Plate, 10 mm Wedge        | 1        |
| P53-104-0212   | BOW & Arrow, Evans Plate, 12 mm Wedge        | 1        |

The complete contents of the Gorilla Plating System can be found in the full system STG, P51-STG-0001

# MR SAFETY INFORMATION

Refer to [www.paragon28.com/ifus](http://www.paragon28.com/ifus) for the complete and most current instructions for use document.

## INSTRUCTIONS FOR USE: BABY GORILLA®/GORILLA® PLATING SYSTEM

Indications, Contraindications, Warnings and Precautions relevant to the Baby Gorilla®/Gorilla® Plating System are contained in the Instructions for Use document of the Baby Gorilla®/Gorilla® Plating System P51-IFU-3001.

| A patient with the Paragon 28® Baby Gorilla®/Gorilla® Plating System may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient. |                                                                                                                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------|
| Name/Identification of Device                                                                                                                                                                      | Paragon 28® Baby Gorilla®/Gorilla® Plating System                                                                       |
| Nominal value(s) of Static Magnetic Field [T]                                                                                                                                                      | 1.5 T or 3 T                                                                                                            |
| Maximum Spacial Field Gradient [T/m and gauss/cm]                                                                                                                                                  | 30 T/m<br>(3000 gauss/cm)                                                                                               |
| RF Excitation                                                                                                                                                                                      | Circularly Polarized (CP)                                                                                               |
| RF Transmit Coil Type                                                                                                                                                                              | Whole body transmit coil,<br>Head RF transmit-receive coil                                                              |
| Operating Mode                                                                                                                                                                                     | Normal Operating Mode                                                                                                   |
| Maximum Whole Body SAR [W/kg]                                                                                                                                                                      | 2.0 W/kg (Normal Operating Mode)                                                                                        |
| Limits on Scan Duration                                                                                                                                                                            | 2.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks) |
| MR Image Artifact                                                                                                                                                                                  | The presence of this implant may produce an image artifact                                                              |
| If information about a specific parameter is not included, there are no conditions associated with that parameter.                                                                                 |                                                                                                                         |

Refer to [www.paragon28.com/ifus](http://www.paragon28.com/ifus) for the complete and most current instructions for use document.

### INSTRUCTIONS FOR USE: MONSTER®/MINI-MONSTER® SCREW SYSTEM

Indications, Contraindications, Warnings and Precautions relevant to the Monster® Screw System are contained in the Instructions for Use document of the Monster® Screw System P20-IFU-3001.

|                                                                                                                    |                                                                                                                                                                                                                                                                                            |  |
|--------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
|                                   | A patient with the Paragon 28® Monster® Screw System may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.                                                                                                         |  |
| <b>Name/Identification of Device</b>                                                                               | <b>Paragon 28® Monster® Screw System</b>                                                                                                                                                                                                                                                   |  |
| Nominal value(s) of Static Magnetic Field [T]                                                                      | 1.5 T or 3 T                                                                                                                                                                                                                                                                               |  |
| Maximum Spacial Field Gradient [T/m and gauss/cm]                                                                  | 30 T/m<br>(3000 gauss/cm)                                                                                                                                                                                                                                                                  |  |
| RF Excitation                                                                                                      | Circularly Polarized (CP)                                                                                                                                                                                                                                                                  |  |
| RF Transmit Coil Type                                                                                              | Whole body transmit coil,<br>Head RF transmit-receive coil                                                                                                                                                                                                                                 |  |
| Maximum Whole Body SAR [W/kg]                                                                                      | 2.0 W/kg (Normal Operating Mode)                                                                                                                                                                                                                                                           |  |
| Limits on Scan Duration                                                                                            | All anatomical regions can be safely scanned under the following conditions:<br><br>2.0 W/kg whole body average SAR for 5 minutes of continuous RF (a sequence or back to back series/scan without breaks) with a 20 minute cooling period between scans for an hour long scanning session |  |
|                                                                                                                    | Scanning of the knees and all anatomy superior to the knees can be safely scanned under the following conditions:<br><br>2.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)                                           |  |
| MR Image Artifact                                                                                                  | The presence of this implant may produce an image artifact of 20 mm                                                                                                                                                                                                                        |  |
| If information about a specific parameter is not included, there are no conditions associated with that parameter. |                                                                                                                                                                                                                                                                                            |  |



## NOTES:

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



## SURGICAL TECHNIQUE GUIDE

### Gorilla® Evans Calcaneal Osteotomy



P51-STG-3005 RevC [2026-06-25]

© Copyright 2026 Zimmer Biomet®. All rights reserved.  
Patents: [www.paragon28.com/patents](http://www.paragon28.com/patents)

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

Paragon 28, Inc.   
14445 Grasslands Dr.  
Englewood, CO 80112 USA  
(855) 786-2828

Australian Sponsor  
Zimmer Biomet Ltd  
(A.C.N: 096 480 992)(A.B.N: 96 096 480 992)  
Level 4, 12 Narabang Way  
BELROSE NSW 2085  
AUSTRALIA  
Phone: (02) 9483 5400

#### DISCLAIMER

Zimmer Biomet does not practice medicine. This technique was developed in conjunction with healthcare professionals. This document is intended for surgeons and is not intended for laypersons. Each surgeon should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training surgeons have received. As with all surgical procedures, the product/technique used in each case will depend on the surgeon's medical judgment as the best treatment for each patient. Results will vary based on health, weight, activity and other variables. Not all patients are candidates for this product and/or procedure. Caution: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only. ©2025 Zimmer Biomet.

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the instructions for use or contact your local representative; visit [www.paragon28.com](http://www.paragon28.com) for additional product information. This material is intended for health care professionals. Distribution to any other recipient is prohibited.

[www.Paragon28.com](http://www.Paragon28.com)