




3D PRINTED ADDITIVE METAL TECHNOLOGY FOR OPTIMAL POROSITY



HIGHLY CROSS-LINKED VITAMIN E INFUSION TO REDUCE OXIDATION, WEAR DEBRIS, AND POTENTIAL FOR OSTEOLYSIS^{1,2}



MULTI-AXIAL ARTICULATION CONFORMS TO NATURAL TIBIOTALAR MOTION



DESIGN RATIONALE

RESEARCH DRIVEN. | SOLUTION FOCUSED.

Exclusively foot & ankle ²⁰
Paragon®

CONTENTS

DESIGN TEAM | Introduction 3

ANALYSIS | Clinical Investigation 4

SOLUTIONS | Principles of Design - Technology Overview 6



3D PRINTED TIBIAL IMPLANTS



VITAMIN E POLYETHYLENE INSERTS



ANATOMICALLY CONTOURED TALAR DOMES

RESEARCH | Ankle Morphology & Kinematics 11

ARTHRITIC TIBIA & TALAR SIZING AND MORPHOLOGY

Evaluation of Arthritic Ankle Joint Morphology

AN EVALUATION OF HEALTHY ANKLE MOTION AND MORPHOLOGY

Different Phases of Simulated Gait Using a Weight-Bearing CBCT Scan

DISTAL TIBIA MORPHOLOGY

Bi-Radial Curvature Morphology of the Healthy Tibiotalar Joint

TALAR MORPHOLOGY

Bi-Radial Curvature Morphology of the Healthy Talus

TIBIOTALAR JOINT KINEMATICS

In Vivo Kinematics of the Healthy Ankle using Weight-Bearing CT



DESIGN SURGEONS



Mark Dalton, MD
Austin, TX



Mark Myerson, MD
Denver, CO



Jeffery Christensen, DPM
Everett, WA



Michael Houghton, MD
Fort Collins, CO


THE PARAGON 28® APEX 3D™ TOTAL ANKLE REPLACEMENT SYSTEM MARKS THE NEXT GENERATION... in total ankle arthroplasty, with advanced technologies and a keen eye towards optionality. The APEX 3D™ System is the newest addition to the Paragon 28® Precision Ankle Solutions portfolio of products and complements the Gorilla® Ankle Fracture Plating System, Silverback® Fusion Plating System, Phantom™ Hindfoot TTC and ActivCore™ Nail Systems.

BASED ON OVER A CENTURY OF COMBINED CLINICAL EXPERIENCE... and cutting-edge biomechanical research, the APEX 3D™ System was designed to address end-stage ankle arthritis and current challenges within the total ankle market including, but not limited to: implant loosening, pathological wear, instability and persistent pain.

THE P28® DESIGN TEAM HAS A GREAT DEAL OF EXPERIENCE... addressing both historic reasons why ankle replacements fail as well as current insufficiencies such as tibial tray radiolucency and talar dome subsidence. It was important for the Design Team and Paragon 28® to investigate clinically reported modes of ankle replacement failure, conduct research to better understand tibiotalar morphology and joint kinematics in order to introduce clinically relevant solutions.

IMPLANT DESIGNS WERE GUIDED BY...

- Newer clinical evidence
- Advanced technologies
- Implant sizing & anatomic footprint based on pre-clinical arthritic tibiotalar morphology studies
- Healthy ankle kinematic weight-bearing CT research
- Extensive R&D testing



MECHANICAL IMPROVEMENTS WERE DESIGNED TO ADDRESS...

- Inadequate initial fixation
- Mechanical loss of fixation
- Biologic loss of fixation

CLINICAL INVESTIGATION

TOTAL ANKLE REPLACEMENT HAS INCREASINGLY BECOME A VIABLE OPTION... for patients who experience end-stage ankle arthritis. However, limitations still exist and based on literature findings, survivorship rates range from 82.6% to 95.8% at 3 to 5 years³⁻⁷ to 84% - 94% at 5 - 12 years respectively.^{55,56,57,58} (Registry data has shown 73% survivorship at 10 years)

END-STAGE ANKLE ARTHRITIS IS A DISABLING CONDITION... that causes severe pain, substantial functional impairment and decreased quality of life. The most common etiologies found in applicable literature were:

- Post-traumatic arthritis^{3,7-34}
- Primary osteoarthritis^{3,8-10,13,15-18,20,22-24,26-29,31-40}



TAR SURGICAL COMPLICATIONS... in the literature were identified as:

- Component Loosening^{3,6,7,12,19,20,27,30-33,44-46}
- Implant Subsidence^{3,5,6,12,19,20,24,27,28,31-34,39,41}
- Persistent Pain^{3,6,12,14,17,19,20,27,30,32,34,38}
- Intraoperative Malleolar Fracture^{3,5,17,18,30,32,38,46}
- Post-operative Malleolar Fracture^{17,19,25,30,32,34,38,39}

HIGH OCCURRENCE COMPLICATIONS... relative to their study populations included:

- Heterotopic ossification (95.6%)⁴⁷
- Impingement (7.5% to 29%)^{6,7,23,25,30,32,39}

MOST COMMON REASONS FOR REVISIONS...

- Talar subsidence, osteolysis, aseptic loosening^{27,28,34,45,46}

DEVICE REVISION COMPLICATIONS LED TO...

- Patient harm such as pain, reoperations, revision to arthrodesis, or below-the-knee amputations.^{6,25-28,31,33,45,48}


CAUSES OF TAR FAILURES⁵⁴...included:

Cause of failure	Total ankle replacement
Aseptic loosening	38
Luxation Instability	8.5
Septic loosening	9.8
Periprosthetic fracture	2
Pathological wear	8
Persistent pain	12
Implant failure	5.3
Technical error	4.6

#1
#2

Table adapted from Barg et al. Total Ankle Replacement, Deutsches Arzteblatt International (2015) 112,177-84.

POTENTIAL LIMITATIONS OF JOINT REPLACEMENT... identified in the literature included:




BONE RESORPTION... is currently considered a natural consequence of total joint arthroplasty due to changes in bone remodeling patterns. Bone remodeling is mediated by the changes in stress distribution caused by implantation. Extensive resorption around the prosthesis may lead to implant loosening and failure.

OSTEOLYSIS... is the result of localized foreign-body reaction to particulate debris generated by cement, metal, and Ultra-High Molecular Weight Polyethylene (UHMWPE). Regarding the etiology, it has been hypothesized that particulate debris generated by the components of a prosthesis migrate into the synovial cavity and the bone-implant interface, where they recruit macrophages and stimulate phagocytic action. The phagocytic action results in the release of cytokines and intercellular mediators (IL-1, 2, PE2) which encourage osteoclastic bone resorption. Clinical and basic research is continuing in order to provide scientific basis for the causes of this phenomenon and the potential ways to reduce its occurrence.

AS A RESULT OF THIS CLINICAL INVESTIGATION, PARAGON 28® SOUGHT TO INTRODUCE... implant design features, precision bone preparation instrumentation and advanced technologies to help mitigate the occurrence of surgical complications.

PRINCIPLES OF DESIGN

PARAGON 28® APEX 3D™ TOTAL ANKLE SYSTEM... has been efficiently streamlined to adapt to surgeon preference, accommodate patient needs and is equipped to address a wide range of arthritic ankle conditions. Implant designs are based on pre-clinical arthritic tibiotalar morphology studies, guided by healthy ankle weight-bearing CT research and extensive clinical investigation findings. The APEX 3D™ System consists of cemented, fixed-bearing anatomically contoured implant components, intuitive instrumentation and precision bone resection guides intended for use in primary or revision surgery for patients with ankle joints damaged by severe rheumatoid, post-traumatic, or degenerative arthritis. Revision surgery is also indicated when sufficient bone stock is present. Components are intended for cemented use only*.




 **APEX 3D™**

*For additional information regarding Indications for Use, Contraindications, Warnings, Precautions, etc. please visit:


<https://www.paragon28.com/ifus/>

PRINCIPLES OF DESIGN

3D PRINTED BONE FACING GEOMETRY



ARC Tibia™ Implant



FLAT Tibia Implant

APEX 3D™ PRINTED TIBIAL COMPONENTS...

- 3D Printed (*Titanium Alloy*)
 - Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial)
- 52 Anatomically contoured low-profile sizing options available in sizes 1 - 6 (Lengths: Short*, Standard & Long)
- Left & right specific components offered in:
 - 2 distinct ARC Tibia™ & FLAT Tibia geometric configurations
- Anterior/posterior lengths, curvatures and medial/lateral widths based on arthritic ankle pre-clinical morphology studies

DESIGNED TO PROVIDE...



- Optimal porosity and rotational stability
- Broad anterior contact surface interface
- Distal tibia bone coverage with off-set posterior curvature to reduce potential for component loosening
- Medial/lateral side-wall surface area for added fixation and rotational stability
- Initial fixation and stabilization with two non-coated vertical pegs positioned slightly posterior to midline where peak bone density is located ⁴⁹

TIBIAL TRAY IMPLANT DESIGN BASED ON... pre-clinical arthritic tibiotalar morphology studies and healthy ankle weight-bearing CT research. 3D printed additive metal manufacturing process was introduced to allow for a low-profile, singular construct with integrated gradient zones of porosity.

Anatomically contoured implants include: 3D printed low-profile style options (ARC Tibia™ & Flat Tibia, 5 mm in thickness) with porous architecture and gradient zones down to solid titanium substrate. The ARC Tibia™ geometric architecture was designed to provide rotational stability and decrease rotational forces reducing potential for component loosening as seen in existing tibial implants. ^{19,20} The basis of the two posterior-biased vertical peg design was introduced to mitigate findings from clinical publications citing posterior lucency and subsidence of existing tibial components. ^{19,20,39} Vertical pegs were also designed to reduce posterior translation during implantation and mitigate potential for posterior gapping compared to designs with anterior peg placement.

Enhanced bone preparation techniques were introduced, featuring corners created from a drill rather than two intersecting saw cuts, resulting in rounded edges intended to help mitigate the risk of tibial stress fracture ^{38,39} and malleolar fracture seen in existing tibial implant designs, ^{19,39} as well as post-operative malleolar fracture seen in other total ankle replacement designs. ^{17,25,30,32,34}

Cross Section of Porous Titanium Architecture



PRINCIPLES OF DESIGN

HIGHLY CROSS-LINKED VITAMIN E INFUSION



Neutral Poly Insert

VITAMIN E TIBIAL INSERT - HIGHLY CROSS-LINKED POLYETHYLENE

- Ultra-high-molecular-weight polyethylene (UHMWPE) + Medical Grade Vitamin E Formula
- 35 Sizing options with universal locking mechanism available in sizes 1 - 5
- 7 primary and revision thicknesses options including 6 - 12mm offered in 1mm increments
- 5 semi-constrained sizing options with bicondylar bearing surface that allow for anatomic conformity

IMPLANT OPTIONALITY & WEAR CHARACTERISTICS...

- Vitamin E infusion introduced to help reduce oxidation, long-term wear and potential for osteolysis that can lead to implant failure.
- Equipped with a universal locking mechanism for additional implant optionality
 - Provides the ability to: upsize by one & downsize without restriction*

*Largest size tibia component, size 6 can be paired with smallest size 1 poly insert.

*Polyethylene size will always follow the selected talar dome

RECENT RESEARCH ON... *In vivo* oxidation has shown that UHMWPE may oxidize, through not only radiation induced free radicals, but also free radicals generated in vivo by cyclic loading and by the reaction of lipids absorbed from the synovial fluid.¹ Vitamin E infusion was introduced to help reduce oxidation, wear debris and potential for osteolysis that can lead to implant failure.^{1,2}

Researchers have suggested that it may be useful to have sufficient Vitamin E activity in UHMWPE to prevent oxidation by thermal degradation during radiation cross-linking and compression molding. In a series of experiments, researchers from Kyoto University demonstrated that UHMWPE sterilized components blended with Vitamin E exhibited resistance to oxidation and fatigue wear, when compared with unstabilized controls.¹

CONTACT AREA ANALYSIS... A 3D CAD model finite element analysis (FEA) was conducted to evaluate and determine contact area.^[TR-19021119]

ANT

Medial



Dorsiflexion Contact Area



Neutral Flexion Contact Area



Plantar Flexion Contact Area

ANT

Lateral


POST

POST


Reflects right-sided limb and joint articulation

PRINCIPLES OF DESIGN

ANATOMICALLY CONTOURED TALAR DOME



Chamfer Talar Implant



Flat-Cut Talar Implant

ANATOMICALLY CONTOURED TALAR DOMES...

- Cobalt-chromium-molybdenum alloy (CoCr)
- Titanium Plasma Spray undercoating (TPS)
- 12 Anatomically contoured sizing options available in sizes 1 narrow (1N) & 1 - 5.
- Left & right specific components offered in two distinct style configurations: Chamfer & Flat
- Anterior/posterior lengths, curvatures, and medial/lateral widths based on - Tibia & Talar Sizing & Morphology Study- Sizes 1N - 5

DESIGNED TO PROVIDE...

- Anatomic fit and optimize talar bone coverage
- Initial fixation and rotational stability
- Intercondylar spacing for coronal plane stability to resist varus/valgus motion
- Multi-axial smooth articulation that mimics natural tibiotalar motion reducing the potential for bony and component impingement
- Tapered medial/lateral shoulders to reduce potential for painful gutter impingement

ANATOMIC RELEVANT LOADING... in healthy patients is critical to understanding interactions of joint articulations. The APEX 3D talar dome implant design considers tibiotalar motion during simulated gait using modern weight-bearing CT technology. The design and enhanced bone preparation techniques allow the implants to be properly impacted, overcoming issues of posterior gapping and potential for component loosening, as noted in current designs.^{12,19,20,31} The chamfer option includes an anterior biased non-coated central fin for ease of insertion and initial stability. The Flat-cut option is equipped with dual non-coated anterior pegs designed for rotational stability, and both options provide posterior medial coverage to help to avoid FHL impingement.


WEIGHT-BEARING CT BASED TALAR MORPHOLOGY RESEARCH⁵⁹

Paragon 28® sought to better understand if the curvature of the medial and lateral aspects of the talus could be more accurately described by dividing the condyles into anterior and posterior regions, creating bi-radial curves.

Paragon 28® found that the healthy ankle joint experiences varying, tri-axial rotation demonstrating that Inman,⁵⁰ Seigler,⁵¹ Barnett and Napier,⁵² were all correct in part with their assessments of ankle joint line axes.


FINDINGS: *Lateral Anterior Radius > Medial Anterior Radius. Medial Posterior Radius > Lateral Posterior Radius*⁵⁹

BIOMECHANICAL IMPLICATIONS: Confirms the ankle does not have a single joint line axis.



IMPLANT SIZING RESEARCH


RESEARCH OVERVIEW Implants were developed based on arthritic ankle morphology scans then compared for fit and coverage to a second subset of subjects to confirm coverage rate. CT scans of arthritic patients for total ankle replacement were segmented and landmarks identified on each bone in the ankle joint: tibia, fibula, and talus. Measurements were taken to determine size and morphology data for the subject population. Statistical analyses were conducted on the measurement data to identify means, outliers, correlations, and measurements that best corresponded with input dimensions for the development of a new total ankle replacement (TAR) implant system. Once identified, the sizing dimension data were used to drive the critical dimensions for a system of TAR tibial trays and talar dome components based on a desired coverage rate of 95%. Implants were developed based on these inputs and compared for fit and coverage to a second subset of subjects to confirm coverage rate. ^[TR-19021117]



THE PURPOSE OF THIS STUDY... was to depict the most appropriate fitting ankle prosthesis on to a range of sizes of subject bones to ensure the design input requirements of this implant system were fulfilled.


RESULTS... at the resection levels, tibia anterior/posterior lengths were taken at the medial, lateral and center to establish optimal bone coverage and size offerings. Similarly, the talar medial-lateral widths at the anterior & posterior aspects were taken.

MORPHOLOGY & KINEMATICS



MEDIAL RADII

LATERAL RADII



BI-RADIAL CURVATURE MORPHOLOGY OF THE HEALTHY TIBIOTALAR JOINT ⁵⁹

Study Goals:

- More accurately describe the curvature of the distal tibia as well as the medial and lateral talar condyles
- Create bi-radial curves for both the medial and lateral sides of the distal tibia and proximal talus
- Compare the regional congruency between the distal tibia and proximal talus

Distal Tibia Findings:

- Anterior Tibial Radii > Anterior Talar Radii
- Lateral Posterior Tibial Radii < Lateral Posterior Talar

Proximal Talus Findings:

- Radius increases from anterior to posterior of both the medial and lateral sides of the talus
- Lateral Anterior Radius > Medial Anterior Radius
- Lateral Posterior Radius < Medial Posterior Radius


IN VIVO KINEMATICS OF THE HEALTHY ANKLE USING WEIGHT-BEARING CT ^[TR-19021110]

Study Goals:

- Measure talar rotation with respect to tibia through three different stages of mid-stance (*early mid-stance, mid-stance, late mid-stance*).

Findings:

- The healthy ankle joint experiences varying, tri-axial rotation



Early Mid-Stance

Mid-Stance

Late Mid-Stance

- 1 Kurtz, S.M., Bracco, P., Costa, L., Oral, E., Muratoglu, O.K. (2016) Vitamin E-Blended UHMWPE Biomaterials. UHMWPE Biomaterials Handbook, 293-325
- 2 G. Rochongar, MD, G. Buia, MD, E. Bourroux, J. Dunet, MD, V. Chapus, MD, and C. Hulet, MD, PhD. (2018) Creep and Wear in Vitamin E-Infused Highly Cross-Linked Polyethylene Cups for Total Hip Arthroplasty A Prospective Randomized Controlled Trial. *Journal of Bone and Joint Surgery, Incorporated*.
- 3 Adams, S. B., Jr., Demetracopoulos, C. A., Queen, R. M., Easley, M. E., DeOrio, J. K., & Nunley, J. A. (2014). Early to mid-term results of fixed-bearing total ankle arthroplasty with a modular intramedullary tibial component. *J Bone Joint Surg Am*, 96(23), 1983-1989. doi:10.2106/jbjs.m.01386
- 4 Alvine, G. F., & Alvine, F. G. (2016). Total Ankle Arthroplasty Using the Agility Stemmed Talar Revisional Component: Three to Eight Year Follow-Up. *S D Med*, 69(4), 151-154.
- 5 Chao, J., Choi, J. H., Grear, B. J., Tenenbaum, S., Bariteau, J. T., & Brodsky, J. W. (2015). Early radiographic and clinical results of Salto total ankle arthroplasty as a fixed-bearing device. *Foot Ankle Surg*, 21(2), 91-96. doi:10.1016/j.fas.2014.09.012
- 6 Harston, A., Lazarides, A. L., Adams, S. B., Jr., DeOrio, J. K., Easley, M. E., & Nunley, J. A., 2nd. (2017). Midterm Outcomes of a Fixed-Bearing Total Ankle Arthroplasty With Deformity Analysis. *Foot Ankle Int*, 38(12), 1295-1300. doi:10.1177/1071100717731853
- 7 Stewart, M. G., Green, C. L., Adams, S. B., Jr., DeOrio, J. K., Easley, M. E., & Nunley, J. A., 2nd. (2017). Midterm Results of the Salto Talaris Total Ankle Arthroplasty. *Foot Ankle Int*, 38(11), 1215-1221. doi:10.1177/1071100717719756
- 8 Ajis, A., Henriquez, H., & Myerson, M. (2013). Postoperative range of motion trends following total ankle arthroplasty. *Foot Ankle Int*, 34(5), 645-656. doi:10.1177/1071100713481433
- 9 Angthong, C. (2015). Health-related quality of life and functional outcomes in ankle arthritis patients based on treating with and without total ankle replacement surgery. *J Med Assoc Thai*, 98 Suppl 2, S131-137.
- 10 Angthong, C., Adams, S. B., Easley, M. E., DeOrio, J. K., & Nunley, J. A., 2nd. (2014). Radiation exposure in total ankle replacement. *Foot Ankle Int*, 35(11), 1131-1136. doi:10.1177/1071100714548062
- 11 Benich, M. R., Ledoux, W. R., Orendurff, M. S., Shofer, J. B., Hansen, S. T., Davitt, J., . . . Sangeorzan, B. J. (2017). Comparison of Treatment Outcomes of Arthrodesis and Two Generations of Ankle Replacement Implants. *J Bone Joint Surg Am*, 99(21), 1792-1800. doi:10.2106/jbjs.16.01471
- 12 Berkowitz, M. J., Clare, M. P., Walling, A. K., & Sanders, R. (2011). Salvage of failed total ankle arthroplasty with fusion using structural allograft and internal fixation. *Foot Ankle Int*, 32(5), S493-502. doi:10.3113/fai.2011.0493
- 13 Bleazey, S. T., Brigido, S. A., & Protzman, N. M. (2013). Perioperative complications of a modular stem fixed-bearing total ankle replacement with intramedullary guidance. *J Foot Ankle Surg*, 52(1), 36-41. doi:10.1053/jjfas.2012.10.009
- 14 Bonnin, M., Gaudot, F., Laurent, J. R., Ellis, S., Colombier, J. A., & Judet, T. (2011). The Salto total ankle arthroplasty: survivorship and analysis of failures at 7 to 11 years. *Clin Orthop Relat Res*, 469(1), 225-236. doi:10.1007/s11999-010-1453-y
- 15 Braito, M., Dammerer, D., Kaufmann, G., Fischler, S., Carollo, J., Reinthaler, A., . . . Biedermann, R. (2014). Are our expectations bigger than the results we achieve? a comparative study analysing potential advantages of ankle arthroplasty over arthrodesis. *Int Orthop*, 38(8), 1647-1653. doi:10.1007/s00264-014-2428-5
- 16 Brigido, S. A., Wobst, G. M., Galli, M. M., Bleazey, S. T., & Protzman, N. M. (2015). Evaluating component migration after modular stem fixed-bearing total ankle replacement. *J Foot Ankle Surg*, 54(3), 326-331. doi:10.1053/jjfas.2014.06.028
- 17 Claridge, R. J., & Sagherian, B. H. (2009). Intermediate term outcome of the agility total ankle arthroplasty. *Foot Ankle Int*, 30(9), 824-835. doi:10.3113/fai.2009.0824
- 18 Clement, R. C., Krynetskiy, E., & Parekh, S. G. (2013). The total ankle arthroplasty learning curve with third-generation implants: a single surgeon's experience. *Foot Ankle Spec*, 6(4), 263-270. doi:10.1177/1938640013493463
- 19 Cody, E. A., Taylor, M. A., Nunley, J. A., 2nd, Parekh, S. G., & DeOrio, J. K. (2019). Increased Early Revision Rate With the INFINITY Total Ankle Prosthesis. *Foot Ankle Int*, 40(1), 9-17. doi:10.1177/1071100718794933
- 20 Penner, M., Davis, W. H., Wing, K., Bemenderfer, T., Waly, F., & Anderson, R. B. (2018). The Infinity Total Ankle System: Early Clinical Results With 2- to 4-Year Follow-up. *Foot Ankle Spec*, 1938640018777601. doi:10.1177/1938640018777601
- 21 Criswell, B., Hunt, K., Kim, T., Chou, L., & Haskell, A. (2016). Association of Short-term Complications With Procedures Through Separate Incisions During Total Ankle Replacement. *Foot Ankle Int*, 37(10), 1060-1064. doi:10.1177/1071100716651964
- 22 Dekker, T. J., Hamid, K. S., Federer, A. E., Steele, J. R., Easley, M. E., Nunley, J. A., & Adams, S. B., Jr. (2018). The Value of Motion: Patient-Reported Outcome Measures Are Correlated With Range of Motion in Total Ankle Replacement. *Foot Ankle Spec*, 11(5), 451-456. doi:10.1177/1938640017750258
- 23 Gross, C. E., Adams, S. B., Easley, M., Nunley, J. A., 2nd, & DeOrio, J. K. (2017). Surgical Treatment of Bony and Soft-Tissue Impingement in Total Ankle Arthroplasty. *Foot Ankle Spec*, 10(1), 37-42. doi:10.1177/1938640016666918
- 24 Taylor, M. A., Green, C. L., Risoli, T. J., DeOrio, J. K., Easley, M. E., Nunley, J. A., & Adams, S. B. (2018). Hindfoot Arthrodesis Screw Position and Trajectory Effect on Talus Subsidence When Performed With Total Ankle Arthroplasty. *Foot Ankle Int*, 1071100718806479. doi:10.1177/1071100718806479
- 25 Criswell, B. J., Douglas, K., Naik, R., & Thomson, A. B. (2012). High revision and reoperation rates using the Agility Total Ankle System. *Clin Orthop Relat Res*, 470(7), 1980-1986. doi:10.1007/s11999-012-2242-6
- 26 Demetracopoulos, C. A., Adams, S. B., Jr., Queen, R. M., DeOrio, J. K., Nunley, J. A., 2nd, & Easley, M. E. (2015). Effect of Age on Outcomes in Total Ankle Arthroplasty. *Foot Ankle Int*, 36(8), 871-880. doi:10.1177/1071100715579717
- 27 DeVries, J. G., Scott, R. T., Berlet, G. C., Hyer, C. F., Lee, T. H., & DeOrio, J. K. (2013). Agility to INBONE: anterior and posterior approaches to the difficult revision total ankle replacement. *Clin Podiatr Med Surg*, 30(1), 81-96. doi:10.1016/j.cpm.2012.08.011
- 28 Ellington, J. K., Gupta, S., & Myerson, M. S. (2013). Management of failures of total ankle replacement with the agility total ankle arthroplasty. *J Bone Joint Surg Am*, 95(23), 2112-2118. doi:10.2106/jbjs.k.00920
- 29 Jiang, J. J., Schipper, O. N., Whyte, N., Koh, J. L., & Toolan, B. C. (2015). Comparison of perioperative complications and hospitalization outcomes after ankle arthrodesis versus total ankle arthroplasty from 2002 to 2011. *Foot Ankle Int*, 36(4), 360-368. doi:10.1177/1071100714558511
- 30 Krause, F. G., Windolf, M., Bora, B., Penner, M. J., Wing, K. J., & Younger, A. S. (2011). Impact of complications in total ankle replacement and ankle arthrodesis analyzed with a validated outcome measurement. *J Bone Joint Surg Am*, 93(9), 830-839. doi:10.2106/jbjs.j.00103

- 31 Lachman, J. R., Ramos, J. A., Adams, S. B., Nunley, J. A., 2nd, Easley, M. E., & DeOrio, J. K. (2019). Patient-Reported Outcomes Before and After Primary and Revision Total Ankle Arthroplasty. *Foot Ankle Int*, 40(1), 34-41. doi:10.1177/1071100718794956
- 32 Pangrazzi, G. J., Baker, E. A., Shaheen, P. J., Okeagu, C. N., & Fortin, P. T. (2018). Single-Surgeon Experience and Complications of a Fixed-Bearing Total Ankle Arthroplasty. *Foot Ankle Int*, 39(1), 46-58. doi:10.1177/1071100717735288
- 33 Raikin, S. M., Sandrowski, K., Kane, J. M., Beck, D., & Winters, B. S. (2017). Midterm Outcome of the Agility Total Ankle Arthroplasty. *Foot Ankle Int*, 38(6), 662-670. doi:10.1177/1071100717701232
- 34 Williams, J. R., Wegner, N. J., Sangeorzan, B. J., & Brage, M. E. (2015). Intraoperative and perioperative complications during revision arthroplasty for salvage of a failed total ankle arthroplasty. *Foot Ankle Int*, 36(2), 135-142. doi:10.1177/1071100714554452
- 35 Coe, M. P., Sutherland, J. M., Penner, M. J., Younger, A., & Wing, K. J. (2015). Minimal clinically important difference and the effect of clinical variables on the ankle osteoarthritis scale in surgically treated end-stage ankle arthritis. *J Bone Joint Surg Am*, 97(10), 818-823. doi:10.2106/jbjs.n.00147
- 36 Coetzee, J. C., Petersen, D., & Stone, R. M. (2017). Comparison of Three Total Ankle Replacement Systems Done at a Single Facility. *Foot Ankle Spec*, 10(1), 20-25. doi:10.1177/1938640016666926
- 37 Hsu, A. R., Davis, W. H., Cohen, B. E., Jones, C. P., Ellington, J. K., & Anderson, R. B. (2015). Radiographic Outcomes of Preoperative CT Scan-Derived Patient-Specific Total Ankle Arthroplasty. *Foot Ankle Int*, 36(10), 1163-1169. doi:10.1177/1071100715585561
- 38 King, A., Bali, N., Kassam, A. A., Hughes, A., Talbot, N., & Sharpe, I. (2018). Early outcomes and radiographic alignment of the Infinity total ankle replacement with a minimum of two year follow-up data. *Foot Ankle Surg*. doi:10.1016/j.fas.2018.11.007
- 39 Saito, G. H., Sanders, A. E., de Cesar Netto, C., O'Malley, M. J., Ellis, S. J., & Demetracopoulos, C. A. (2018). Short-Term Complications, Reoperations, and Radiographic Outcomes of a New Fixed-Bearing Total Ankle Arthroplasty. *Foot Ankle Int*, 39(7), 787-794. doi:10.1177/1071100718764107
- 40 Stavrakis, A. I., & SooHoo, N. F. (2016). Trends in Complication Rates Following Ankle Arthrodesis and Total Ankle Replacement. *J Bone Joint Surg Am*, 98(17), 1453-1458. doi:10.2106/jbjs.15.01341
- 41 Bonnin, M. P., Laurent, J. R., & Casillas, M. (2009). Ankle function and sports activity after total ankle arthroplasty. *Foot Ankle Int*, 30(10), 933-944. doi:10.3113/fai.2009.0933
- 42 Lefrancois, T., Younger, A., Wing, K., Penner, M. J., Dryden, P., Wong, H., ... Glazebrook, M. (2017). A Prospective Study of Four Total Ankle Arthroplasty Implants by Non-Designer Investigators. *J Bone Joint Surg Am*, 99(4), 342-348. doi:10.2106/jbjs.16.00097
- 43 Merrill, R. K., Ferrandino, R. M., Hoffman, R., Ndu, A., & Shaffer, G. W. (2018). Comparing 30-day all-cause readmission rates between tibiotalar fusion and total ankle replacement. *Foot Ankle Surg*. doi:10.1016/j.fas.2017.12.011
- 44 Brigido, S. A., Mulhern, J. L., Wobst, G. M., & Protzman, N. M. (2015). Preoperative and Postoperative Range of Motion: A Retrospective Comparison of Two Total Ankle Replacement Systems. *J Foot Ankle Surg*, 54(5), 809-814. doi:10.1053/jjfas.2014.12.036
- 45 Myerson, M. S., Shariff, R., & Zonno, A. J. (2014). The management of infection following total ankle replacement: demographics and treatment. *Foot Ankle Int*, 35(9), 855-862. doi:10.1177/1071100714543643
- 46 Roukis, T. S., Iceman, K., & Elliott, A. D. (2016). Intraoperative Radiation Exposure During Revision Total Ankle Replacement. *J Foot Ankle Surg*, 55(4), 732-737. doi:10.1053/jjfas.2016.01.039
- 47 Anghong, C., Adams, S. B., Easley, M. E., & Nunley, J. A., 2nd. (2014). Heterotopic ossification following total ankle replacement: clinical significance and factors affecting its formation. *J Surg Orthop Adv*, 23(4), 203-213.
- 48 Patton, D., Kiewiet, N., & Brage, M. (2015). Infected total ankle arthroplasty: risk factors and treatment options. *Foot Ankle Int*, 36(6), 626-634. doi:10.1177/1071100714568869
- 49 Hvid, I. et al. (1985) Trabecular Bone Strength Profiles at the Ankle Joint. *Clinical Orthopaedics and Related Research*, 306-312.
- 50 Inman VT. *The Joints of the Ankle*. Williams & Wilkins, Baltimore, 1976
- 51 Siegler S, et al. New Observations on the morphology of the talar dome and its relationship to ankle kinematics. *Clinical Biomechanics* October 2013
- 52 Barnett CH, Napier JR. The Axis of Rothation at the Ankle Joint in Man: its Influence upon the form of the Talus and the Mobility of the Fibula. *J Anat* 1952;86(91):1-9
- 53 Myerson M, et al. Bi-Radial Curvature Morphology of the Healthy Talus. In: Chicago, IL: American Orthopaedic Foot & Ankle Society Annual Meeting (e-poster); 2019.
- 54 Barg, A., Wimmer, M. D., Wiewiorski, M., Wirtz, D. C., Pagenstert, G. I., & Valderrabano, V. (2015). Total ankle replacement. *Dtsch Arztebl Int*, 112(11), 177-184. doi:10.3238/arztebl.2015.0177
- 55 Zaidi R, Cro S, Gurusamy K, Siva N, Macgregor A, Henricson A, Goldberg A. The outcome of total ankle replacement: a systematic review and meta-analysis. *Bone Joint J*. 2013 Nov;95-B(11):1500-7. doi:10.1302/0301-620X.95B11.31633. PMID: 24151270.
- 56 Barg, Alexej MD; Zwicky, Lukas MSc; Knupp, Markus MD; Henninger, Heath B. PhD; Hintermann, Beat MD HINTEGRA Total Ankle Replacement: Survivorship Analysis in 684 Patients, *The Journal of Bone & Joint Surgery: July 3, 2013 - Volume 95 - Issue 13 - p 1175-1183* doi:10.2106/JBJS.L.01234
- 57 Daniels TR, Mayich DJ, Penner MJ. Intermediate to Long-Term Outcomes of Total Ankle Replacement with the Scandinavian Total Ankle Replacement (STAR). *J Bone Joint Surg Am*. 2015 Jun 3;97(11):895-903. doi: 10.2106/JBJS.N.01077. PMID: 26041850.
- 58 Palanca A, Mann RA, Mann JA, Haskell A. Scandinavian Total Ankle Replacement: 15-Year Follow-up. *Foot & Ankle International*. 2018;39(2):135-142. doi:10.1177/1071100717738747
- 59 Anderle MR, Obert RM, Paxson RD, Brinker LZ, Clancy JT. A bi-radial approach to define the sagittal geometry of the healthy ankle. *Foot Ankle Surg*. 2020 Nov 4;S1268-7731(20)30226-5. doi:10.1016/j.fas.2020.10.011. Epub ahead of print. PMID: 33214032.

A series of horizontal dashed lines providing a template for writing notes.



P10-SDR-0001 Rev. A

™Trademarks and ®Registered Marks of Paragon 28®, Inc.
© Copyright 2021 Paragon 28®, Inc. All rights reserved.
Patents: www.paragon28.com/patents

Paragon 28, Inc. 
14445 Grasslands Dr.
Englewood, CO 80112
(855) 786-2828

PATENTED, DESIGNED & EXCLUSIVELY
DISTRIBUTED BY

Exclusively foot & ankle **28**
Paragon®

www.Paragon28.com