



V92-FC+
PRODUCT INFORMATION

Viable Allogeneic Bone Scaffold



HARNESSING NATURE & TECHNOLOGY WITH EASE

YOU ASKED FOR AN **EASY-TO-USE**, **STREAMLINED PREPARATION PROCESS**, AND WE HEARD YOU!

V92-FC+ allograft is the next generation solution for bone formation to support a variety of potential clinical applications. The allograft is delivered in an easy-to-use syringe with minimal preparation time of under 15 minutes. Paragon 28[®] preserves the native bone cells in a DMSO-free (free of dimethyl sulfoxide) cryoprotectant, requiring no rinsing or decanting — just thaw and use!

V92-FC+ PROVIDES THE THREE KEY ELEMENTS IDEAL FOR BONE FORMATION

- **An osteoconductive three-dimensional scaffold with cortical and cancellous components**
- **A demineralized cortical bone scaffold**
 - ▶ Demineralized cortical bone has been identified to have osteoinductive potential¹
- **Viable endogenous bone cells to support osteogenic healing processes**

BONE SCAFFOLD DELIVERS OSTEOCONDUCTIVE AND OSTEOINDUCTIVE POTENTIAL

V92-FC+ provides an osteoconductive bone scaffold composed of mineralized cancellous bone along with demineralized cortical fibers. Bone fibers offer superior osteoconductivity when compared to powder due to the increased ability for cells to migrate along fibers, creating “cellular highways” for bone formation.³ In contrast, particulate-based demineralized bone matrices (DBMs) have gaps between the particles that osteoblasts cannot always bridge across.³ The demineralized cortical fibers are supplemented with cancellous chips to deliver a 100% human-derived product that mimics the particulate structure of native bone.

OPERATING ROOM EASE OF USE

- **Packaged in easy-to-use syringe**
- **No rinsing or decanting steps required — native bone cells are preserved in a DMSO-free cryoprotectant**
- **Average cell viability exceeds 92% post-thaw²**
- **Average of 1.5 million viable cells per cc of allograft²**
- **Four-hour working window for implantation after thaw without loss of cell viability**



CELLS PROTECTED BY PROPRIETARY CRYOPROTECTANT

- **Protective coating preserves allograft and prevents crack propagation and membrane lysis²**
- **Retains over 92% cell viability after thaw²**
- **Non-cytotoxic, non-DMSO**

- ▶ Reduces concerns about cytotoxicity and negative effects on cell differentiation^{4,5,6}
- ▶ Does not require rinsing or decanting




FIGURE 1

Cells protected with proprietary cryoprotectant to prevent crystalline damage (frozen)
*image captured by SEM




FIGURE 2

Viable endogenous cells are preserved with the use of proprietary cryoprotectant (frozen)




FIGURE 3

2.5% DMSO-compromised cells showing reduced viability (thawed)

Proper preservation of cellular allografts requires strict adherence to recovery and processing protocols. To manufacture V92-FC+, viable endogenous bone cells are collected from the donor and preserved with the use of a next-generation DMSO-free cryoprotectant, which uses an extracellular protective coating on the cell to prevent crack propagation and membrane lysis² (Figure 1). Industry standard DMSO penetrates the cell and prevents crystal formation from within (Figure 3). At room temperature, DMSO-based cryoprotectants raise concerns about cytotoxicity and negative effects on cell differentiation.^{4,5,6}


Paragon 28[®]'s patented and proprietary cryoprotectant is a differentiated technology that is applied to a number of products in our portfolio. The technology provides our products with distinct advantages over DMSO-based cryoprotectant technology used in competitive products. As an added bonus, because our cryoprotectant is DMSO-free, it does not require the multiple rinsing and decanting steps of DMSO-based cryoprotectants. Rinsing steps can diminish both the cell viability and the inherent regenerative properties of allografts.

Our cryoprotectant provides a surgical procedure advantage over other cryoprotectants containing DMSO. Allografts treated with our cryoprotectant experience minimal cell loss and retain, on average, over 80% cell viability after thaw² (Figure 2). Our cryoprotectant also allows for usage up to four hours after thawing and V92-FC+ allografts can be stored for up to one year at or below -65°C.

The bone cells are endogenous to the cancellous bone, remaining attached throughout the donor tissue processing event. Strict donor criteria and quality control processes verify a viable cell population for osteogenic supplementation as a viable structural allograft.

IT ALL ADDS UP

- **Improved storage container streamlines preparation: thaw product in provided syringe and use**
- **A natural, 100% tissue scaffold of demineralized cortical bone fibers coupled with chips rich with endogenous bone cells provides an optimal microenvironment for osteogenesis and excellent handling**
- **A proprietary DMSO-free cryoprotectant protects and allows for consistent delivery of viable allograft to the patient**
- **A viable cell population for osteogenic supplementation as a viable structural allograft**



ORDERING INFORMATION

PART NUMBER	DESCRIPTION
P01-V92-0102	V92-FC+, BONE MATRIX, 1.0CC
P01-V92-0252	V92-FC+, BONE MATRIX, 2.5CC
P01-V92-0502	V92-FC+, BONE MATRIX, 5.0CC
P01-V92-1002	V92-FC+, BONE MATRIX, 10.0CC



V92FC+-BRO-01 Rev B
2022-09-12

™Trademarks and ®Registered Marks of Paragon 28®, Inc.
© Copyright 2022 Paragon 28®, Inc. All rights reserved.
Patents: www.paragon28.com/patents

Distributed by:
Paragon 28, Inc.
14445 Grasslands Dr.
Englewood, CO 80112 USA
(855) 786-2828

Paragon 28 Medical Devices Trading Limited
First Floor Block 7 Beckett Way
Park West Business Park
Dublin 12, D12 X884, Ireland
+353 (0) 1588 0350

For the contraindications, potential complications and adverse reactions, warnings and precautions associated with this device, please refer to the device specific instructions for use at <http://www.paragon28.com/ifus>

Processed By:
VIVEX Biologics, Inc.
2430 NW 116th Street, Miami, FL 33167
(888) 684-7783 | vivex.com | customercare@vivex.com

1. Gruskin, E. et al., Demineralized bone matrix in bone repair: history and use. *Advanced Drug Delivery Reviews*, 2012. 64:1063-1077
2. Data on file at VIVEX Biologics, Inc.
3. Martin GJ Jr, Boden SD, Titus L, Scarborough NL, "New formulations of demineralized bone matrix as a more effective graft alternative in experimental posterolateral lumbar spine arthrodesis.", *Spine*. 1999 Apr 1;24(7):637-45.
4. Best, Benjamin. P. Cryoprotectant Toxicity: Facts, Issues, and Questions. *Rejuvenation Research*, 2015. Vol. 18, No. 5.
5. Renzi, S., et al., Mesenchymal stromal cell cryopreservation. *Biopreservation and Biobanking*, 2012. 10(3): p. 276-281.
6. Asghar, W., et al., Preserving human cells for regenerative, reproductive, and transfusion medicine. *Biotechnology Journal*, 2014. 9: p. 895-903.