



Gorilla® Lapidus Medial Wall Plating System

PRODUCT INFORMATION



LAPIDUS PLATE FEATURES

- Medial wall fixation allows for improved resistance to bending compared to dorsal and dorsal medial plates due to the placement of the plate on its side¹
- Ramped proximal portion of plate designed to help reduce irritation to the tibialis anterior
- Patented plate shape with plantar arm helps to reduce plantar gapping (excludes Armless Standard plate)
- Medial wall plate curvature aids in limiting adductory forces on 1st metatarsal by supporting corrected position during distal screw insertion
- Anatomically contoured to 1st metatarsal and medial cuneiform

PRECISION™ GUIDE LAPIDUS

- Patent pending guide allows for placement of cross-screw across the arthrodesis while maintaining on-axis placement of plate screws
- Provides four trajectories of guide wire paths for variations among patient anatomy
- Precision™ Guide Lapidus allows for 3.5 mm, 4.0 mm or 4.5 mm cannulated cross-screw

	Standard Armless	Standard	Medial Wall Step-Off	Graft Spanning
Size	Standard	Standard	Step off sizes in 1 mm increments from 1-5 mm	Small (5 mm Graft) Medium (8 mm/10 mm Graft) Large (12 mm/14 mm Graft)
Plate Thickness	1.3 mm	1.3 mm	1.4 mm	1.6 mm
Total Plates	2	2	10	6

LAPIDUS PLATES (20 Total Plates)

Gorilla® Lapidus Plating System

ADDITIONAL INFORMATION

PRESERVE™ LAPIDUS GRAFT

- Patented shape features both dorsal to plantar and medial to lateral taper allowing for biplanar correction
- Donor harvest site is density matched specific to Lapidus indication for strength demands and blood flow requirements
- Aseptically processed without gamma irradiation or hydrogen peroxide to help preserve the native mechanical advantages of human bone and the osteoinductivity of the environment in which the graft is being implanted





Lapidus Trial Sizers

12 mm

10 mm

8 mm

5 mm

- Allows surgeon to determine ideal graft size by demonstrating correction

LAPIDUS NIPPER



- Sharp Toothed Jaw
- Patent pending instrument specifically designed to aid in removal of the two bone fragments created after sagital saw cartilage resection of the 1st TMT Joint
- · Osteotome jaw designed to aid in completion of saw cut
- · Sharp toothed jaw helps to release remaining soft tissue attachments
- Osteotome jaw and sharp toothed jaw clamp together allowing for extraction of fragment with less disruption of surrounding tissue
- Long jaws designed to grasp around entire bone fragment from dorsal to plantar

LAPD-01 Rev D

[™]Trademarks and [®] Registered Trademarks of Paragon 28, Inc. © Copyright 2020 Paragon 28, Inc. All rights reserved. Patents: www.paragon28.com/index.php/patents

For the contraindications, potential complications and adverse reactions, warnings and precautions associated with this device, please refer to the device specific instructions for use at http://www.paragon28.com/fus

Paragon 28, Inc. 14445 Grasslands Dr. Englewood CO80112 (855) 786-2828



Paragon28 Medical Devices Trading Limited First FloorBlock7 Beckett Way Park WestBusiness Park Dublin12 D12 X884 Ireland +353 (0) 1588 0350



Lapidus Nipper inserted into the 1st TMT joint space

