

Exclusively foot & ankle ²⁰
Paragon®



The mission of Paragon 28®

Our mission is to strategically build a company around the core principles that **drive innovation and quality**. Relentlessly working to advance the science behind foot and ankle surgery, Paragon 28 is passionate about and committed to:



Blending different surgical philosophies from various global thought leaders to develop biomechanically and clinically relevant surgical solutions.



A customer service based, dedicated, and highly-trained distribution network to the foot and ankle market.



Meeting the needs of increasing pricing pressures and reimbursement while continuing to develop cutting-edge surgical implants.



Creating complete surgeon-centric systems, specialty instruments, and next generation implants to solve real world issues faced by foot and ankle surgeons.



Improving clinical outcomes through meaningful and unbiased research. We are ready to embrace research data regardless of marketing impact.

Paragon 28® History

In 2010, Paragon 28® was established as an orthopedic foot and ankle company. The name “Paragon 28” was chosen to show that we are exclusively a foot and ankle company, with the “28” representing the number of bones in the foot. We will remain true to that vision. Paragon 28 was started as a small, family-based company and we have kept those core ideals as we have grown.

Our first product was the Monster® Screw System, followed by the Gorilla® Plating System. As we continue to add more products, we maintain a fine level of detail to every feature of every product we design. The goal is to give options to the foot and ankle surgeon, knowing there is more than one way to achieve a great result. We have listened to surgeons and will continue to do so in order to understand and address their specific needs.



SOFT TISSUE


TenoTac™ 2.0



Grappler® Interference Screw System




R3LEASE™ Stabilization System




R3ACT™ Stabilization System




Grappler® Suture Anchor System—All Suture



Grappler® Suture Anchor System—Titanium



Grappler® Knotless Anchor System





Grappler® Suture Anchor System— PEEK



Mister Tendon™ Harvester System



Paratrooper™ Plantar Plate Repair System



BIOLOGICS

V92*^{*}




V92-FC+^{TM*}




Bone Marrow Aspirate KitTM



Bone Graft Harvest SystemTM




Paraderm[®] Dermal Matrix*



Pro3[®]-P Amniotic Placental Membrane*



Pro3[®]-C Amniotic Placental Membrane*



MgNumTM BVF*




BEASTTM Injectable Putty Demineralized Bone Matrix*




BEAST+ Demineralized Bone Matrix*



BEASTTM Osteologic Sponge*



PRESERVETM Bone Wedges*




FOREFOOT

HammerTube® System



Baby Gorilla® Akin Plates




Phantom® Metatarsal Shortening System



JAWS® Nitinol Staple System




Monster® BITE Snap-Off Screws




Gorilla® MTP Plates



PRESERVE® MTP Graft*



AVITRAC®**



FJ2000™ Power Console System**



Wedge


Corta

Recta



FOREFOOT (Cont.)


PRECISION® Jones
Fracture Screw System



Baby Gorilla® 5th Metatarsal
Plating System



Gorilla® BOW Plates




PROMO™ Proximal Rotational
Metatarsal Osteotomy



Mini-Monster® Solid Screws
2.7/3.5/4.0



Mini-Monster® Canulated Screws 2.0/2.5/3.0/3.5/4.0



Precision® MIS Bunion System


Chamfer Head
Sits flush to bone

Sharp Tip and Sharp Threads
Four sharp tips are self-drilling

Forward Cutting Flutes
Creates thread pattern during
screw insertion and allows the
screw to be self-tapping



Modified Hexalobe Head
Driver tip sits in one position and
allows for orientation of screw




MIDFOOT


Gorilla® NC Fusion
Plating System



Gorilla® Cotton Plates




Gorilla® Universal Plates



PRESERVE™ Lapidus Grafts




Phantom® Small Bone
Intramedullary Nail





Phantom® MIS Nail System



Gorilla® Lapidus Plating System



Lapidus Cut Guide System




Bun-Yo-Matic™ Lapidus Clamp

MIDFOOT (Cont.)


JAWS™ Great White
Staple System



Gorilla® Medial Column
Straddle Plates



Gorilla® Medial Column
Plating System



Gorilla® Lateral Column
Fusion Plating System



Titan 3-D™
Cotton Metal Wedges



PRESERVE™ Cotton Grafts




Gorilla® Lisfranc Slated T Plates




Gorilla® LisFranc Dual Ray Plates



Gorilla® Central Column
Plating System




Baby Gorilla® Navicular Plates



Baby Gorilla® Zig Zag Plates



Baby Gorilla® Cuboid Plates



Baby Gorilla® Y-Plates




HINDFOOT


Silverback™ Mini-Open Ankle Fusion Plating System




Silverback™ Span Plating System



Silverback™ Straddle Plating System



Silverback™ Posterior Ankle Fusion Plating System



Phantom® ActivCore Nail System



Phantom® TTC/TC Nail System




Phantom® Hindfoot Nail Precision® Guide



HINDFOOT (Cont.)

Silverback™ Lateral Ankle Fusion Plating System




Silverback™ Anterior Ankle Fusion Plating System




Gorilla® BOW & ARROW Evans Plate




Titan 3-D™ Evans Metal Wedge




Gorilla® Tibia Plates




Gorilla® Ankle Fracture Plates




Gorilla® Calc Fracture Plates




Gorilla® Teddy Bear Plates



Baby Gorilla® Talar Neck Plates




PRESERVE™ STJDA Grafts



PRESERVE™ Cal-Cuboid Grafts



PRESERVE™ Evans Graft and Plates




Gorilla® Pilon Plates




Gorilla® Anatomic Fibular Plates




Gorilla® Hook Plates



Gorilla® Calc Slide Plates



Monster® Hindfoot Screws
4.5/5.5/7.0



APEX 3D™ Total Ankle Replacement




MAVEN™ Patient Specific Instrumentation




HINDFOOT (Cont.)

FasTrac™ Distal Tibia Alignment System




Patient Specific Talus Spacer*



JAWS™ Great White Staple System



Joust™ Beaming Screw System



Monkey Rings™ External Fixation System



Monkey Bars™ Pin to Bar External Fixation System



*Humanitarian Device. Authorized by Federal law for use in the treatment of avascular necrosis of the ankle joint in adult patients. The effectiveness of this device for this use has not been demonstrated

Innovative Products
QUALITY SOLUTIONS

PRODUCT PORTFOLIO

Exclusively foot & ankle **28**
Paragon[®]

www.Paragon28.com

P28SS-01 Rev G
2024-04-26

Paragon 28[®], Inc. 
14445 Grasslands Dr.
Englewood, CO 80112 USA
(855) 786-2828

Paragon 28[®] Medical Devices Trading Limited
First Floor Block 7 Beckett Way
Park West Business Park
Dublin 12
D12 X884
Ireland
+353 (0) 1588 0350

For the contraindications, potential complications and adverse reactions, warnings and precautions associated with this device, please refer to the device specific instructions for use at <http://www.paragon28.com/ifus>