


JuggerKnot[®] All-Suture Anchor Portfolio

Ordering Information Guide



OsseoCoat™ JuggerKnot All-Suture Anchors

Ordering Information



A JuggerKnot Soft Anchor with OsseoCoat Technology is identical to a JuggerKnot Soft Anchor except for the presence of bioceramics embedded in the anchor portion of the implant. The hydroxyapatite and bioceramic particles, embedded in the surface of the suture anchor material, provide a surface coating comprised of bioactive materials shown to have osteoconductive properties.^{1,2}


The OC JuggerKnot Anchors utilize the existing and corresponding instrumentation within the JuggerKnot portfolio. Please refer to the following pages for all instrumentation.

1.45 mm OC JuggerKnot Rigid Implants

Product	Drill Size	Insert Length	Suture Size	# Sutures	Sliding vs. Non-Sliding	Configuration and Colors	Needle	Qty	Part Number
	1.45 mm	Rigid	1.5 mm BB	1	Sliding	Blue/Black	MO-6	Single	110048390

2.9 mm OC JuggerKnot Rigid Implants

Product	Drill Size	Insert Length	Suture Size	# Sutures	Sliding vs. Non-Sliding	Configuration and Colors	Needle	Qty	Part Number
	2.9 mm	Rigid	#2 MB	2	Sliding	White/Blue, Black/White	MO-6	Single	110048372
	2.9 mm	Rigid	1.5mm BB	2	Sliding	White/Blue, Black/White	MO-6	Single	110048374
	2.9 mm	Rigid	1.5mm BB	2	Non-Sliding	Blue/Black SL, Black NS	MO-6	Single	110048378
	2.9 mm	Rigid	1.5mm BB	3	Sliding	White/Black, Blue/White, Black/White	MO-6	Single	110048376



OsseoCoat coating
on JuggerKnot implant sleeve

1. Dewi AH, Ana ID. The use of hydroxyapatite bone substitute grafting for alveolar ridge preservation, sinus augmentation, and periodontal bone defect: A systematic review. Heliyon. 2018 Nov 2;4(10):e00884. doi: 10.1016/j.heliyon.2018.e00884. PMID: 30417149; PMCID: PMC6218667.
2. Cheah CW, Al-Namnam NM, Lau MN, Lim GS, Raman R, Fairbairn P, Ngeow WC. Synthetic Material for Bone, Periodontal, and Dental Tissue Regeneration: Where Are We Now, and Where Are We Heading Next? Materials (Basel). 2021 Oct 15;14(20):6123. doi: 10.3390/ma14206123. PMID: 34683712; PMCID: PMC8537464.

OsseoCoat Juggerknot has not been clinically evaluated for osteoconductivity results.

JuggerKnot® All-Suture Anchors

Ordering Information

Short and Rigid JuggerKnot Systems

1.0 mm JuggerKnot Mini Implants

Product	Drill Size	Insert Length	Suture Size	# Sutures	Sliding vs. Non-Sliding	Configuration and Colors	Needle	Qty	Part Number
	1.0 mm	Mini	2-0 MB	1	Sliding	White/Blue	RB-1	Single*	912076
	1.0 mm	Mini	2-0 MB	1	Sliding	White/Blue	RB-1	10 Pack*	912080
	1.0 mm	Mini	3-0 MB	1	Sliding	White/Blue	RB-1	Single*	912082
	1.0 mm	Mini	3-0 MB	1	Sliding	White/Blue	RB-1	10 Pack*	912084

1.0 mm JuggerKnot Mini Instrumentation

Product	Description	Part Number
	1.0 mm Mini Step Drill Bit	912077

1.4 mm JuggerKnot Short Implants

Product	Drill Size	Insert Length	Suture Size	# Sutures	Sliding vs. Non-Sliding	Configuration and Colors	Needle	Qty	Part Number
	1.4 mm	Short	#1 MB	1	Sliding	Blue/White	SH	Single	912068
	1.4 mm	Short	#1 MB	1	Sliding	Blue/White	SH	10 Pack	912069
	1.4 mm	Short	1.0 mm BB	1	Sliding	White/Blue	CT-2	Single	CM-99114BN
	1.4 mm	Short	1.0 mm BB	1	Sliding	White/Blue	CT-2	10 Pack	CM-99114BNT

1.4 mm JuggerKnot Short Instrumentation

Product	Description	Part Number
	1.4 mm Short Drill Bit	912071
	1.4 mm Short Disposable Kit (Guide, Drill Bit)	912073
	1.4 mm Short Reusable Guide	912072

*Includes Drill Bit

Short and Rigid JuggerKnot Systems (cont.)

1.45 mm JuggerKnot Rigid Implants

Product	Drill Size	Insertor Length	Suture Size	# Sutures	Sliding vs. Non-Sliding	Configuration and Colors	Needle	Qty	Part Number
	1.45 mm	Rigid	#1 MB	1	Sliding	Blue/White	SH	Single*	110005307
	1.45 mm	Rigid	#1 MB	1	Sliding	Blue/White	SH	10 Pack*	110005319
	1.45 mm	Rigid	#2 MB	1	Sliding	Blue/White	MO-6	Single*	110005315
	1.45 mm	Rigid	#2 MB	1	Sliding	Blue/White	MO-6	10 Pack*	110005320
	1.45 mm	Rigid	1.5 mm BB 1.5 mm BB	1	Sliding	Blue/Black	MO-6	Single	CM-99145BRN
	1.45 mm	Rigid		1	Sliding	Blue/Black	MO-6	10 Pack	110046069

1.45 mm JuggerKnot Rigid Instrumentation

Product	Description	Part Number
	1.45 mm Rigid Drill Bits (5 pack)	110005323

2.9 mm JuggerKnot Rigid Implants

Product	Drill Size	Insertor Length	Suture Size	# Sutures	Sliding vs. Non-Sliding	Configuration and Colors	Needle	Qty	Part Number
	2.9 mm	Rigid	#2 MB	2	Sliding	Blue/White, Blue	FSL	Single	110005093
	2.9 mm	Rigid	#2 MB	2	Sliding	Blue/White, Blue	FSL	10 Pack	110005241
	2.9 mm	Rigid	#2 MB	2	Sliding	Blue/White, Blue	MO-6	Single	110005096
	2.9 mm	Rigid	#2 MB	2	Sliding	Blue/White, Blue	MO-6	10 Pack	110005242
	2.9 mm	Rigid	1.5 mm BB	2	Sliding	White/Blue, Black/White	MO-6	Single	CM-99129BN
	2.9 mm	Rigid	1.5 mm BB	2	Sliding	White/Blue, Black/White	MO-6	10 Pack	CM-99129BNT
	2.9 mm	Rigid	1.5 mm BB	3	Sliding	White/Black, Blue/White, Black/White	MO-6	Single	110046088
	2.9 mm	Rigid	1.5 mm BB	3	Sliding	White/Black, Blue/White, Black/White	MO-6	10 Pack	110046089
	2.9 mm	Rigid	1.5 mm BB	2	Non-Sliding	Blue/Black SL, Black NS	MO-6	Single	CM-99229BN
	2.9 mm	Rigid	1.5 mm BB	2	Non-Sliding	Blue/Black SL, Black NS	MO-6	10 Pack	CM-99229BNT

*Includes Drill Bit

Short and Rigid JuggerKnot Systems (cont.)

2.9 mm JuggerKnot Short Instrumentation

Product	Description	Part Number
	2.9 mm Disposable Drill Bit Short	110005306
	2.9 mm Reusable Short Guide	110016951

JuggerKnot Implant Needle Chart

Product	RB-1	SH	CT-2	MO-6	FSL
Needle Geometry	1/2 Circle	1/2 Circle	1/2 Circle	1/2 Circle	3/8 Circle
Needle Point	Taper Point	Taper Point	Taper Point	Taper Point	Reverse Cutting
Needle Length	17 mm	26 mm	26 mm	26 mm	30 mm
Needle Outer Diameter	0.61 mm	0.86 mm	0.99 mm	1.28 mm	1.28 mm

Products within this system are under the design control of various legal manufacturers. Refer to the product labeling of each device for the legal manufacturer. All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. This material is intended for health care professionals. For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the instructions for use or contact your local representative; visit www.zimmerbiomet.com for additional product information.

©2023, 2024, 2026 Zimmer Biomet



Legal Manufacturer

Riverpoint Medical
825 NE 25th Avenue
Portland, OR 97232
USA



Legal Manufacturer

Biomet Sports Medicine
56 E. Bell Drive, P.O. Box 587
Warsaw, IN 46581-0587
USA

Sports Medicine Customer Service:

1.800.348.2759
customerservice@zimmerbiomet.com

www.zimmerbiomet.com