

GENIE™ & GENIESP™ KNOTLESS IMPLANTS

"Soft Tissue Repair Granted"



The GENIE & GENIESP Biocomposite and PEEK Anchors offer a **Reproducible** and **Versatile** Knotless solution with secure soft tissue to bone fixation, self-punching ease, and a fully **Biocomposite option** accommodating various soft tissue repairs.

Dual pitch thread for Cancellous and Cortical fixation, available in multiple sizes and in three different materials:

- 1. All Biocomposite:** Anchor body and suture eyelet are composed of bioabsorbable biocomposite material.
- 2. All PEEK:** Anchor body and suture eyelet are composed of PEEK self-punching tip material.
- 3. Hybrid:** Biocomposite anchor body and a PEEK self-punching tip eyelet.



- Complete biocomposite threaded knotless implant with a non-self-punching biocomposite tip
- Vented anchor body
- Biocomposite threaded knotless implant with PEEK self-punching tip
- Vented anchor body
- PEEK threaded knotless implant with PEEK self-punching tip
- Vented anchor body

GENIE KNOTLESS IMPLANTS

Part Number	Description
375BCP	3.75 MM GENIE BIOCOMPOSITE ANCHOR W/PEEK TIP
375PP	3.75 MM GENIE PEEK ANCHOR W/PEEK TIP
475BCBC	4.75 MM GENIE BIOCOMPOSITE ANCHOR W/BIOCOMPOSITE TIP
550BCBC	5.5 MM GENIE BIOCOMPOSITE ANCHOR W/BIOCOMPOSITE TIP
475BCPSP	4.75 MM GENIESP BIOCOMPOSITE ANCHOR W/PEEK SELF-PUNCHING TIP
550BCPSP	5.5 MM GENIESP BIOCOMPOSITE ANCHOR W/PEEK SELF-PUNCHING TIP
650BCPSP	6.5 MM GENIESP BIOCOMPOSITE ANCHOR W/PEEK SELF-PUNCHING TIP
475PPSP	4.75 MM GENIESP PEEK ANCHOR W/PEEK SELF-PUNCHING TIP
550PPSP	5.5 MM GENIESP PEEK ANCHOR W/PEEK SELF-PUNCHING TIP
650PPSP	6.5 MM GENIESP PEEK ANCHOR W/PEEK SELF-PUNCHING TIP

GENIE KNOTLESS INSTRUMENTATION

Part Number	Description
375TAP	3.75 MM GENIE BONE TAP
475AWL	4.75 MM GENIE BONE AWL
475TAP	4.75 MM GENIE BONE TAP
375BKT	5.5 MM GENIE BONE TAP
375PKT	GENIE BONE DRILL GUIDE

GENIE KNOTLESS DISPOSABLE

Part Number	Description
475AD	4.75 MM DISPOSABLE AWL
375DRILL	3.75 MM GENIE BONE DRILL
475DRILL	4.75 MM GENIE BONE DRILL
550DRILL	5.5 MM GENIE BONE DRILL
650DRILL	6.5 MM GENIE BONE DRILL

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

For US distribution.

Intended for Healthcare Professionals. For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the Instructions for use or contact your local representative; visit www.zimmerbiomet.com for additional product information.

©2025 Zimmer Biomet



Legal Manufacturer
Osteonic Co., Ltd
505-2, 505-3, 1021, 1206, 1207 Ho
38, Digital-ro, 29-gil, Guru-gu
Seoul, Republic of Korea

www.osteonic.com

3781.2-US-en Issue Date 2025-05